


**COSMOSIL**



High Performance  
Liquid Chromatography

## CORPORATE PROFILE

Nacalai Tesque dates back to 1846 when the company's founder Mansuke Nakarai opened Nakarai Mansuke Shoten, Ltd., an apothecary selling traditional Japanese and Chinese medicines. In 1958, this company's reagent department became an independent company, Nakarai Chemicals, Ltd.

The company has since dedicated itself to expanding its corporate base and has strived to be an enterprise that our customers always rely on, while taking pride in its contribution to scientific and industrial development.


Making the most of this 30-year history and as a step toward the future, Nakarai Chemicals changed its corporate name to Nacalai Tesque, Inc. in 1988. At Nacalai Tesque, we have fostered a corporate commitment to the pursuit of reliable quality and the creation of products of real value, while serving as a vital link between humanity and science. Centering around research chemicals, the fields of our activities include fine chemicals, diagnostics and related laboratory equipment and supplies.

The pace of scientific and technological progress in every industrial field is rapidly accelerating, and all business partners and affiliates are required to provide even more diversified and advanced expertise.

It is our corporate policy to strive for our lofty ideals for excellence while respecting our long history and tradition. We consider it our mission to maintain close contact with our customers by offering reliable quality in all our products, information and services, and by making full use of the knowledge and experience of our staff.

### The 30<sup>th</sup> Anniversary of COSMOSIL HPLC Column

In 2010, Nacalai Tesque reached another important milestone in the company's history with the 30<sup>th</sup> anniversary of COSMOSIL brand HPLC column. For the past 30 years, Nacalai Tesque has provided innovative and reliable chromatography products to the analytical industry. Our goal has not changed since the first COSMOSIL HPLC column was developed in 1980; it is to focus on customer needs and make customer's separations as successful as possible by offering the highest level of products and technical support. Keeping this goal in mind, we look forward to serving the industry for the next 30 years and beyond.



CORPORATE PROFILE .....	1
CONTENTS .....	2-3
GENERAL INFORMATION .....	4

## I. COSMOSIL HPLC COLUMN

1. COSMOSIL/COSMOGEL .....	6
2. Column selection guide .....	7
3. COSMOSIL silica packing materials .....	8-9
4. Performance guarantee .....	10-11
1) Quality guarantee of packing materials .....	10-11
2) Quality guarantee of COSMOSIL packed columns .....	11
5. Reversed phase chromatography .....	12-35
1) Octadecyl types .....	12-19
C <sub>18</sub> -MS-II .....	14-15
C <sub>18</sub> -AR-II .....	16-17
C <sub>18</sub> -PAQ .....	18-19
2) Special columns .....	20-29
Cholester .....	21-23
πNAP .....	24-25
PYE .....	26-27
NPE .....	28
PBB-R .....	29
3) Other columns .....	30-32
CN-MS .....	31
C <sub>22</sub> -AR-II, C <sub>8</sub> -MS, C <sub>4</sub> -MS, TMS-MS, PE-MS .....	32
4) Silica based preparative columns .....	33
15C <sub>18</sub> -MS-II, 15C <sub>18</sub> -AR-II, 15C <sub>18</sub> -PAQ .....	33
5) Ultra-Fast Liquid Chromatography (UFLC) .....	34-35
2.5C <sub>18</sub> -MS-II, 2.5Cholester .....	34-35
6. Normal phase chromatography .....	36-37
SL-II .....	36-37
7. Hydrophilic interaction chromatography .....	38-39
HILIC .....	38-39
8. Saccharide analysis .....	40-43
Sugar-D .....	41-43
NH <sub>2</sub> -MS .....	43
9. Protein separation wide pore columns .....	44-53
1) Reversed phase chromatography .....	44-47
Protein-R .....	44-45
C <sub>18</sub> -AR-300, C <sub>8</sub> -AR-300, C <sub>4</sub> -AR-300, Ph-AR-300 .....	46-47
2) Gel filtration chromatography .....	48-49
Diol-120-II, Diol-300-II .....	48-49
3) Ion-exchange chromatography .....	50-52
DEAE, QA, CM, SP .....	50-52
4) Hydrophobic interaction chromatography .....	53
HIC .....	53
10. Special columns for fullerenes .....	54-59
Buckyprep .....	56
Buckyprep-M .....	57
PBB .....	58
PYE .....	59
NPE .....	59
11. Special columns for carbon nanotubes .....	60
CNT-300, CNT-1000, CNT-2000 .....	60
12. Conventional columns versus high performance columns .....	61

## II. LIQUID CHROMATOGRAPHY RELATED PRODUCT

1. Packing materials for column chromatography .....	64-71
C <sub>18</sub> -OPN .....	65-67
C <sub>18</sub> -PREP .....	68
SL-II-PREP .....	69
Silica Gel (spherical • neutral) .....	70
Silica Gel (for column chromatograph) .....	71

2. Liquid chromatography related products (for mobile phase) .....	72
Ion-pair reagents .....	72-73
3. Liquid chromatography related products (for pretreatment) .....	74-75
Cosmonice filter .....	74-75
Cosmospin filter .....	74-75
COSMOSIL HPLC accessories .....	76

## III. APPLICATION DATA

1. COSMOSIL Chromatogram Index .....	78
2. Application data of substances in Japanese Pharmacopoeia, 15 <sup>th</sup> version .....	79
3. COSMOSIL Application Data .....	80-149
1) Drugs .....	80-112
2) Crude Drugs .....	113-117
3) Natural Compounds .....	118-120
4) Pesticides .....	121-122
5) Food Additives .....	123-124
6) Vitamins .....	125-127
7) Metabolites .....	128-129
8) Carbohydrates .....	130-132
9) Lipids .....	133-134
10) Nucleic Acid Related Substances .....	135
11) Amino Acids, Peptides and Proteins .....	136-142
12) The others .....	143-149
4. Reference list .....	150-151

## IV. TECHNICAL NOTE

1. Selectivity of packing materials in reversed phase liquid chromatography .....	154-160
2. Preparation of mobile phase for HPLC .....	161-163
3. Sample pretreatment for HPLC .....	164-167
4. Baseline noise in gradient elution .....	168
5. Troubleshooting for increased pressure .....	169-171
6. Effect of guard column .....	172
7. Troubleshooting for normal phase chromatography .....	173-174
8. Inner diameter of column (scale down and scale up) .....	175

## V. INDEX

1. Sample name .....	178-187
2. Packing material's name .....	188-202



## General information

### General ordering information

When placing an order with us or making an enquiry, please contact our International Business Development Group or your local distributor. Please clearly identify the product in question when submitting your enquiry.

The speed of innovation is accelerating. We always have brand new or improved columns not listed here. There are also many other products Nacalai Tesque can supply. Therefore we urge you to make enquiries.

### Product description and availability

Please note that the product specifications are subjected to changes and the manufacturing of some product may be stopped. Please consult the table on page 11 for cross-reference information on old products and their newer and better equivalents.

### Column identification

At the end of each section, the COSMOSIL and COSMOGEL packed columns are listed in a way that the particle size, stationary phase, column size of the packing material can be easily determined.

38019-81 COSMOSIL 5 C<sub>18</sub>-MS-II 4.6 mm I.D. x 150 mm

(1)        (2)        (3)        (4)        (5)

When placing an order, please clearly indicate the product number (1), product name (2), particle size (3), type of stationary phase (4) and column size (5).

### Warranty claims

The manufacturer will replace defective columns if notified within 2 weeks of receipt of the product by the customer under the following conditions:

- 1) Column abnormalities are due to accidents in shipping or rough handling.
- 2) The number of effective plates of the column is considerably lower than the minimum guaranteed theoretical plate number documented in the inspection report that accompanies each column.

Please contact the International Business Development Section of Nacalai Tesque or your local distributor for additional information.

### Terms and conditions of sale

Terms are subject to conditions set forth by the authorized Nacalai Tesque distributors in each country.

### Not for clinical use

Nacalai Tesque products are not intended for clinical use. While clinical applications may be shown, these products are not validated for clinical use.

1. COSMOSIL/COSMOGEL .....	6
2. Column selection guide .....	7
3. COSMOSIL silica packing materials .....	8
4. Performance guarantee .....	10
1) Quality guarantee of packing materials .....	10
2) Quality guarantee of COSMOSIL packed columns .....	11
5. Reversed phase chromatography .....	12
1) Octadecyl types .....	12
C <sub>18</sub> -MS-II .....	14
C <sub>18</sub> -AR-II .....	16
C <sub>18</sub> -PAQ.....	18
2) Special columns.....	20
Cholester .....	21
πNAP .....	24
PYE .....	26
NPE .....	28
PBB-R .....	29
3) Other columns .....	30
CN-MS .....	31
C <sub>22</sub> -AR-II, C <sub>8</sub> -MS, C <sub>4</sub> -MS, TMS-MS, PE-MS .....	32
4) Silica based preparative columns .....	33
15C <sub>18</sub> -MS-II, 15C <sub>18</sub> -AR-II, 15C <sub>18</sub> -PAQ .....	33
5) Ultra-Fast Liquid Chromatography (UFLC) .....	34
2.5C <sub>18</sub> -MS-II, 2.5Cholester.....	34
6. Normal phase chromatography .....	36
SL-II .....	36
7. Hydrophilic interaction chromatography .....	38
HILIC .....	38
8. Saccharide analysis .....	40
Sugar-D.....	41
NH <sub>2</sub> -MS.....	43
9. Protein separation wide pore columns .....	44
1) Reversed phase chromatography .....	44
Protein-R .....	44
C <sub>18</sub> -AR-300, C <sub>8</sub> -AR-300, C <sub>4</sub> -AR-300, Ph-AR-300 .....	46
2) Gel filtration chromatography .....	48
Diol-120-II, Diol-300-II .....	48
3) Ion-exchange chromatography .....	50
DEAE, QA, CM, SP.....	50
4) Hydrophobic interaction chromatography .....	53
HIC.....	53
10. Special columns for fullerenes .....	54
Buckyprep.....	56
Buckyprep-M .....	57
PBB .....	58
PYE .....	59
NPE .....	59
11. Special columns for carbon nanotubes .....	60
CNT-300, CNT-1000, CNT-2000 .....	60
12. Conventional columns versus high performance columns .....	61

# 1. COSMOSIL/COSMOGEL


## General description of the COSMOSIL/COSMOGEL packing materials

Sample	Separation mode	Packing material	Stationary phase	Special features and applications	USP category	Page	
Organic compounds (low M.W.)	Reversed phase	C <sub>18</sub> -MS-II	Octadecyl group	Multi-purpose C <sub>18</sub> column. Monofunctional silylation on ultra-pure silica gel for separation of the widest range of compounds	L1	14	
		C <sub>18</sub> -AR-II		Multi-purpose C <sub>18</sub> column using ultrapure silica gel. Features strong acid resistance and suitable for a wide range of separation	L1	16	
		C <sub>18</sub> -PAQ		Reversed phase column, compatible with 100% water based mobile phases	L1	18	
		Cholester	Cholesteryl group	Usable under the same condition as C <sub>18</sub> . Unique rigid cholesteryl structure improves separation		21	
		πNAP	Naphthylethyl group	Stronger π-π interaction than phenyl column		24	
		PYE	Pyrenylethyl group	The most powerful π-π interaction		26	
		NPE	Nitrophenylethyl group	Separation utilizing π-π interaction and Dipole-dipole interaction		28	
		PBB-R	Pentabromobenzyl group	Separation utilizing dispersion force		29	
		CN-MS	Cyanopropyl group	Great reproducibility using isocratic elution mode	L10	31	
		C <sub>22</sub> -AR-II	Docosyl group	Alkyl chain columns except C <sub>18</sub> column			32
		C <sub>8</sub> -MS	Octyl group		L7		
		C <sub>4</sub> -MS	Butyl group		L26		
		TMS-MS	Trimethyl group		L13		
		PE-MS	Phenylethyl group	π-π interaction	L11		
	Normal phase	SL-II	**	Normal Phase chromatography with nonpolar organic solvents	L3	36	
	Hydrophilic interaction	HILIC	Triazole	Retains highly polar compounds that would not be retained in C <sub>18</sub> column		38	
Monosaccharides Oligosaccharides	Hydrophilic interaction	Sugar-D	Secondary/Tertiary amine	A novel stationary phase for saccharide separation. Extended column life and increased stability. Alternative to aminopropyl type		41	
		NH <sub>2</sub> -MS	Aminopropyl group	Primary amine bonded column		43	
Proteins	Reversed phase	Protein-R	Octadecyl group	The most suitable reverse phased column for proteins		44	
		C <sub>18</sub> -AR-300	Octadecyl group	Wide pore type reversed phase columns with high acid resistance recommended for the separation of proteins, polypeptides, nucleic acids and other large molecules	L1	46	
		C <sub>8</sub> -AR-300	Octyl group		L7		
		C <sub>4</sub> -AR-300	Butyl group		L26		
		Ph-AR-300	Phenyl group		L11		
	Gel permeation	Diol-120-II	Diol group	Silica-based gel filtration column for high speed separation of proteins and water soluble polymer	L20	48	
		Diol-300-II					
	Ion-exchange	DEAE	Diethylaminoethyl type	Weak anion-exchange		50	
		QA	Quarternary ammonium type	Strong anion-exchange	L23		
		CM	Carboxymethyl type	Weak cation-exchange			
SP		Sulfopropyl type	Strong cation-exchange				
	Hydrophobic interaction	HIC	**	Hydrophobic interaction chromatography column for protein separation		53	
Fullerenes	**	Buckyprep	Pyrenylpropyl group	Standard column for fullerenes separation		56	
		Buckyprep-M	Phenothiazinyl group	Designed to separate metallofullerenes		57	
		PBB	Pentabromobenzyl group	Designed for the preparative separation of fullerenes using carbon disulphide, o-dichlorobenzene and toluene		58	
		NPE	Nitrophenylethyl group	Separation of derivatized fullerenes		59	
		PYE	Pyrenylethyl group	Separation of fullerenes and structural isomer			
Carbon nanotubes	Gel permeation	CNT-300 NEW	Hydrophilic group (neutral)	Separation of soluble carbon nanotubes		60	
		CNT-1000 NEW					
		CNT-2000 NEW					

For old type columns, please refer to page 61.

# 2. Column selection guide

## COSMOSIL columns selection guide



## COSMOSIL columns selection guide

### Organic compounds (low M.W.)

Octadecyl group bonded column (C<sub>18</sub>, ODS) are recommended as first-choice columns for separations of organic compounds (low M.W.). If there is not enough separation or no retention using COSMOSIL C<sub>18</sub> columns, COSMOSIL series offer many kinds of specialty columns.

- Medicines
- Crude drugs
- Natural compounds
- Pesticides
- Food additives
- Vitamins
- Lipids etc.

### Saccharides

- COSMOSIL Sugar-D is recommended for the separation of monosaccharides and oligosaccharides as a first-choice column.
- For the separation of sugar derivatives, COSMOSIL C<sub>18</sub>-PAQ is suitable as well.

### Proteins

- Please select based on the separation mode. Please refer to page 44.

### Fullerenes

- COSMOSIL Buckyprep is most suitable for the separation of fullerenes.

# 3. COSMOSIL silica packing material

## Introduction

Superior HPLC columns can be produced only with excellent packing materials and superb packing technique. COSMOSIL columns are well known for their high efficiency and high-resolution separations. Based on spherical, totally porous silica, COSMOSIL columns provide enhanced chemical and mechanical stability as well as very high surface coverage.

The selection of the C<sub>18</sub> chemistries available enables the chromatographer to tailor separation to special applications. The ultra pure silica based MS-II series with widely extended pH range are developed for improved separation of basic compounds. The C<sub>18</sub> AR-II phase provides increased acid resistance. Four unique bonded chemistries are available for COSMOSIL specialty columns: Cholester, πNAP, PYE and HILIC. These specialty columns may improve the separation compared with conventional columns. Five highly effective phases for fullerene separation are also available: Buckyprep, Buckyprep-M, PBB, PYE and NPE. COSMOGEL packing materials are non-silica based and provide superior in ion exchange columns.






Figure. Microscopic photograph of the silica gel

## Packing material and view showing a frame format of column



## Raw material silica gel

COSMOSIL is based on ultra pure porous spherical silica gel (purity: 99.99% or higher). Low-purity silica gel contains metal impurity which may cause interference in the separation, especially for metal coordination compounds.



### Metal coordination compounds

The compounds, which have 2 or more hydroxyl groups, carboxyl groups, carbonyl groups and/or amino groups in the adjacent positions, can complex with metal impurity, which results in peak tailing.




Table. Metal content percentage (ppm)

Packing material	Metal content (ppm)						Total
	Al	Ca	Fe	Mg	Na	Ti	
Ultra pure silica gel	4	10	11	3	12	<1	41
Normal silica gel	32	450	26	88	56	134	786

## Stationary phase construction

While C<sub>18</sub> columns are most widely used in reversed phase HPLC, it is important to distinguish between two very different bonded phase formats. Monomeric type C<sub>18</sub> format incorporates the bonding of the C<sub>18</sub> alkyl chain to a single silica atom on the silica gel backbone. Monomeric type columns such as the COSMOSIL C<sub>18</sub>-MS-II and the MS series have excellent synthesis reproducibility, very good lot-to-lot reproducibility and short mobile phase equilibration times. On the other hand, the polymeric C<sub>18</sub> format incorporates a tri-functional silylation procedure whereby the octadecyl group is bonded to 2 or 3 silica atoms on the silica gel backbone. This increases silylation results in far greater column stability particularly in acidic mobile phase conditions. Stereo recognition capability is also greater than that of the monofunctional silylation type C<sub>18</sub> columns. The polymeric format is offered in the AR-II and the entire AR-300 series of COSMOSIL columns. Please refer to product descriptions and application chromatograms for selection guidance.





Figure. Diagrams of different stationary phase constructions (before end-capping treatment)

## End-capping treatment

The silanols (Si-OH groups) on the silica surface provided bonding site for stationary phases. However, part of the silanol groups remain uncapped as residual silanol groups even after the end-capping treatment, they cause peak tailing for basic compounds. COSMOSIL packing materials for reversed phase chromatography are of near-perfectly end-capped residual silanol groups.

- R-NH<sup>+</sup> (Basic compounds)
- Si-O<sup>-</sup> (Free silanol group)

Basic compound can form ionic bonds with residual free silanols. The ionic bonding causes peak tailing of basic compounds if a silica based column is not perfectly end-capped.



## Synthesis reproducibility

By using strictly selected silica gel and constant synthesis conditions, the chemically bonded type column retains a variance of the capacity factor (k') between synthetic lots of within ±10% and a variance of the separation factor (α) of within ± 5%. The figures below show in graphic form the lot inspection results of synthesized packing material (COSMOSIL 5C<sub>18</sub>-MS-II). Figure 1 shows the variance of stationary phase (octadecyl group) introduced volume which is the basic indicator of the quality of the packing material. Figure 2 shows the end-capping efficiency of the packing material. The variance among the lots is reduced to the minimum in the COSMOSIL packed columns.





Figure 1. Variance of the combining volume between silica gel and C<sub>18</sub>

Figure 2. Variance of end capping efficiency of the packing material

# 4. Performance guarantee

## 1) Quality guarantee of packing materials

The strict quality control system of Nacalai Tesque supports the customers with an individual "Inspection Report" which accompanies each and every COSMOSIL and COSMOGEL Packed Column (except guard columns) and an additional "Certificate of Analysis" for the COSMOSIL 5C<sub>18</sub>-MS-II and 5C<sub>18</sub>-AR-II (4.6 mm I.D. x 150 mm and 4.6 mm I.D. x 250 mm).


### Validated columns

Product name	Product number	Column size
COSMOSIL 5C <sub>18</sub> -MS-II	38019-81	4.6 mm I.D. x 150 mm
	38020-41	4.6 mm I.D. x 250 mm
COSMOSIL 5C <sub>18</sub> -AR-II	38144-31	4.6 mm I.D. x 150 mm
	38145-21	4.6 mm I.D. x 250 mm
COSMOSIL Cholester	05976-61	4.6 mm I.D. x 150 mm
	05977-51	4.6 mm I.D. x 250 mm
COSMOSIL HILIC	07056-51	4.6 mm I.D. x 150 mm
	07057-41	4.6 mm I.D. x 250 mm

### COSMOSIL certificate of analysis


Validate terms of the physical properties of the silica gel, the carbon content, polar selectivity, hydrophobicity, silanol capacity, steric selectivity, inactive degree to basic and chelating compounds.

Certificate of Analysis		MS-II
COSMOSIL 5C <sub>18</sub> -MS-II		GEL Lot No. 29
Base silicagel material	Specification	Results
Median Particle Size		
50% cum vol [μm]	4.3-4.6	4.5
Surface Area [m <sup>2</sup> /g]	320-350	321
Pore Volume [ml/g]	0.9-1.1	0.9
Median Pore Diameter [nm]	11.0-13.0	11.8
Carbon content [%]	15.5-17.5	15.6
Atomic Emission [ppm]		
Al	≤ 5	1.5
Fe	≤ 20	8.3
Ti	≤ 0.5	0.1
Na	≤ 20	1.7
Chromatographic Results		
TEST A		
α [k'(Methyl Benzoate)/k'(Benzene)]	≤ 0.80	0.73
α [k'(Toluene)/k'(Benzene)]	≥ 1.60	1.76
TEST B		
α [k'(Caffeine)/k'(Phenol)]	≤ 0.44	0.40
TEST C		
α [k'(N-Acetyl Procainamide) /k'(Benzyl Alcohol)]	≤ 0.60	0.6
α [N(N'-Acetyl Procainamide) /N(Benzyl Alcohol)]	≥ 0.30	0.38
TEST D		
α [k'(8-Quinolol) /k'(Propylbenzene)]	≤ 0.14	0.09
α [N(8-Quinolol) /N(Propylbenzene)]	≥ 0.15	0.17
Nacalai Tesque Inc. Kyoto, Japan		
Approved: Quality Control Dept.		
Name: <i>K. Iwaguchi</i>	Date: 2009.8.21	00-3




- 3 lot numbers
- 3 lot numbers are available.

Please contact us. (info.intl@nacalai.com)



## 2) Quality guarantee of COSMOSIL packed columns

Inspection report contains data of number of theoretical plates (N), peak asymmetry (s), capacity factor (k') and separation factor (α) and pressure.



Inspection report


# 5. Reversed phase chromatography

## 1) Octadecyl types

### Introduction

The reversed phase HPLC column is most commonly used because of the high theoretical plate number, excellent separation characteristics, reproducibility, affordable cost and ease of use. Columns packed with the octadecyl group bonded type silica gel (C<sub>18</sub>, ODS) are the most widely employed. We offer three types of octadecyl group bonded columns: COSMOSIL C<sub>18</sub>-MS-II, C<sub>18</sub>-AR-II and C<sub>18</sub>-PAQ, each of which has a different separation property.


### Material characteristics

Packing material	C <sub>18</sub> -MS-II	C <sub>18</sub> -AR-II	C <sub>18</sub> -PAQ
Silica gel	High purity porous spherical silica		
Average particle size	3, 5, 15 μm*	3, 5, 15 μm	5, 15 μm
Average pore size	approx. 120 Å		
Specific surface area	approx. 300 m <sup>2</sup> /g		
Stationary phase			
	Octadecyl group		
Bonding type	Monomeric	Polymeric	
Main interaction	Hydrophobic interaction		
End-capping treatment	Near-perfect treatment		
Carbon content	approx. 16%	approx. 17%	approx. 11%
pH range	2 ~ 10**	1.5 ~ 7.5**	2 ~ 7.5
Feature	This phase is recommended for most applications but particularly effective for basic compounds.	This phase is recommended for separations requiring acidic mobile phase conditions. It also shows superior molecular shape selectivity to monomeric type C <sub>18</sub> columns.	This phase is designed to offer superior retention of polar compounds and excellent reproducibility in highly aqueous mobile phases, even in 100% aqueous.

\*For 2.5C<sub>18</sub>-MS-II, please refer to page 34.

\*\*Optimum pH range of columns based on silica gel is between 2 and 7.5.

### Difference of separation property



### Column selection based on application data

We prepare following application data to help you select separation conditions.

#### • COSMOSIL Chromatogram Index

More than 6,100 single compound elution profiles with full chromatographic condition description are available. They are not only an incredible help for chromatographers, but also can be used as references in choosing conditions for similar compounds. These data are available at our web site: <http://www.nacalai.com>

#### • COSMOSIL Application Data


COSMOSIL Application Data is now available on our website. The online version includes more than 1,000 application data using COSMOSIL columns. The online data are searchable by name of sample and column. If you have any questions regarding the application data or separations of compounds not listed here, please feel free to e-mail us at [info.intl@nacalai.com](mailto:info.intl@nacalai.com).

#### • Application data of substances in Japanese Pharmacopoeia, 15<sup>th</sup> version (246 data)

We prepare data of drugs using three kinds of C<sub>18</sub> columns that are specified in HPLC analysis in application data of substances in Japanese Pharmacopoeia, 15<sup>th</sup> version. The data are available at our web site. <http://www.nacalai.co.jp/en/cosmosil/TheJP15.htm>, or type "COSMOSIL Japanese Pharmacopoeia" at a search site.


### Column selection by mobile phase

- If a mobile phase is determined, use the following chart to select an appropriate COSMOSIL column.
- Refer to application data above for choosing a mobile phase of new analysis.
- Adjustment of pH is required for dissociative compounds.
- Generally acidic mobile phase is suitable for acidic compounds, and neutral mobile phase is suitable for basic compounds.
- If you are not sure about the mobile phase, try C<sub>18</sub>-MS-II first.




# C<sub>18</sub>-MS-II

The COSMOSIL 5C<sub>18</sub>-MS-II is a well-balanced column with better basic performance such as sharper peaks of basic compounds and chelating compounds, and high theoretical plate number. It is the most popular HPLC column because we produce very consistent products and minimize variation from lot-to-lot. Furthermore we provide abundant application data with the column, so it will help you to choose an analysis condition of your sample.



## Analysis of basic compounds and metal coordination compounds

The COSMOSIL 5C<sub>18</sub>-MS-II column, taking advantage of a new end-capping treatment, can replace the original COSMOSIL C<sub>18</sub> (ODS) column. A new end-capping treatment with polar groups for "shield effect" has significantly improved peak shape for basic compounds. Ultra pure silica gel with low trace-metal content is used for COSMOSIL columns; thus the columns provide excellent peak shapes for chelating compounds.




## Validation


The strict quality control system of Nacalai Tesque supports the customers with an individual "Inspection Report" which accompanies each and every COSMOSIL and COSMOGEL Packed Column (except guard columns) and an additional "Certificate of Analysis" for the COSMOSIL 5C<sub>18</sub>-MS-II (4.6 mm I.D. x 150 mm and 4.6 mm I.D. x 250 mm). For more informations, please refer to page 10.

## Application data

### Parasympatholytic agents



### Analgesic antipyretic drugs



## Ordering information

### Analytical / Preparative column (Particle size: 5 μm)

#### COSMOSIL 5C<sub>18</sub>-MS-II Packed Column

Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
1.0× 50	02824-31	4.6×100	38018-91
1.0×150	02896-01	4.6×150	38019-81
2.0× 30	05876-71	4.6×150 3 lots set*	09397-73
2.0× 50	04355-21	4.6×250	38020-41
2.0×100	05597-31	6.0×150	38021-31
2.0×150	38025-91	6.0×250	38022-21
2.0×250	05761-61	10× 50	05789-21
3.0×100	05458-51	10×150	34355-91
3.0×150	34245-31	10×250	38023-11
3.0×250	34254-11	20×150	05091-41
4.6× 30	34341-61	20×250	38024-01
4.6× 50	38017-01	28×250	05760-71

#### COSMOSIL 5C<sub>18</sub>-MS-II Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	38014-31
4.6×10 Cartridge**	38015-89
10×20	38016-11
20×20	05790-81
20×50	34371-71
28×50	34347-01

\* For 4.6×150 3 lots set, please refer to page 11.  
 \*\* 3 cartridges included, needs a holder.  
 Please refer to page 76.

### Preparative column (Particle size: 15 μm)

#### COSMOSIL 15C<sub>18</sub>-MS-II Packed Column

Column size I.D. x length (mm)	Product number
28×250	34525-61
50×250	05886-41
50×500	34531-71

#### COSMOSIL 15C<sub>18</sub>-MS-II Guard Column

Column size I.D. x length (mm)	Product number
28×50	05885-51
50×50	34527-41

### Fast LC column (Particle size: 3 μm)

#### COSMOSIL 3C<sub>18</sub>-MS-II Packed Column

Column size I.D. x length (mm)	Product number
2.0× 50	05514-01
4.6× 10	38065-71
4.6× 50	38066-61
4.6×100	38067-51

For more informations, please refer to page 33 for 15C<sub>18</sub>-MS-II and page 34 for 2.5C<sub>18</sub>-MS-II.  
 For flow rate and device of semi-micro columns, or preparative columns, please refer to page 175.

# C<sub>18</sub>-AR-II

The COSMOSIL 5C<sub>18</sub>-AR-II packed column features a polymeric type of C<sub>18</sub> reversed phase material. The column employs an ultra-pure silica gel low in metal impurities and has near-perfect end-capping. In addition, it has stronger acid resistance than the COSMOSIL 5C<sub>18</sub>-AR. The COSMOSIL 5C<sub>18</sub>-AR-II column is especially effective for the separation of chelating compounds as well as both acidic and basic compounds.

## Acid resistance

The acidic resistance of COSMOSIL 5C<sub>18</sub>-AR-II is much improved compared with commercially available monomeric type octadecyl stationary phases. It retains high performance even in case of acidic mobile phases commonly used to separate acidic compounds and peptides.




Figure.  
Decomposition test in 0.1% Trifluoroacetic acid solution at 60°C.  
Capacity factor(k') = Naphthalene,  
Mobile phase: 70% Methanol


## Application data

### Peptides

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
Column size: 4.6mm I.D.-150mm  
Mobile phase: A; 0.05%TFA-H<sub>2</sub>O  
B; 0.05%TFA-Acetonitrile  
B conc. 10→40% 20min Linear gradient  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220nm

Sample: 1; Oxytocin (0.8μg)  
2; Angiotensin II(Human) (0.8μg)  
3; Angiotensin I(Human) (0.8μg)  
4; Substance P (0.8μg)




NACALAI TESQUE, INC  
AP-0351

### Organic acids

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
Column size: 4.6mm I.D.-150mm  
Mobile phase: Methanol/ 20mmol/l Phosphoric Acid = 20/80  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254nm, 0.5AUFS

Sample: 1; Gallic Acid (0.63μg)  
2; Protocatechuic Acid (0.63μg)  
3; Gentisic Acid (0.63μg)  
4; Phthalic Acid (0.63μg)




NACALAI TESQUE, INC  
AP-0159

### Salicylic acid esters

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
Column size: 4.6mm I.D.-150mm  
Mobile phase: Acetonitrile/ H<sub>2</sub>O = 70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254nm, 0.2AUFS

Sample: 1; Methyl Salicylate (2.3μg)  
2; Ethyl Salicylate (2.6μg)  
3; Propyl Salicylate (2.3μg)  
4; Butyl Salicylate (2.6μg)




NACALAI TESQUE, INC  
AP-0165

### Parabens

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
Column size: 4.6mm I.D.-150mm  
Mobile phase: Acetonitrile/ H<sub>2</sub>O = 50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254nm, 0.12AUFS


Sample: 1; Methyl *p*-Hydroxybenzoate (0.125μg)  
2; Ethyl *p*-Hydroxybenzoate (0.125μg)  
3; Propyl *p*-Hydroxybenzoate (0.125μg)  
4; Butyl *p*-Hydroxybenzoate (0.125μg)



NACALAI TESQUE, INC  
AP-0099

## LC/MS Application data

### Identification of herbal medicine constituents by LC/MS



## Ordering information

### Analytical / Preparative column (Particle size: 5 μm)

#### COSMOSIL 5C<sub>18</sub>-AR-II Packed Column

Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
1.0× 50	02955-21	4.6×100	38143-41
1.0×150	02951-61	4.6×150	38144-31
2.0× 30	05098-71	4.6×150 3 lots set*	09396-83
2.0× 50	34400-81	4.6×250	38145-21
2.0×100	34469-11	6.0×150	38146-11
2.0×150	37992-51	6.0×250	38147-01
2.0×250	05272-71	10× 50	05369-21
3.0×100	05791-71	10×150	34350-41
3.0×150	38028-61	10×250	38149-81
3.0×250	38029-51	20×150	34316-01
4.6× 30	05877-61	20×250	38150-41
4.6× 50	38142-51	28×250	34362-91

#### COSMOSIL 5C<sub>18</sub>-AR-II Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	38141-61
4.6×10 Cartridge**	38008-89
10×20	38148-91
20×20	34458-51
20×50	34479-81
28×50	34363-81

\* For 4.6×150 3 lots set, please refer to page 11.  
\*\* 3 cartridges included, needs a holder.  
Please refer to page 76.

### Preparative column (Particle size: 15 μm)

#### COSMOSIL 15C<sub>18</sub>-AR-II Packed Column

Column size I.D. x length (mm)	Product number
28×250	37978-51
50×250	38058-71
50×500	05884-61

#### COSMOSIL 15C<sub>18</sub>-AR-II Guard Column

Column size I.D. x length (mm)	Product number
28×50	38030-11
50×50	38057-81

### Fast LC column (Particle size: 3 μm)

#### COSMOSIL 3C<sub>18</sub>-AR-II Packed Column

Column size I.D. x length (mm)	Product number
2.0× 50	05478-91
4.6× 10	38068-41
4.6× 50	38069-31
4.6×100	38070-91


For more information on 15C<sub>18</sub>-AR-II, please refer to page 33.

For flow rate and device of semi-micro columns, or preparative columns, please refer to page 175.

# C<sub>18</sub>-PAQ

The COSMOSIL 5C<sub>18</sub>-PAQ maintains stable retention time even in 100% aqueous mobile phases. The new polymeric linking style gives this column a strong acidic resistance so that it is compatible with mobile phases of acidic pH that can permanently damage conventional octadecyl stationary phases. This type is especially good for separation of hydrophilic compounds.


## Stable performance under 100% aqueous conditions




The figure shows the change of retention time for thymine with 100% aqueous mobile phase (20 mmol/l phosphate buffer pH7). The sample was analyzed 4 times (1 hour). Flow of mobile phase was then stopped for 1 hour. The sample was analyzed under the same condition again after 1 hour. The conventional C<sub>18</sub> column showed change of retention time, but COSMOSIL 5C<sub>18</sub>-PAQ maintained stable retention time.

## Application data


### • Organic acids




### • Nucleobases and Nucleosides




### • dNTPs




### • Catecholamines



### • Ascorbic acid



### • 2-Phenylethyl glycosides



## Ordering information

### • Analytical / Preparative column (Particle size: 5 µm)

#### COSMOSIL 5C<sub>18</sub>-PAQ Packed Column

Column size I.D. x length (mm)	Product number
1.0× 50	05792-61
1.0×150	05793-51
2.0× 30	05878-51
2.0× 50	05794-41
2.0×100	05470-71
2.0×150	34449-71
2.0×250	05795-31
3.0×100	05796-21
3.0×150	05797-11
3.0×250	05798-01
4.6× 30	05879-41
4.6× 50	34451-21

#### COSMOSIL 5C<sub>18</sub>-PAQ Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	02484-91
10×20	34457-61
20×20	05803-11
20×50	05804-01
28×50	34455-81

### • Preparative column (Particle size: 15 µm)

#### COSMOSIL 15C<sub>18</sub>-PAQ Packed Column

Column size I.D. x length (mm)	Product number
28×250	05888-21
50×250	05890-71
50×500	05891-61

#### COSMOSIL 15C<sub>18</sub>-PAQ Guard Column

Column size I.D. x length (mm)	Product number
28×50	05887-31
50×50	05889-11






For more information on 15C<sub>18</sub>-PAQ, please refer to page 33.  
 For flow rate and device of semi-micro columns, or preparative columns, please refer to page 175.

## 2) Special columns

### Introduction


Reversed phase HPLC columns have been widely used because of their superior resolution, high theoretical plate number and ease of use. Since hydrophobic interaction is the dominant separation mechanism in reversed phase chromatography, conventional stationary phases such as C<sub>18</sub> and C<sub>8</sub> do not offer optimum selectivity for compounds with similar hydrophobicity. COSMOSIL offers a broad selection of columns with unique stationary phases for separation of these difficult analytes. These columns offer improved separation of structurally similar compounds that are difficult to analyze with a C<sub>18</sub> type column.

### Material characteristics

Packing material	Cholester	$\pi$ NAP	PYE	NPE	PBB-R
Silica gel	High purity porous spherical silica				
Average particle size	5 $\mu$ m*	5 $\mu$ m			
Average pore size	approx. 120 Å				
Specific surface area	approx. 300 m <sup>2</sup> /g				
Stationary phase					
	Cholesteryl group	Naphtylethyl group	Pyrenylethyl group	Nitrophenylethyl group	Pentabromobenzyl group
Bonding type	Monomeric type				
Main interaction	-Hydrophobic interaction -Molecular shape selectivity	-Hydrophobic interaction - $\pi$ - $\pi$ interaction	-Hydrophobic interaction - $\pi$ - $\pi$ interaction -Dispersion force -Charge-transfer interaction	-Hydrophobic interaction - $\pi$ - $\pi$ interaction -Dipole-dipole interaction	-Hydrophobic interaction -Dispersion force
End-capping treatment	Near-perfect treatment				
Carbon content	approx. 20%	approx. 11%	approx. 18%	approx. 9%	approx. 8%
Feature	-Usable under condition the same as C <sub>18</sub> -High molecular sharp selectivity	-Stronger $\pi$ - $\pi$ interaction than phenyl column	-Strongest $\pi$ - $\pi$ interaction	-Dipole-dipole interaction	-Dispersion force interaction

\* For 2.5  $\mu$ m particle size, please refer to page 34.

### Column selection guide




## Cholester

COSMOSIL Cholester is a reversed phase HPLC column with cholesteryl groups bonded silica packing material, which provides equivalent hydrophobicity like traditional alkyl group bonded silica packing materials (C<sub>18</sub>, C<sub>30</sub>). However, Cholester offers strong molecular shape selectivity for hydrophobic compounds to yield unique and reproducible separation patterns following the same analytical conditions used with other ODS columns.


### Hydrophobic interaction

Figure shows the comparison of hydrophobic interactions with competitive C<sub>18</sub> columns. Cholester provides the same hydrophobicity as alkyl group bonded types (C<sub>18</sub>, C<sub>30</sub>). It is not necessary to change the analytical conditions when replacing C<sub>18</sub> or C<sub>30</sub> columns with Cholester.



### Molecular shape selectivity



The stationary phase of Cholester has very rigid structures and can distinguish different molecular shapes. Cholester offers improved separation for structurally similar compounds that are difficult to analyze with alkyl group bonded materials (C<sub>18</sub> and C<sub>30</sub>). As in the following example Cholester retains planar Triphenylene longer than stereoscopic o-Terphenyl.



### Comparison of molecular shape selectivity


Column: 5C<sub>18</sub>-MS-II      Cholester

Column size: 4.6mm I.D.-150mm  
Mobile phase: Methanol/H<sub>2</sub>O = 90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254nm  
Sample: 1: o-Terphenyl (0.1  $\mu$ g)  
2: Triphenylene (0.01  $\mu$ g)




## Improvement in separation

COSMOSIL Cholester provides enhanced selectivity over traditional C<sub>18</sub> columns and offers greater performance in separating isomers or other closely related compounds. COSMOSIL Cholester is ideal for method development and serves as an excellent alternative to traditional C<sub>18</sub> columns. The figure below shows analytical data of chloroacetophenone isomers. These isomers are difficult to separate with C<sub>18</sub> and C<sub>30</sub>, but they are well resolved by COSMOSIL Cholester.



## Efficiency of preparative separation

The figure below shows the comparison of efficiency of preparative separation with a C<sub>18</sub> column. Both columns show good separation. However, sample loading capacity for preparative separations can be affected by a slight difference in separation ability. COSMOSIL Cholester can load 4 times of sample volume compared with C<sub>18</sub> columns.




## Application data

### • Catechins

#### COSMOSIL Application Data


Column: 4.6mm I.D. x 150mm  
 Mobile phase: A: Acetonitrile/ 20mmol/l Phosphate buffer(pH2.5) = 10/90  
 B: Acetonitrile/ 20mmol/l Phosphate buffer(pH2.5) = 30/70  
 B conc. 0→100% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm  
 Sample: 1; Gallic acid(GC) (0.80mg/ml)  
 2; Caffeine (0.08mg/ml)  
 3; Epigallocatechin(EGC) (0.80mg/ml)  
 4; Catechin (C) (0.40mg/ml)  
 5; Epicatechin(EC) (0.40mg/ml)  
 6; Epigallocatechin gallate(EGCG) (0.20mg/ml)  
 7; Gallic acid gallate(GCG) (0.40mg/ml)  
 8; Epicatechin gallate(EGC) (0.20mg/ml)  
 9; Catechin gallate(CG) (0.20mg/ml)  
 Injection Vol. 1.0μl



### • Flavones

#### COSMOSIL Application Data

Column: 4.6mm I.D. x 150mm  
 Mobile phase: A: Acetonitrile/ 20mmol/l Phosphate buffer(pH2.5) = 20/80  
 B: Acetonitrile/ 20mmol/l Phosphate buffer(pH2.5) = 70/30  
 B conc. 0→100% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm  
 Sample: 1; Fisetin (0.25μg)  
 2; Myricetin (0.20μg)  
 3; 7,8-Dihydroxyflavone (0.05μg)  
 4; Luteolin (0.20μg)  
 5; Quercetin (0.20μg)  
 6; 7-Hydroxyflavone (0.10μg)      10; Chrysin (0.05μg)  
 7; Baicalein (0.05μg)      11; 6-Methoxyflavone (0.05μg)  
 8; 6-Hydroxyflavone (0.05μg)      12; 3-Hydroxyflavone (0.25μg)  
 9; Flavone (0.05μg)      13; 5-Hydroxyflavone (0.05μg)



## Ordering information

### • Analytical / Preparative column (Particle size: 5 μm)

#### COSMOSIL Cholester Packed Column

Column size I.D. x length (mm)	Product number
1.0×150	05968-71
1.0×250	05969-61
2.0× 30	08565-51
2.0× 50	06352-91
2.0×100	06948-01
2.0×150	05971-11
2.0×250	05972-01
3.0×150	05973-91
3.0×250	05974-81

Column size I.D. x length (mm)	Product number
4.6×150	05976-61
4.6×150 3 lots set※	07970-03
4.6×250	05977-51
10×150	08011-91
10×250	05979-31
20×150	06088-71
20×250	05982-71
28×250	05985-41

※ For 4.6×150 3 lots set, please refer to page 11.


※ For 2.5 μm particle size, please refer to page 34

#### COSMOSIL Cholester Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	05975-71
10×20	05978-41
20×20	05980-91
20×50	05981-81
28×50	05983-61

COSMOSIL πNAP is a reversed phase HPLC column with naphthylethyl group bonded silica packing material. The naphthylethyl group is composed of two fused aromatic rings and forms strong π-π interactions with unsaturated compounds. This column offers improved separation of compounds such as positional isomers that are difficult to analyze with alkyl group bonded materials.

Comparison of π-π interactions




COSMOSIL πNAP shows stronger π-π interactions than phenyl columns.


Figure. Comparison of π-π interaction

Application data


Tolunitriles




Adrenal Cortical Hormones



Berberines



Isoflavones



Ordering information

Analytical / Preparative column (Particle size: 5 μm)

COSMOSIL πNAP Packed Column		COSMOSIL πNAP Guard Column	
Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
1.0×150	08076-61	3.0×250	08081-81
1.0×250	08077-51	4.6×150	08085-41
2.0× 30	08566-41	4.6×250	08086-31
2.0× 50	08567-31	10×150	08088-11
2.0×100	08299-51	10×250	08089-01
2.0×150	08078-41	20×150	08092-41
2.0×250	08079-31	20×250	08093-31
3.0×150	08080-91	28×250	08095-11

COSMOSIL PYE column is a reversed phase column with 2-(1-pyrenyl) ethyl groups bonded silica packing material. This column utilizes  $\pi$ - $\pi$  interactions originating from the planar pyrene ring structure to separate structural isomers.

Comparison  $\pi$ - $\pi$  interaction






Figure. Comparison of  $\pi$ - $\pi$  interactions

Application data


• Sterols




• Phthalonitriles



• Tocopherols



• Capsaicins



Attention

1. Methanol is recommended as a mobile phase for COSMOSIL PYE column. Acetonitrile is not recommended because it has many  $\pi$  electrons and interferes  $\pi$ - $\pi$  interactions between a sample and the stationary phase.
2. The stationary phase of COSMOSIL PYE, nitrophenyl group, has a large UV absorption. When the stationary phase detaches from silica gel and elutes, even a slight quantity can be detected and causes baseline noise. In such a case, wash the column with tetrahydrofuran. Detachment of a small amount of the stationary phase does not deteriorate a column's separation ability.
3. COSMOSIL PYE column is not suitable for gradient analysis.

Ordering information

• Analytical / Preparative column (Particle size: 5  $\mu$ m)

COSMOSIL 5PYE Packed Column

Column size I.D. x length (mm)	Product number
1.0×150	02851-71
2.0×150	38042-61
2.0×250	34450-31

COSMOSIL 5PYE Guard Column


Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
4.6×150	37837-91	4.6×10	37903-11
4.6×250	37989-11	10×20	38041-71
10×250	37996-11	20×20	05867-91
20×250	38044-41	20×50	34475-21

# NPE

COSMOSIL NPE column is a reversed phase column with nitrophenylethyl groups bonded silica packing material. This column provides unique retention characteristics, slightly different from the COSMOSIL PYE column, utilizing both dipole-dipole and  $\pi$ - $\pi$  interactions.

## Selectivity for dipole-dipole interactions

COSMOSIL NPE strongly retains 1,8-dinitronaphthalene because of the strong dipole formed by the two nitro groups positioned on the same side of naphthalene.



## Attention

1. Methanol is recommended as a mobile phase for COSMOSIL NPE column. Acetonitrile is not recommended because it has many  $\pi$  electrons and interferes  $\pi$ - $\pi$  interactions between a sample and the stationary phase.
2. The stationary phase of COSMOSIL NPE, nitrophenyl group, has a large UV absorption. When the stationary phase detaches from silica gel and elutes, even a slight quantity can be detected and causes baseline noise. In such a case, wash the column with tetrahydrofuran. Detachment of a small amount of the stationary phase does not deteriorate a column's separation ability.
3. COSMOSIL NPE column is not suitable for gradient analysis.

## Ordering information

- Analytical / Preparative column (Particle size: 5  $\mu$ m)

### COSMOSIL 5NPE Packed Column

Column size I.D. x length (mm)	Product number
1.0×150	05897-01
2.0×150	34328-51
2.0×250	34379-91

### COSMOSIL 5NPE Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	37902-21
4.6×250	37990-71
10×20	38045-31
20×20	05868-81
20×50	05869-71


# PBB-R

COSMOSIL PBB-R is a reversed phase column with pentabromobenzyl groups bonded silica packing material. This column provides unique selectivity for structurally similar compounds utilizing the dispersion force interaction. The dispersion force interaction of COSMOSIL PBB-R makes it useful for separation of structural isomers differing only by a double bond.

## Application data


- Halogen exchange reaction products

COSMOSIL PBB-R strongly retains iodine atom which has a large dispersion force, than bromine atom. So it can separate halogen exchange reaction products that are difficult to analyze with C<sub>18</sub> column.



- Surfactant agents

C<sub>18</sub> column can not separate Triton X-100 mixture, because (-OCH<sub>2</sub>CH<sub>2</sub>-) group has little hydrophobicity. However, COSMOSIL PBB-R can separate them because it distinguishes difference in the dispersion force, which depends on its molecular weight.



## Ordering information

- Analytical / Preparative column (Particle size: 5  $\mu$ m)

### COSMOSIL 5PBB-R Packed Column

Column size I.D. x length (mm)	Product number
1.0×150	05899-81
2.0×150	05900-31
2.0×250	05904-91

### COSMOSIL 5PBB-R Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	05697-21
4.6×250	05698-11
10×20	05699-01
20×20	05700-51
4.6×10	05704-11
10×20	05721-81
20×20	05911-91
20×50	05722-71




### 3) Other columns

#### Introduction


COSMOSIL alkyl type columns also include stationary phases with packing materials of C<sub>22</sub> (Dococyl group), C<sub>8</sub> (Octyl group), C<sub>4</sub> (Butyl group) and TMS (Trimethyl group). The order of retention force by hydrophobicity of each packing material is C<sub>22</sub>>C<sub>8</sub>>C<sub>4</sub>>TMS. The columns having lower hydrophobicity than C<sub>18</sub> or C<sub>22</sub> are effective for separation of high hydrophobic compounds and compounds with big difference in hydrophobic. Hydrophobicity of C<sub>22</sub> is about the same as C<sub>18</sub>. However, stereoselectivity of C<sub>22</sub> is higher than C<sub>18</sub>, and so C<sub>22</sub> may provide better separation.

The COSMOSIL PE-MS columns (phenylethyl group) and the COSMOSIL CN-MS columns (cyanopropyl group) provide a secondary separation mode ( $\pi$ - $\pi$  interaction). These columns are recommended when the other alkyl chain stationary phases do not offer optimum selectivity for structurally similar compounds.

#### Material characteristics

Packing material	C <sub>22</sub> -AR-II	C <sub>8</sub> -MS	C <sub>4</sub> -MS	TMS-MS	PE-MS	CN-MS
Silica gel	High purity porous spherical silica					
Average particle size	5 $\mu$ m					
Average pore size	approx. 120 Å					
Specific surface area	approx. 300 m <sup>2</sup> /g					
Stationary phase						
	Dococyl group	Octyl group	Butyl group	Trimethyl group	Phenylethyl group	Cyanopropyl group
Bonding type	Polymeric	Monomeric				
Main interaction	Hydrophobic interaction				Hydrophobic interaction $\pi$ - $\pi$ interaction	
End-capping treatment	Near-perfect treatment					
Carbon content	approx. 19%	approx. 10%	approx. 7%	approx. 5%	approx. 10%	approx. 7%


#### Different of separation characteristic



## CN-MS

Gradient elution is commonly used for the samples containing both polar and non-polar compounds. However, gradient elution may cause reproducibility problem depending on the gradient mixer and pump, and need an equilibration time for each analysis. COSMOSIL 5CN-MS offers rapid analysis and great reproducibility using isocratic elution mode.

#### Rapid analysis of samples containing both polar and non-polar compounds




#### Application data

##### • Acetoaminophen

#### COSMOSIL Application Data

Column: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 50mmol/l KH<sub>2</sub>PO<sub>4</sub>  
 (pH4.7 with NaOH) = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV225nm

Sample: 1; *p*-Aminophenol Hydrochloride (0.2 $\mu$ g)  
 2; Acetaminophen (0.2 $\mu$ g)  
 3; 4'-Acetoxyacetanilide (0.2 $\mu$ g)



#### Ordering information

##### • Analytical / Preparative column (Particle size: 5 $\mu$ m)

##### COSMOSIL 5CN-MS Packed Column

Column size I.D. x length (mm)	Product number
4.6× 50	38233-61
4.6×100	38234-51
4.6×150	38235-41
4.6×250	38236-31

Column size I.D. x length (mm)	Product number
6.0×150	38237-21
6.0×250	38238-11
10×250	38239-01
20×250	38240-61

##### COSMOSIL 5CN-MS Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	38231-81
10×20	38232-71

# C<sub>22</sub>-AR-II, C<sub>8</sub>-MS, C<sub>4</sub>-MS, TMS-MS, PE-MS

## Ordering information

- Analytical / Preparative column (Particle size: 5 μm)

### COSMOSIL 5C<sub>22</sub>-AR-II Packed Column

Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
4.6× 50	05848-41	6.0×150	05850-91
4.6×100	05849-31	6.0×250	05851-81
4.6×150	04598-51	10×250	04969-91
4.6×250	04599-41	20×250	05183-41

### COSMOSIL 5C<sub>22</sub>-AR-II Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	04881-21
10×20	05554-81

### COSMOSIL 5C<sub>8</sub>-MS Packed Column

Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
4.6× 50	38153-11	6.0×150	38157-71
4.6×100	38154-01	6.0×250	38158-61
4.6×150	38155-91	10×250	38159-51
4.6×250	38156-81	20×250	38160-11

### COSMOSIL 5C<sub>8</sub>-MS Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	38151-31
10×20	38152-21

### COSMOSIL 5C<sub>4</sub>-MS Packed Column

Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
4.6× 50	38163-81	6.0×150	38167-41
4.6×100	38164-71	6.0×250	38168-31
4.6×150	38165-61	10×250	38169-21
4.6×250	38166-51	20×250	38170-81

### COSMOSIL 5C<sub>4</sub>-MS Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	38161-01
10×20	38162-91

### COSMOSIL 5TMS-MS Packed Column

Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
4.6× 50	38173-51	6.0×150	38177-11
4.6×100	38174-41	6.0×250	38178-01
4.6×150	38175-31	10×250	38179-91
4.6×250	38176-21	20×250	38180-51

### COSMOSIL 5TMS-MS Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	38171-71
10×20	38172-61

### COSMOSIL 5PE-MS Packed Column

Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
4.6× 50	38183-21	6.0×150	38187-81
4.6×100	38184-11	6.0×250	38188-71
4.6×150	38185-01	10×250	38189-61
4.6×250	38186-91	20×250	38190-21

### COSMOSIL 5PE-MS Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	38181-41
10×20	38182-31

## 4) Silica based preparative columns

# 15C<sub>18</sub>-MS-II, 15C<sub>18</sub>-AR-II, 15C<sub>18</sub>-PAQ

COSMOSIL series is available in 10 mm I.D. and 20 mm I.D. for semi-preparative applications and 28 mm I.D. and 50 mm I.D. for preparative scales. For column sizes and packing materials not listed below, contact either your local distributor or the manufacturer directly.


## Material characteristics

Packing material	15C <sub>18</sub> -MS-II	15C <sub>18</sub> -AR-II	15C <sub>18</sub> -PAQ
Silica gel	High purity porous spherical silica		
Average particle size	15 μm		
Average pore size	approx. 120 Å		
Specific surface area	approx. 300 m <sup>2</sup> /g		
Stationary phase	Octadecyl group (please refer to page 12.)		
Bonding type	Monomeric	Polymeric	
Main interaction	Hydrophobic interaction		
End-capping treatment	Near-perfect treatment		
pH range	2 ~ 10	1.5 ~ 7.5	2 ~ 7.5
Carbon content	approx. 16%	approx. 17%	approx. 11%
Feature	This phase is recommended for most of applications but particularly effective for basic organic compounds.	This phase is recommended for the separations requiring acidic mobile phase conditions. It also shows superior molecular shape selectivity to monomeric type C <sub>18</sub> columns.	This phase is designed to offer superior retention of polar compounds and excellent reproducibility in highly aqueous mobile phases, even in 100% aqueous.

## Application data

- Preparative separation using 50 mm I.D. column

Carbazole is extracted from anthracene oil (coal tar) and required high purity because it is often used for analytical applications. Following is the preparative separation of carbazole using a 50 mm I.D. COSMOSIL 15C<sub>18</sub>-AR-II.



Please refer to TECHNICAL NOTE 8, Inner diameter of column (scale down and scale up) at page 175.

## Ordering information

Please refer to page 15 for 15C<sub>18</sub>-MS-II, page 17 for 15C<sub>18</sub>-AR-II and page 19 for 15C<sub>18</sub>-PAQ.

## 5) Ultra-Fast Liquid Chromatography (UFLC)

# 2.5C<sub>18</sub>-MS-II, 2.5Cholester

Ultra-Fast Liquid Chromatography (UFLC) columns filled with sub-2 μm particles have become widely used. Smaller particle columns generate higher pressures, while provide higher theoretical plate number. Therefore the column length should be shorter, and lead to lower resolution. COSMOSIL 2.5Cholester columns provide enhanced selectivity over traditional C<sub>18</sub> materials, and greater performance in separating isomers or other closely related compounds.


### Material characteristics

Packing material	C <sub>18</sub> -MS-II	Cholester
Silica gel	High purity porous spherical silica	
Average particle size	2.5 μm	
Average pore size	approx. 130 Å	
Specific surface area	approx. 330 m <sup>2</sup> /g	
Stationary phase	Octadecyl group	Cholesteryl group
Bonding type	Monomeric	
Main interaction	Hydrophobic interaction	Hydrophobic interaction Molecular shape selectivity
End-capping treatment	Near-perfect treatment	
Carbon content	approx. 18%	approx. 21%
Features	◦First choice of reversed phase column	◦The same mobile phase as C <sub>18</sub> ◦High molecular shape selectivity

### Ultra-Fast Liquid Chromatography


Very fast and efficient separation can be achieved using 2.5 μm particles.

Note: Ultra high pressure liquid chromatography system or some modification of HPLC system is required for UFLC analysis.



### Comparison of analytical pressure


COSMOSIL 2.5Cholester and 2.5MS-II are reversed phase columns with 2.5 μm particle size of silica packing material. It can be used under lower pressure than competitor's 2 μm columns.




### Application data

#### • Saikosaponins

COSMOSIL 2.5Cholester offers improved resolution for compounds difficult to analyze with C<sub>18</sub> without changing analytical condition.



#### • Catechins



### Ordering information

#### • Analytical column (Particle size: 2.5 μm)

##### COSMOSIL 2.5C<sub>18</sub>-MS-II Packed Column

Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
2.0× 50	08994-31	3.0× 50	08997-01
2.0× 75	08995-21	3.0× 75	08998-91
2.0×100	08996-11	3.0×100	08999-81

##### COSMOSIL 2.5Cholester Packed Column

Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
2.0× 50	09000-01	3.0× 50	09049-91
2.0× 75	09047-11	3.0× 75	09050-51
2.0×100	09048-01	3.0×100	09051-41

# 6. Normal phase chromatography


## SL-II

Ultra-pure silica gel of more than 99.99% purity is used for the COSMOSIL SL-II packed column series. This column provides improved separation and reproducibility for compounds with carbonyl or phenol hydroxyl groups, which are often problematic to separate using conventional silica gel columns.

### Material characteristics


Packing material	SL-II
Silica gel	High purity porous spherical silica
Average particle size	3.5, 15 μm
Average pore size	approx. 120 Å
Specific surface area	approx. 300 m <sup>2</sup> /g
Feature	•High purity silica gel (>99.99%) with specially treatment •Suitable for preparative separation

### Comparison with old type




### Comparison with medium-pressure column

COSMOSIL SL-II offers sharper peak compared with packing materials for medium-pressure liquid chromatography and open chromatography.



### Scaling up from analytical to preparative separation



Please refer to TECHNICAL NOTE 8, Inner diameter of column (scale down and scale up) at page 175

### Ordering information

- Analytical / Preparative column (Particle size: 5 μm)

#### COSMOSIL 5SL-II Packed Column

Column size I.D. x length (mm)	Product number
4.6× 50	37999-81
4.6×100	38000-01
4.6×150	38001-91
4.6×250	38002-81

#### COSMOSIL 5SL-II Guard Column

Column size I.D. x length (mm)	Product number
6.0×150	38003-71
6.0×250	38004-61
10×250	38005-51
20×250	38006-41
28×250	34358-61

Column size I.D. x length (mm)	Product number
4.6×10	37997-01
10×20	37998-91
20×20	05874-91
20×50	05875-81
28×50	34359-51

- Preparative column (Particle size : 15 μm)

#### COSMOSIL 15SL-II Packed Column

Column size I.D. x length (mm)	Product number
28×250	05893-41
50×250	05895-21
50×500	05896-11

#### COSMOSIL 15SL-II Guard Column

Column size I.D. x length (mm)	Product number
28×50	05892-51
50×50	05894-31

- Columns of high number of theoretical plates (Particle size : 3 μm)

#### COSMOSIL 3SL-II Packed Column


Column size I.D. x length (mm)	Product number
4.6× 10	38059-61
4.6× 50	38060-21
4.6×100	38061-11

# 7. Hydrophilic interaction chromatography

## HILIC


COSMOSIL HILIC is a new column for hydrophilic interaction chromatography with Triazole bonded silica packing material. The hydrophilic interaction chromatography is a variation of normal phase chromatography where a polar stationary phase is used with a mobile phase which contains a high concentration of water miscible organic solvent and a low concentration of aqueous eluent. The main retention mechanism is the partitioning of the polar analytes between the polar stationary and the non-polar mobile phase. As it is also called "aqueous normal phase", the elution order is similar to that of normal phase and the sample elution is in the order of increasing hydrophilicity. Without using ion-pair reagent COSMOSIL HILIC retains highly polar analytes that would not be retained in reversed phase chromatography. It also shows a weak anion-exchange mechanism with the positively charged stationary phase, thus acidic compound is strongly retained.

### Material characteristics

Packing material	HILIC
Silica gel	High purity porous spherical silica
Average particle size	5 μm
Average pore size	approx. 120 Å
Specific surface area	approx. 300 m <sup>2</sup> /g
Stationary phase	 Triazole
Interactions	Hydrophilic interaction, anion exchange
Object substance	Hydrophilic compounds, acidic compounds
Feature	Suitable for non-retaining by C <sub>18</sub>


### Comparison with C<sub>18</sub>

COSMOSIL HILIC can separate glycine and glycylglycine without ion-pair reagent. Although C<sub>18</sub> column can separate them with ion-pair reagents, there are some disadvantages such as column equilibration, preparation of mobile phase and column deterioration.




### Application data

#### Melamine related compounds





#### BCAA (amino acid branched-chain)



### Optimizing analytical conditions

COSMOSIL HILIC Chromatogram Index, which includes 154 chromatograms using COSMOSIL HILIC, is now available online at <http://www.nacalai.co.jp/en/cosmosil/>. This index is useful for optimizing analytical conditions for hydrophilic interaction chromatography.



### Ordering information

#### Analytical / Preparative column (Particle size: 5 μm)

##### COSMOSIL HILIC Packed Column

Column size I.D. x length (mm)	Product number
1.0×150	07869-11
1.0×250	07870-71
2.0× 30	08568-21
2.0× 50	07052-91
2.0×100	08569-11
2.0×150	07054-71
2.0×250	07489-91
3.0×150	07871-61
3.0×250	07872-51

Column size I.D. x length (mm)	Product number
4.6×150	07056-51
4.6×150 3 lots set*	09385-23
4.6×250	07057-41
10×250	07059-21
20×250	07060-81
28×250	07875-21

\* for 4.6×150 3 lots set, please refer to page 11.

##### COSMOSIL HILIC Guard Column


Column size I.D. x length (mm)	Product number
4.6×10	07055-61
10×20	07058-31
20×20	07854-91
20×50	07873-41
28×50	07874-31

# 8. Saccharide analysis

## Introduction


Saccharides are not retained on standard C<sub>18</sub> columns because of the low hydrophobicity of compounds. COSMOSIL Sugar-D and NH<sub>2</sub>-MS are specifically designed for separation of saccharides. COSMOSIL C<sub>18</sub>-PAQ is recommended for hydrophobic glycosides or saccharide derivatives.

## Material characteristics

Packing material	Sugar-D	NH <sub>2</sub> -MS
Silica gel	High purity spherical silica	
Average particle size	5 μm	
Average pore size	-	approx. 120 Å
Specific surface area	-	approx. 300 m <sup>2</sup> /g
Stationary phase	-	
Bonding type	Secondary/Tertiary amine	Aminopropyl
Object substance	Monosaccharides, oligosaccharides	
End capping treatment	-	Near-perfect treatment
Carbon content	-	approx. 4 %
Feature	<ul style="list-style-type: none"> <li>First choice of saccharide analysis</li> <li>High durability</li> <li>Good quantitative analysis</li> </ul>	<ul style="list-style-type: none"> <li>Different selectivity from Sugar-D</li> </ul>

## Comparison of retention


The conventional aminopropyl column is slightly more retentive than Sugar-D. The retention time of Sugar-D can be adjusted by increasing the concentration of acetonitrile in the mobile phase by 5%-10%.



# Sugar-D


Conventionally aminopropyl bonded stationary phases are used for liquid chromatographic analysis of mono- and oligosaccharides. General shortcomings of the conventional aminopropyl bonded phases are tailing and adsorption of certain saccharides and the general low durability (short active life) of these columns. These problems are addressed and solved by the novel COSMOSIL Sugar-D, resulting in better (sharper) separation and much improved durability.

## Durability



## Adsorption characteristics

Certain types of saccharides such as arabinose or galactose are partially or temporarily adsorbed on conventional aminopropyl stationary phases causing tailing or no elution at all. COSMOSIL Sugar-D provides superior separation and high recovery for these saccharides.




## Application data

### • Mono- and Oligosaccharides

#### COSMOSIL Application Data

Column: Sugar-D  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile / H<sub>2</sub>O = 75/25  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample: 1; Rhamnose (10μg)  
 2; Fructose (10μg)  
 3; Glucose (10μg)  
 4; Sucrose (10μg)  
 5; Maltose (10μg)  
 6; Raffinose (10μg)




### • Polyols

#### COSMOSIL Application Data

Column: Sugar-D  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile / H<sub>2</sub>O = 75/25  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample: 1; Glycerol (10μg)  
 2; meso-Erythritol (10μg)  
 3; Xylitol (10μg)  
 4; Glucitol (10μg)  
 5; Maltitol (10μg)  
 6; Inositol (10μg)





# 9. Protein separation wide pore columns

## 1) Reversed phase chromatography

### Protein-R


COSMOSIL Protein-R is a reversed phased HPLC column designed specifically for peptide and protein separations.

#### Material characteristics


Packing material	Protein-R
Silica gel	High purity spherical porous silica
Average particle size	5 μm
Average pore size	approx. 300 Å
Specific surface area	approx. 150 m <sup>2</sup> /g
Stationary phase	Octadecyl group
Bonding type	Polymeric
Main Interactions	Hydrophobic interaction
End capping treatment	Near-perfect treatment
Feature	•High recovery rate •Acid-resistant

#### Comparison of separation

Protein-R shows sharper peaks for proteins than conventional C<sub>18</sub> wide pore columns.




#### Recovery rate



The figure left shows recovery rates for proteins using different columns. Protein-R shows a higher recovery rate than C<sub>4</sub>-300 and a much higher recovery rate than C<sub>18</sub>-300.

#### Durability


The figure below shows durability against acidic mobile phase of various columns. Protein-R shows a higher acid durability than C<sub>4</sub>-300.



Degradation test with 0.1%-Trifluoroacetic Acid at 60°C  
 (K): Napthalene in the mobile phase (Methanol : Water = 50 : 50 )

#### Application of peptide separation

5C<sub>18</sub>-MS-II (pore size 120 Å) shows better separation of low molecular weight proteins, but Protein-R shows better separation of high-molecular weight proteins.



#### Ordering information

• Analytical / Preparative column (Particle size: 5 μm)

COSMOSIL Protein-R Packed Column

Column size I.D. x length (mm)	Product number
2.0×150	06514-71
4.6× 50	06525-31
4.6×150	06526-21
4.6×250	06527-11





COSMOSIL Protein-R Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	06518-31
10×20	06528-01
20×20	08692-81


# C<sub>18</sub>-AR-300, C<sub>8</sub>-AR-300, C<sub>4</sub>-AR-300, Ph-AR-300

COSMOSIL offers a variety of stationary phases with wide-pore silica gel material for separations of polypeptides and proteins.

## Material characteristics

Packing material	5C <sub>18</sub> -AR-300	5C <sub>8</sub> -AR-300	5C <sub>4</sub> -AR-300	5Ph-AR-300
Silica gel	High purity spherical porous silica			
Average particle size	5 μm			
Average pore size	approx. 300 Å			
Specific surface area	approx. 150 m <sup>2</sup> /g			
Stationary phase				
	Octadecyl group	Octyl group	Butyl group	Phenyl group
Bonding type	Polymeric			
Main Interaction	Hydrophobic interaction			Hydrophobic interaction π-π interaction
End capping treatment	Near-perfect treatment			
Carbon content	approx. 12%	approx. 7%	approx. 6%	approx. 7%

## Comparison of separation



## Ordering information

• Analytical / Preparative column (Particle size: 5 μm)

### COSMOSIL 5C<sub>18</sub>-AR-300 Packed Column

Column size I.D. x length (mm)	Product number
4.6× 50	37911-01
4.6×150	37913-81
4.6×250	37914-71

Column size I.D. x length (mm)	Product number
10×150	37917-41
10×250	37918-31
20×150	37919-21
20×250	37920-81

### COSMOSIL 5C<sub>18</sub>-AR-300 Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	37910-11
10×20	37965-11

### COSMOSIL 5C<sub>8</sub>-AR-300 Packed Column

Column size I.D. x length (mm)	Product number
4.6× 50	37951-81
4.6×150	37953-61
4.6×250	37954-51

Column size I.D. x length (mm)	Product number
10×150	34345-21
10×250	34247-11
20×150	05861-51
20×250	34364-71

### COSMOSIL 5C<sub>8</sub>-AR-300 Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	37950-91
10×20	34464-61

### COSMOSIL 5C<sub>4</sub>-AR-300 Packed Column

Column size I.D. x length (mm)	Product number
4.6× 50	37956-31
4.6×150	37958-11
4.6×250	37959-01

Column size I.D. x length (mm)	Product number
10×150	34249-91
10×250	38047-11
20×150	34477-01
20×250	38048-01

### COSMOSIL 5C<sub>4</sub>-AR-300 Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	37955-41
10×20	05862-41

### COSMOSIL 5Ph-AR-300 Packed Column

Column size I.D. x length (mm)	Product number
4.6× 50	37961-51
4.6×150	37963-31
4.6×250	37964-21

Column size I.D. x length (mm)	Product number
10×150	05865-11
10×250	34267-51
20×150	05866-01
20×250	34468-21

### COSMOSIL 5Ph-AR-300 Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	37960-61
10×20	34268-41

## 2) Gel filtration chromatography


# Diol-120-II, Diol-300-II

COSMOSIL Diol-120-II and Diol-300-II gel filtration columns are ideal for the size-based separation of proteins and other water soluble polymers. The separation MW range is 5,000 - 700,000 daltons for proteins and 300 - 300,000 daltons for water soluble polymers when Diol-120-II and Diol-300-II are used in series.

### Material characteristics

Packing material	5Diol-120-II	5Diol-300-II
Silica gel	High purity spherical porous silica	
Average particle size	5 μm	
Average pore size	approx. 120 Å	approx. 300 Å
Specific surface area	Diol group	
Object substance	Proteins, water soluble polymers	
Flow rate	0.5-1.0 (ml/min)	
Selection of pore size (Protein)	5,000-100,000 Da	10,000-700,000 Da
Selection of pore size (Water soluble polyme)	300-30,000 Da	500-300,000 Da


### Calibration curve of proteins



Column: COSMOSIL 5Diol-II 7.5mm I.D. x 600mm  
 Mobile phase: 20mmol/l Phosphate Buffer (pH 7.0) + 100mmol/l Na<sub>2</sub>SO<sub>4</sub>  
 Flow rate: 1.0ml/min  
 Temperature: 30°C


Sample	M.W.	Sample	M.W.
Thyroglobulin	660,000	Peroxidase	40,000
Catalase	250,000	Carbonic Anhydrase	30,000
Glucose Oxidase	160,000	α-Chymotrypsinogen A	25,700
Uricase	128,000	α-Chymotrypsin	25,200
Choline Oxidase	95,000	Trypsinogen	24,000
Transferrin	85,000	Trypsin (bovine)	23,300
Conalbumin	77,500	Myoglobin	17,000
Malate Dehydrogenase	70,000	Lysozyme	14,300
α-Glucosidase	68,500	Ribonuclease A	13,700
Albumin (BSA)	66,000	Cytochrome C	12,400
α-Amylase	52,500	Aprotinin	6,500
Fetuin	48,000	Gly-Gly	132
Albumin (Ovalbumin)	45,000		

### Comparison of separation property




## Application data


### • E.Coli Proteins




### • Water Soluble Polystyrene



### • Polyethylene Glycol



### • Pullulan



## Ordering information

### • Analytical column (Particle size: 5 μm)

COSMOSIL 5Diol-120-II Packed Column		COSMOSIL 5Diol-120-II Guard Column	
Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
7.5×300	38050-51	7.5×50	38049-91
7.5×600	38051-41		

COSMOSIL 5Diol-300-II Packed Column		COSMOSIL 5Diol-300-II Guard Column	
Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
7.5×300	38053-21	7.5×50	38052-31
7.5×600	38054-11		

### 3) Ion-exchange chromatography

## DEAE, QA, CM, SP

The packing materials for COSMOGEL ion-exchange glass columns are based on hydrophilic polymethacrylate 10 µm particles with a 1000 Å pore size. COSMOGEL packed columns are available with DEAE, a weak anion exchanger; QA, a strong anion exchanger; CM, a weak cation exchanger; and SP, a strong cation exchanger. The availability of four different ion exchangers provides chromatographers with the flexibility of column selection based on charge differences of samples.

#### Material characteristics

Packing material	DEAE	QA	CM	SP
Type	Diethylaminoethyl type Weak anion exchange	Quaternary ammonium type Strong anion exchange	Carboxymethyl type Weak cation exchange	Sulfopropyl type Strong cation exchange
Gel	Totally porous spherical hydrophilic polymer			
Average particle size	10 µm			
Average pore size	approx. 1,000 Å			
Functional group	N <sup>+</sup> H(C <sub>2</sub> H <sub>5</sub> ) <sub>2</sub>	N <sup>+</sup> H(CH <sub>3</sub> ) <sub>3</sub>	COO <sup>-</sup>	SO <sub>3</sub> <sup>-</sup>
Counter ion	Cl <sup>-</sup>	Cl <sup>-</sup>	Na <sup>+</sup>	Na <sup>+</sup>
Capacity (meq/g)	0.6	0.4	0.3	0.4
50% Ionization pH	10.8	11.0	5.7	2.6
pH range	<11	Whole area	>4	Whole area
Flow rate (appropriate)	7.5 mm I.D. / 8.0 mm I.D. ; 0.5 ~ 1.0 ml/min		20 mm I.D. ; 4.0 ~ 6.0 ml/min	
Flow rate (maximum)	7.5 mm I.D. / 8.0 mm I.D. ; 1.5 ml/min		20 mm I.D. ; 8.0 ml/min	
Pressure (maximum)	1.5 MPa			
Temperature	10 ~ 50°C			

#### Collection rate of proteins

##### Anion exchange columns

Protein	Collection rate(%)	
	DEAE	QA
Catalase	93	92
BSA	98	97
Ovalbumin	100	103
β-Lactoglobulin	95	102
Myoglobin	97	101
Transferrin	105	102

Sample: 0.02 M Tris-HCl buffer (pH 8.2)  
Elution: Tris-HCl buffer (pH 8.2) + 0.5 M NaCl

##### Cation exchange columns

Pretein	Collection rate(%)	
	CM	SP
γ-Globulin	97	96
α-Chymotrypsinogen	96	103
Myoglobin	99	95
Cytochrome C	98	94

Sample: 0.02 M Tris-HCl buffer (pH 6.5)  
Elution: Tris-HCl buffer (pH 6.5) + 0.5 M NaCl


### Difference between anion exchange columns and cation exchange columns

- Separation of proteins on anion exchange columns

The higher the negative charge, the longer the sample is retained on an anion exchange column. As shown below, with a weak alkaline mobile phase, the lower the isoelectric point, the longer the sample is retained.


#### Separation on anion exchange columns

Column: 7.5mm I.D.-75mm  
 Column size: 7.5mm I.D.-75mm  
 Mobile phase: A: 20mmol/l Tris-HCl buffer (pH 8.2)  
 B: 0.5mol/l NaCl in buffer A  
 B conc. 0 → 100% (30min)  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm  
 Sample: Lysoxtyase (pI 5)



#### Separation on cation exchange columns

Column: 7.5mm I.D.-75mm  
 Column size: 7.5mm I.D.-75mm  
 Mobile phase: A: 20mmol/l Acetate buffer (pH 4.5)  
 B: 0.5mol/l Na<sub>2</sub>SO<sub>4</sub> in buffer A  
 B conc. 0 → 100% (30min)  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm  
 Sample: Lysoxtyase (pI 5)




### Sample loading capacity and resolution

Up to 1 mg of sample can be well separated on an 8 mm I.D. column.

If the resolution is high enough, 10 mg of sample can be separated.

#### Sample loading capacity and resolution

Column: DEAE  
 Column size: 8.0mm I.D.-75mm  
 Mobile phase: A: 20mmol/l Tris-HCl buffer (pH 8.2)  
 B: 0.5mol/l NaCl in buffer A  
 B conc. 0 → 100% (30min)  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm, 0.16AUFS  
 Sample: 1; Myoglobin (20µg)  
 2; Transferrin (50µg)  
 3; Trypsin Inhibitor (70µg)



## Selection of the mobile phase

Generally, anion-exchange columns are operated with the mobile phase pH at least one point higher than the isoelectric point (pI) of samples, while cation-exchange columns are operated with the mobile phase pH at least one point lower than the pI. The elution force of bivalent ions such as Bis-tris HCl is stronger than univalent ions such as Tris HCl.

Table. Buffer type and pH

Anion exchange (DEAE, QA)	pH	Cation exchange (CM, SP)
	4.0	Formic acid buffer
Piperazine buffer	5.0	Acetic acid buffer
Bis-Tris buffer	6.5	Phosphoric acid buffer
Tris buffer	8.0	HEPES buffer
Monoethanolamine buffer	9.5	

The initial mobile Phase (A) is 20-50 mmol/l of one of the above mentioned buffer solutions and the final mobile phase (B) is the mobile phase (A) with an addition of 20-600 mmol/l of salt.

## Selection of salts

High concentration of salts, generally NaCl, is used in elution buffers. When stronger elution buffer is needed, CaCl<sub>2</sub> or MgCl<sub>2</sub> can be used for DEAE and QA columns, and Na<sub>2</sub>SO<sub>4</sub> can be used for CM and SP columns.

## Selection of organic solvents

Water miscible solvents can be used when the elution is not strong enough. COSMOGEL columns can be used with up to 20% of water miscible solvents such as acetonitrile and 2-propanol.

## Ordering information

- Analytical column (Particle size: 10 μm)

### COSMOGEL Stainless Packed Column

Product name	Column size I.D. x length (mm)	Product number
COSMOGEL DEAE Stainless Packed Column	7.5×75	43371-91
COSMOGEL QA Stainless Packed Column		43373-71
COSMOGEL CM Stainless Packed Column		43375-51
COSMOGEL SP Stainless Packed Column		43377-31

### COSMOGEL Glass Packed Column

Product name	Column size I.D. x length (mm)	Product number
COSMOGEL DEAE Glass Packed Column	8.0×75	37845-81
COSMOGEL QA Glass Packed Column		37846-71
COSMOGEL CM Glass Packed Column		37844-91
COSMOGEL SP Glass Packed Column		37847-61

## 4) Hydrophobic interaction chromatography

# HIC

COSMOSIL 5HIC is designed for one step desalting and separation of proteins. Hydrophobic interaction chromatography (HIC) is an effective method for purification and separation of proteins (especially enzymes) based on differences in their surface hydrophobicity. Since this method does not use organic solvents like reversed phase chromatography, there is only a little loss in enzyme activity and the tertiary structure of proteins.


## Material characteristics

Packing material	HIC
Silica gel	High purity spherical porous silica
Average particle size	5 μm
Average pore size	approx. 300Å
Specific surface area	approx. 150m <sup>2</sup> /g
Main interaction	Hydrophobic interaction


## Application data

A buffer with high salt concentration, usually 1-2 mol/l of (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, is used as an initial mobile phase for adsorption of samples to a weakly hydrophobic stationary phase. The elution is done with a decreasing salt gradient.

- Separation of protein standards



- Separation of β-Glucosidase



## Ordering information

- Analytical column (Particle size: 5 μm)

### COSMOSIL 5HIC Packed Column






Column size I.D. x length (mm)	Product number
4.6×50	04263-21

# 10. Special column for fullerenes


## Introduction

Separation of fullerenes, especially preparative-scale separation, on conventional HPLC columns is always problematic due to the low solubility and small recovery of fullerenes. COSMOSIL offers a variety of columns designed for preparative-scale separation of fullerenes including higher fullerenes, metallofullerenes and fullerene derivatives.

## Material characteristics

Packing material	Buckyprep	Buckyprep-M	PBB	PYE	NPE
Silica gel	High purity porous spherical silica				
Average particle size	5 μm				
Average pore size	approx. 120Å				
Specific surface area	approx. 300 m <sup>2</sup> /g				
Stationary phase					
	Pyrenylpropyl group	Phenothiazinyl group	Pentabromobenzyl group	Pyrenylethyl group	Nitrophenylethyl group
Bonding type	Monomeric				
End capping treatment	Near-perfect	None	Near-perfect treatment		
Carbon content	approx. 17%	approx. 13%	approx. 8%	approx. 18%	approx. 9%
Feature	Standard column for fullerene separation.	Designed to separate metallofullerenes.	Designed for preparative separation of C <sub>60</sub> , C <sub>70</sub> .	Separation of fullerene and structural isomers.	Separation of fullerene derivatives

## Comparison of retention in toluene




## Solubility and boiling point of each solvent for C<sub>60</sub>

Solubility and boiling point of each solvent for C<sub>60</sub>

Solvent	mg/ml	b.p.(°C)
Methanol	0.001	64.5
Acetonitrile	0.018	81.8
n-Hexane	0.046	68.7
Toluene	3.2	111
Chlorobenzene*	7.0	132
Carbon disulfide	12	46.3
o-Dichlorobenzene*	27	180
1,2,4-Trichlorobenzene	21.3	213

\*: R.S.Ruoff, et al., J.Phys.Chem.,97,3379(1993)

## Comparison of separation property



## Suggested solvents

Solvent	Feature
Chlorobenzene	Stronger eluent than toluene. Recommended for higher fullerenes.
o-Dichlorobenzene	Stronger eluent than chlorobenzene.
1,2,4-Trichlorobenzene	Strongest eluent. It can be used as a washing solvent for higher fullerenes. To wash a column, inject 3 ml of 1,2,4-trichlorobenzene to a 4.6 mm I.D. × 250 mm column and 50 ml to a 20 mm I.D. × 250 mm column after every operation.
n-Hexane	Weak eluent. Recommended for weakly retained fullerenes.
Acetonitrile	Weak eluent. Recommended for weakly retained fullerenes.


Note : Use them after filtration or distillation, if they are not for HPLC.

# Buckyprep


COSMOSIL Buckyprep is a pyrenylpropyl group bonded silica based column specifically designed for fullerene separation. The unparalleled separation capabilities have enabled COSMOSIL Buckyprep to become the world benchmark of HPLC column for fullerene separation. COSMOSIL Buckyprep retains fullerenes very strongly with a mobile phase of 100% toluene and exceeds the injection volume of a standard C<sub>18</sub> column by a factor of 35. Therefore, preparative-scale separation can be obtained with a 250 mm × 4.6 mm I.D. analytical column.

## Application data

### Higher fullerenes




### Detivatized fullerene



## Fullerene Chromatogram Index

A comprehensive index with more than 100 chromatograms for fullerene separation is available from Nacalai Tesque, Inc. Please feel free to contact us for a copy.



## Ordering information

### Analytical / Preparative column (Particle size: 5 µm)



COSMOSIL Buckyprep Packed Column		COSMOSIL Buckyprep Guard Column	
Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
4.6×250	37977-61	4.6×10	37983-71
10×250	37981-91	10×20	37984-61
20×250	37982-81	20×50	34374-41
28×250	34346-11	28×50	05871-21

# Buckyprep-M

COSMOSIL Buckyprep-M is a phenothiazinyl group bonded silica based column specifically designed for metallofullerene separation. Metallofullerenes are retained more strongly than other fullerenes on this column. COSMOSIL Buckyprep-M is also effective for the separation of higher fullerenes and fullerene derivatives.

## Application data

### Metallo fullerenes



## Ordering information


### Analytical / Preparative column (Particle size: 5 µm)

COSMOSIL Buckyprep-M Packed Column		COSMOSIL Buckyprep-M Guard Column	
Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
4.6×250	04138-71	4.6×10	04139-61
10×250	04141-11	10×20	04140-21
20×250	04142-01	20×50	34474-31
28×250	05873-01	28×50	05872-11

# PBB


COSMOSIL PBB is a pentabromobenzyl group bonded silica based column useful for preparative-scale separation of fullerenes. It can be used with *o*-Dichlorobenzene, which has greater solubility for fullerenes than toluene. The loading capacity of COSMOSIL PBB for C<sub>60</sub> and C<sub>70</sub> can be three times greater than COSMOSIL Buckyrep.

## Separation of fullerenes with different mobile phases



## Preparative-scale separation

The loading capacity of COSMOSIL PBB for C<sub>60</sub> and C<sub>70</sub> can be three times greater than COSMOSIL Buckyrep.



## Ordering information


• Analytical / Preparative column (Particle size: 5 μm)

COSMOSIL 5PBB Packed Column		COSMOSIL 5PBB Guard Column	
Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
4.6×250	37980-01	4.6×10	37987-31
10×250	37985-51	10×20	37988-21
20×250	37986-41	20×50	34375-31

# PYE • NPE

## Application data

• PCBM, PCBB



## Ordering information

• Analytical / Preparative column (Particle size: 5 μm)

COSMOSIL 5PYE Packed Column		COSMOSIL 5PYE Guard Column	
Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
4.6×250	37989-11	4.6×10	37903-11
10×250	37996-11	10×20	38041-71
20×250	38044-41	20×50	34475-21
28×250	34300-91		

COSMOSIL 5NPE Packed Column		COSMOSIL 5NPE Guard Column	
Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
4.6×150	37902-21	4.6×10	37904-01
4.6×250	37990-71	10×20	38045-31
10×250	05469-11	20×50	05869-71
20×250	38046-21		

# 11. Special columns for carbon nanotubes

## CNT-300, CNT-1000, CNT-2000

COSMOSIL CNT series are ideal for separation of soluble carbon nanotubes based on sizes. COSMOSIL CNT series are packed with hydrophilic group-bonded silica packing material. The columns are specially designed to avoid adsorption of carbon nanotubes to silica support and thus ensure high resolution and maximum recovery of carbon nanotubes. COSMOSIL CNT series are available in three different pore sizes, 300 Å, 1000 Å and 2000 Å, respectively.


### Material characteristics

Packing material	CNT-300	CNT-1000	CNT-2000
Silica gel	High purity porous spherical silica		
Average particle size	5 µm		
Average pore size	approx. 300 Å	approx. 1000 Å	approx. 2000 Å
Stationary phase	Hydrophilic group (neutral)		
pH range	2.0-7.5		
Pressure	15 MPa and below		

### Application data

#### • Carbon nanotubes

COSMOSIL CNT columns offered improved separation for DNA wrapped carbon nanotubes by connecting three columns with different pore sizes.



### Ordering information

#### • Analytical column (Particle size: 5 µm)

COSMOSIL CNT-300 Packed Column		COSMOSIL CNT-300 Guard Column	
Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
7.5×300	09195-71	7.5×50	09194-81

COSMOSIL CNT-1000 Packed Column		COSMOSIL CNT-1000 Guard Column	
Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
7.5×300	09197-51	7.5×50	09196-61

COSMOSIL CNT-2000 Packed Column		COSMOSIL CNT-2000 Guard Column	
Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
7.5×300	09199-31	7.5×50	09198-41

# 12. Conventional columns versus high performance columns

### Introduction

A period of more than 30 years has passed since the first COSMOSIL 5C<sub>18</sub> columns were developed and offered for sale. Continuous technical improvement has made many of these columns obsolete and not of the highest quality and performance available any more. However, many long-term users continue to employ these older conventional columns for routine analysis and quality control. Nevertheless, the manufacture of these older columns will eventually cease and we strongly urge customers to replace the conventional columns with their higher performance equivalents outlined in the table below. For additional information, contact the manufacturer or your local distributor directly.

Conventional columns (old)		High performance columns (new)
5C <sub>18</sub> -AR	→	5C <sub>18</sub> -AR-II
5C <sub>18</sub>	→	5C <sub>18</sub> -MS-II
5C <sub>18</sub> -MS	→	5C <sub>18</sub> -MS-II
5C <sub>18</sub> -P	→	5C <sub>18</sub> -PAQ
5C <sub>18</sub> -P-MS	→	5C <sub>18</sub> -PAQ
5C <sub>8</sub>	→	5C <sub>8</sub> -MS
5TMS	→	5TMS-MS
5PE	→	5PE-MS
5CN-R	→	5CN-MS
5NH <sub>2</sub>	→	5NH <sub>2</sub> -MS
5C <sub>18</sub> -300	→	5C <sub>18</sub> -AR-300
5C <sub>8</sub> -300	→	5C <sub>8</sub> -AR-300
5C <sub>4</sub> -300	→	5C <sub>4</sub> -AR-300
5SL	→	5SL-II

### Ordering information

Product name	Column size I.D. x length (mm)	Product number
COSMOSIL 5C <sub>18</sub> Packed Column	4.6×150	39047-81
	4.6×250	39265-21
COSMOSIL 5C <sub>18</sub> -MS Packed Column	4.6×150	37971-21
	4.6×250	37972-11
COSMOSIL 5C <sub>18</sub> -AR Packed Column	4.6×150	37861-61
	4.6×250	37862-51
COSMOSIL 5C <sub>18</sub> -P Packed Column	4.6×150	39103-31
	4.6×250	39280-11
COSMOSIL 5C <sub>18</sub> -P-MS Packed Column	4.6×150	37995-21
	4.6×250	37994-31
COSMOSIL 5C <sub>8</sub> Packed Column	4.6×150	39042-31
	4.6×250	39260-71
COSMOSIL 5TMS Packed Column	4.6×150	39057-51
	4.6×250	39275-91
COSMOSIL 5CN-R Packed Column	4.6×150	39114-91
	4.6×250	39285-61
COSMOSIL 5NH <sub>2</sub> Packed Column	4.6×150	39150-11
	4.6×250	39290-81
COSMOSIL 5C <sub>18</sub> -300 Packed Column	4.6×150	39607-41
	4.6×150	39037-11
COSMOSIL 5SL Packed Column	4.6×150	39037-11
	4.6×250	39255-51


For more information on other columns, please feel free to contact us.



# LIQUID CHROMATOGRAPHY RELATED PRODUCT



1. Packing materials for column chromatography .....	64
C <sub>18</sub> -OPN .....	65
C <sub>18</sub> -PREP .....	68
SL-II-PREP .....	69
Silica gel (spherical • neutral) .....	70
Silica gel (for column chromatograph) .....	71
2. Liquid chromatography related products (for mobile phase) .....	72
Ion-pair reagents .....	72
3. Liquid chromatography related products (for pretreatment) .....	74
Cosmonice filter .....	74
Cosmospin filter .....	74
COSMOSIL HPLC accessories .....	76



# 1. Packing materials for column chromatography

## Introduction


Open column chromatography is an excellent and easy technique for large-scale preparation and purification at low cost. COSMOSIL offers both normal and reversed phase packing materials based on totally porous spherical silica, which provides higher separation, less pressure and higher reproducibility than irregular silica.

## Material characteristics


Packing material	C <sub>18</sub> -OPN	C <sub>18</sub> -PREP	SL-II-PREP	Silica gel 60(neutral)
Silica gel	High purity porous spherical silica			
Average particle size	75 · 140 μm	40 · 75 · 140 μm	75 · 140 μm	
Average pore size	approx. 120 Å			approx. 60 Å
Specific surface area	approx. 300 m <sup>2</sup> /g			approx. 500 m <sup>2</sup> /g
Stationary phase	Octadecyl group		None	
Carbon content	—	approx. 19%	0%	
End-capping treatment	Treated	None	—	
Useful range	Open column chromatography / Flash column chromatography			
	Reversed phase chromatography		Normal phase chromatography	

For more informations on other silica gel, please refer to page 71.

## Selection guide (Reversed phase)



## Selection guide (Normal phase)



## C<sub>18</sub>-OPN

Conventional reversed phase C<sub>18</sub> packing materials are restricted to about 30-50% water in the mobile phase. The COSMOSIL C<sub>18</sub>-OPN is a new "Water-Wet" C<sub>18</sub> packing material developed for reversed phase open column chromatography. The C<sub>18</sub>-OPN material can be used in 100% aqueous effluents.

## Characteristic

The external surface of the C<sub>18</sub>-OPN gel is coated with hydrophilic group to increase wettability of the gel, and octadecyl group is bonded in the pore of the gel. This physical characteristic of the gel makes the reversed phase open column chromatography possible with 100 % water.




Figure 1. Structure of C<sub>18</sub>-OPN



Figure 2. Packing material in water

Left : C<sub>18</sub>OPN provides good resolution  
Right : C<sub>18</sub>PREP float up

## Application data

### • Separation of Theobromine and Theophylline




Figure shows the sample are clearly separated by reversed open column chromatography with 70% of water.

## Influence of particle size





Table. Comparison between 75 μm and 140 μm particle size silica.

## Flow rate




Since reversed phase chromatography generally employs high viscosity solvents such as water and methanol, the flow rate is lower than that of normal phase chromatography. The flow rate of reversed phase depends on the mobile phase composition. Figure indicates that the flow rate of the COSMOSIL 140C<sub>18</sub>-OPN (140 μm in particle size) is about 2.5 times higher than that of the COSMOSIL 75C<sub>18</sub>-OPN.

Figure. Concentration of methanol against flow rate  
 Column size: 10 mm I.D. x 180 mm bed height (gravitational liquid flow)

## Comparison of normal phase

- Separation of *p*-Cresol and *p*-Ethylphenol by normal and reversed phase mode.

Since the structural difference between *p*-Cresol and *p*-Ethylphenol is only one methylene group, it is difficult to separate such samples under normal phase condition. On the other hand, the samples are clearly separated under reversed phase condition with COSMOSIL C<sub>18</sub>-OPN packing material.



## Column size and required amount of packing material

Table. Column size and required amount of C<sub>18</sub>-OPN packing material

Column I.D. (mm)	Bed height (mm)	Amount of C <sub>18</sub> -OPN(g)
10	150	4
	250	7
20	150	17
	250	28
30	150	38
	250	63

## Reproducibility and washing methods

Wash the COSMOSIL C<sub>18</sub>-OPN packing material with tetrahydrofuran, chloroform or other solvents to remove the impurities. This packing material has excellent reproducibility and can be used repeatedly.

### “CAUTION”

Do not wash with basic solvents of pH 7 or more which will dissolve the silica gel or pH 2 or less which will cleave the C<sub>18</sub> stationary phase. Dry the packing material at 50°C or less. See end of this chapter for packing method.

## Ordering information

- COSMOSIL C<sub>18</sub>-OPN

Product name	Average particle size	Product number	PKG size
COSMOSIL 75C <sub>18</sub> -OPN	75 μm	37842-66	100 g
		37842-95	500 g
		37842-11	1 kg
COSMOSIL 140C <sub>18</sub> -OPN	140 μm	37878-16	100 g
		37878-45	500 g
		37878-61	1 kg

# C<sub>18</sub>-PREP

The large particle size C<sub>18</sub> bulk materials are widely used for lab to process scale purifications. COSMOSIL offers three different particle sizes of C<sub>18</sub> packing materials.

## Particle size, flow rate and theoretical plate number

Because reversed phase chromatography employs effluents of high viscosity such as methanol and water, the flow rate is lower than that of normal phase chromatography, which uses effluents of low viscosity such as hexane and ethyl acetate.




Figure 1. Flow rate against theoretical plate number  
Column size: 20 mm I.D. x 300 mm




Figure 2. Concentration of methanol against flow rate  
Column size: 10 mm I.D. x 180 mm bed height (gravitational liquid flow)


## Application data

### • Vitamin E

#### COSMOSIL Application Data

Column: 40C<sub>18</sub>-PREP  
Column size: 20mm I.D.-300mm (Closed column)  
Mobile phase: Methanol  
Flow rate: 9.9 ml/min  
Temperature: Room temperature  
Detection: UV280nm

Sample: 1; DL- $\alpha$ -Tocopherol (5mg)  
2; DL- $\alpha$ -Tocopherol Acetate (5mg)




### • Natural compounds

#### COSMOSIL Application Data

Column: 40C<sub>18</sub>-PREP  
Column size: 20mm I.D.-300mm (Closed column)  
Mobile phase: Methanol/ 0.05%TFA-H<sub>2</sub>O = 70/30  
Flow rate: 9.9 ml/min  
Temperature: Room temperature  
Detection: UV254nm

Sample: 1; Baicalin (40 $\mu$ g)  
2; Baicalein (120 $\mu$ g)  
3; Wogonin (40 $\mu$ g)



## Ordering information

### • COSMOSIL C<sub>18</sub>-PREP

Product name	Average particle size	Product number	PKG size
COSMOSIL 40C <sub>18</sub> -PREP	40 $\mu$ m	37932-86	100 g
		37932-15	500 g
		37932-31	1 kg
COSMOSIL 75C <sub>18</sub> -PREP	75 $\mu$ m	37933-76	100 g
		37933-05	500 g
		37933-21	1 kg
COSMOSIL 140C <sub>18</sub> -PREP	140 $\mu$ m	37934-66	100 g
		37934-95	500 g
		37934-11	1 kg

# SL-II-PREP

COSMOSIL SL-II-PREP is ultra pure silica gel packing material more than 99.99% purity. COSMOSIL SL-II-PREP provides improved separation and reproducibility for compounds with carbonyl or phenol hydroxyl groups, which are often problematic on conventional silica gel materials.

\*All chromatograms shown below are obtained with silica gel packed into stainless steel columns.

## Performance for chelating compounds


Highly purified silica gel of COSMOSIL SL-II-PREP enables separation of metal coordination compounds without adsorption.

### • Metal Coordination Compounds

#### COSMOSIL Application Data

Column: 10mm I.D.-250mm  
Mobile phase: Hexane/Ethanol = 95/5  
Flow rate: 5.0 ml/min  
Temperature: 30°C  
Detection: UV254nm

Sample: 1; Quinizarin  
2; p-Nitrobenzyl Alcohol




### • Organic Acid and Amide

#### COSMOSIL Application Data

Column: 10mm I.D.-250mm  
Mobile phase: Hexane/Ethanol = 90/10  
Flow rate: 5.0 ml/min  
Temperature: 30°C  
Detection: UV254nm

Sample: 1; Salicylic Acid  
2; Salicylamide



## Ordering information

### • COSMOSIL SL-II

Product name	Average particle size	Product number	PKG size
COSMOSIL 75SL-II-PREP	75 $\mu$ m	38012-64	100 g
		38012-35	500 g
		38012-51	1 kg
COSMOSIL 140SL-II-PREP	140 $\mu$ m	38013-54	100 g
		38013-25	500 g
		38013-41	1 kg



# Silica gel (spherical • neutral)

Since conventional silica gels are weakly acidic, some pH sensitive compounds may be decomposed during the purification by column chromatography with the acidic silica gels. The pH of Silica gel 60 (spherical • neutral) is adjusted to nearly neutral for the separation of not only pH sensitive compounds but also new compounds that the physical properties are still unknown.

\*All chromatograms shown below are obtained with silica gel packed into stainless and steel columns.

## Comparison with conventional silica gel

### • Purification of Acetal



## Ordering information

### • Silica gel (spherical • neutral)

Product name	Average particle size	Product number	PKG size
Silica gel 60 (spherical • neutral) for column chromatograph	75 µm	30511-64	100 g
		30511-35	500 g
		30511-51	1 kg
		30511-06	5 kg
		30511-22	25 kg
	140 µm	30518-94	100 g
		30518-65	500 g
		30518-81	1 kg
		30518-52	25 kg

# Silica gel (for column chromatograph)

## Ordering information

### • Silica gel (spherical)

Product name	Particle size	Pore size	Grade	Product number	PKG size
Silica Gel 60, spherical	approx. 70 ~ 230 mesh	60 Å	SP	30731-71	1 kg
				30731-42	25 kg
	approx. 150 ~ 325 mesh		SP	30733-51	1 kg
				30733-22	25 kg
Silica Gel 120, spherical	approx. 70 ~ 230 mesh	120 Å	SP	30734-41	1 kg
	approx. 150 ~ 325 mesh			30735-31	1 kg

### • Silica gel (irregular)

Product name	Particle size	Pore size	Grade	Product number	PKG size
Silica Gel 60	approx. 70 ~ 230 mesh	60 Å	SP	30724-55	500 g
				30724-71	1 kg
				30724-84	5 kg
				30724-42	25 kg
	approx. 230 ~ 400 mesh		SP	30721-85	500 g
				30721-01	1 kg
				30721-14	5 kg
				30721-72	25 kg

## 2. Liquid chromatography related products (for mobile phase)


### Ion-pair reagents

The use of ion pair reagents as mobile phase additives extends the applicability of reversed phase HPLC. Ionic or highly polar compounds are difficult to analyze by reversed phase using only organic solvent and buffer solution because of the short retention time. Ion pair reagents are strong hydrophobic ions which form neutral ion pairs with oppositely charged samples molecules, making the efficient ODS columns amenable to separate ionic or highly polar samples. Nacalai Tesque offers a broad range of ion pair reagents for pharmaceutical compounds and other highly polar materials.

#### General use of ion-pair reagents in the mobile phase

When using ion pair reagents, ample time should be allowed for establishing equilibrium and for cleaning the column.


When using ion pair reagents with an alkyl chain of C<sub>10</sub> or shorter, it typically takes 20 minutes for establishing equilibrium and 30 minutes for cleaning. It may take more than 1 hour to clean the column when using ion pair reagents with an alkyl chain longer than C<sub>10</sub>. Therefore, it is highly recommended to prepare a column for exclusive use with ion pair reagents.



\*For HILIC column, please refer to page 38.

#### Separation adjustment

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV270nm  
 Sample: 1; L-Noradrenaline, 2; L-Adrenaline, 3; L-DOPA, 4; Dopamine Hydrochloride




#### Application data

- Low-molecular-weight unsaturated carboxylic acids
- Amino acids

##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 5mmol/l Tetra-n-butylammonium Phosphate, 20mmol/l Phosphate buffer(pH7) = 10 /90  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV230nm, 0.5AUFS


Sample: 1; Acrylic Acid (2.0µg), 2; Crotonic Acid (2.0µg), 3; Tiglic Acid (2.0µg), 4; Sorbic Acid (2.0µg)



##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 5mmol/l Sodium 1-Hexanesulfonate, 20mmol/l Phosphate buffer(pH3) = 10/90  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV215nm, 0.5AUFS

Sample: 1; Histidine (1.0µg), 2; Tyrosine (3.0µg), 3; Phenylalanine (1.0µg), 4; Tryptophan (0.5µg)



#### Ordering information

- for Basic samples

Product name	R:	Grade	Product number	PKG size
Sodium 1 - Butanesulfonate	C <sub>4</sub> H <sub>9</sub> -	SP	31331-94	5 g
Sodium 1 - Pentanesulfonate	C <sub>5</sub> H <sub>11</sub> -	SP	31730-64	5 g
			31730-22	25 g
Sodium 1 - Hexanesulfonate	C <sub>6</sub> H <sub>13</sub> -	SP	31529-24	5 g
			31529-82	25 g
Sodium 1 - Heptanesulfonate	C <sub>7</sub> H <sub>15</sub> -	SP	31528-34	5 g
			31528-92	25 g
Sodium 1 - Octanesulfonate	C <sub>8</sub> H <sub>17</sub> -	SP	31729-04	5 g
			31729-62	25 g
Sodium 1 - Nonanesulfonate	C <sub>9</sub> H <sub>19</sub> -	SP	31626-44	5 g
Sodium 1 - Decanesulfonate	C <sub>10</sub> H <sub>21</sub> -	SP	31429-34	5 g
Sodium 1 - Undecanesulfonate	C <sub>11</sub> H <sub>23</sub> -	SP	32030-04	5 g
Sodium 1 - Dodecane sulfonate	C <sub>12</sub> H <sub>25</sub> -	SP	31426-64	5 g
Sodium Lauryl Sulfate	**	SP	31623-32	25 g

#### 0.5M Solution

Sodium 1 - Butanesulfonate	C <sub>4</sub> H <sub>9</sub> -	SP	31332-84	5×10 ml
Sodium 1 - Hexanesulfonate	C <sub>6</sub> H <sub>13</sub> -	SP	31532-64	10 ml
			31532-06	5×10 ml
Sodium 1 - Octanesulfonate	C <sub>8</sub> H <sub>17</sub> -	SP	31733-34	10 ml
			31733-76	5×10 ml

- for Acid samples

Product name	X-	Grade	Product number	PKG size
Tetra - n - butylammonium Bromide	-Br	SP	32824-72	25 g
Tetra - n - butylammonium Chloride	-Cl	EP	32935-51	1 g
			32935-64	5 g
			32935-22	25 g
Tetra - n - butylammonium Hydrogensulfate	-HSO <sub>4</sub>	GR	32924-62	25 g
Tetra - n - butylammonium Iodide	-I	SP	32905-54	5 g
			32905-12	25 g
Tetra - n - butylammonium Perchlorate	-ClO <sub>4</sub>	SP	32906-44	5 g
			32906-02	25 g
Tetra - n - butylammonium Phosphate	-H <sub>2</sub> PO <sub>4</sub>	SP	32929-54	5 g

#### 0.5M Solution

Tetra - n - butylammonium Phosphate	-H <sub>2</sub> PO <sub>4</sub>	SP	32926-26	10 ml
			32926-84	5×10 ml

### 3. Liquid chromatography related products (for pretreatment)

## Cosmonice filter

Injection of samples containing particulates (microparticles, precipitates, colloid substances) will clog HPLC columns, shorten injector life, and result in extensive maintenance on pumps. Cosmonice filters are used to remove particulates from samples and prolong the life of HPLC system components. There are two types of Cosmonice filters as stated below.



#### W series (Aqueous solution)

W series are installed with low adsorption hydrophilic durapore filter (polyvinylidenedifluoride, PVDF). W series can be used with both aqueous and organic solvents. They are best suited for prefiltration of protein and other biological samples.

#### S series (Organic solvents)

S series are installed Teflon filter (polytetrafluoroethylene, PTFE) with strong resistance to organic solvents, acids and alkalis. They are best suited for prefiltration of samples with aggressive organic solvents such as chloroform and tetrahydrofuran.

Please refer to TECHNICAL NOTE 3, Sample pretreatment for HPLC at page 164.

#### Ordering information

##### • Cosmonice filter

Product name	Diameter (mm)	Pore size (µm)	Process volume	Hold-up volume	Product number	PKG size
Cosmonice Filter W (Aqueous)	4	0.45	less than 1 ml	< 10 µl	06543-04	100 pkg
	13	0.45	0.5 - 10 ml	< 30 µl	06544-94	100 pkg
	25	0.45	3 - 50 ml	< 100 µl	06545-84	50 pkg
Cosmonice Filter S (Solvent)	4	0.45	less than 1 ml	< 10 µl	06541-24	100 pkg
	13	0.45	0.5 - 10 ml	< 30 µl	06542-14	100 pkg

[Connection] Inlet: luer-lock, Outlet: luer-slip, Connectable needles

## Cosmospin fileter

Cosmospin filters are used to remove fine particles and precipitates from samples by centrifugation. They utilize omnipore hydrophilic PTFE membrane filter, which has a wide range of chemical resistance. Cosmospin filters are the best choice for HPLC sample filtration. Two pore sizes, G (0.2 µm) and H (0.45 µm), are available.



Please refer to TECHNICAL NOTE 3, Sample pretreatment for HPLC at page 164.

#### Ordering information

##### • Cosmospin filter

Product name	Pore size (µm)	Maximum sample volume	Hold-up volume	Maximum centrifugal force	Rotor size (fixed-angle)	Filtration area	Color	Product number	PKG size
Cosmospin Filter G	0.2	0.4 ml	5 µl	5,000× g	1.5 ml	0.2 cm <sup>2</sup>	brown	06549-44	100 pkg
Cosmospin Filter H	0.45	0.4 ml	5 µl	5,000× g	1.5 ml	0.2 cm <sup>2</sup>	white	06540-34	100 pkg

Dimension: Diameter 10.6 mm x Length 45 mm Membrane: Omnipore Hydrophilic PTFE Sample reservoir and collection tube: Polypropylene

### Chemical compatibility

Solvent	Cosmonice W series	Cosmonice S series	Cosmospin	Solvent	Cosmonice W series	Cosmonice S series	Cosmospin
Acetic acid, 98%	+	+	+	Hydrogen gas	+	+	+
Acetone	-	+	+	Hydrogen peroxide (3%)	+	+	
Acetonitrile	+	+	+	Hydraulic oil (5606)	+	+	+
Ammonia solution (6N)	+	+	+	Hypo (photo)	+	+	+
Ammonium hydroxide (conc.)	+	+	-	Isopropyl acetate	+	+	+
Amyl alcohol	+	+	+	Isopropyl alcohol	+	+	+
Benzene	+	+	-	Kerosene	+	+	+
Benzyl alcohol	+	+	-	Methanol	+	+	+
Boric acid	+		+	Methyl ethyl ketone	-	+	+
Butyl acetate		+		Methyl isobutyl ketone	+	+	-
Carbon tetrachloride	+	+	+	2-Methyl-1-propanol	+	+	+
Chloroform	+	+	+	Nitric acid (6N)	+	+	
Cyclohexanone	-	+	-	Nitrobenzene	+	+	-
Dichloromethane	+	+	-	Ozone gas	-	+	-
Dimethylacetamide	-	+	+	Paraldehyde		+	
Dimethylformamide	+	+	+	Pentane	+	+	-
Dimethylsulfoxide	-	+	-	Petroleum ether	+	+	
Dioxane	+	+	+	Phenol (water saturation)	+	+	-
DMSO	-	+	-	Phosphate buffer solution	+		+
Ethers	+	+	+	2-Propanol	+	+	+
Ethyl acetate	+	+	+	Pyridine	-	+	+
Ethyl alcohol	+	+	+	Seawater	+	+	+
Ethyl cello solve	+	+	+	Silicone oils	+	+	+
Ethylene glycol	+	+	+	Sodium hydroxide (conc.)	+	+	+
Formamide	+	+	+	Sulfuric acid (6N)		+	
Freon, TF or PCA solvent	+	+	+	Toluene	+	+	-
Gasoline	+	+	+	THF	-	+	-
Glycerine (Glycerol)	+	+	+	Trichloroacetic acid	+	+	+
Helium gas		+	+	Trichloroethane	+	+	-
Hexane	+	+	-	Trichloroethylene	+	+	-
Hydrochloride (6N)	+	+	+	TFA	+	+	-
Hydrofluoric acid	-	+	-	Xylene	+	+	+

+ : Recommended, - : Not recommended, (blank) : Not data available

# COSMOSIL HPLC accessories

I. COSMOSIL HPLC column

II. Liquid chromatography related product

III. Application data

IV. Technical note

V. Index

## Ordering information

### COSMOSIL Guard Cartridge Holder

Product name	Product number	PKG size
COSMOSIL Guard Cartridge Holder	38009-79	1 PKG



### COSMOSIL Column Prefilter

Product name	Product number	PKG size
COSMOSIL Column Prefilter	39361-19	1 PKG



### COSMOSIL Column Spare Filter for Prefilter

Product name	Product number	PKG size
COSMOSIL Column Spare Filter for Prefilter	39539-09	2 PKG



### COSMOSIL Column Connecting Tube

Product name	Product number	PKG size
COSMOSIL Column Connecting Tube	37843-69	1 PKG



# APPLICATION DATA

1. COSMOSIL Chromatogram Index .....	78
2. Application data of substances in Japanese Pharmacopoeia, 15 <sup>th</sup> version .....	79
3. COSMOSIL Application Data .....	80
1) Drugs .....	80
2) Crude Drugs .....	113
3) Natural Compounds .....	118
4) Pesticides .....	121
5) Food Additives .....	123
6) Vitamins .....	125
7) Metabolites .....	128
8) Carbohydrates .....	130
9) Lipids .....	133
10) Nucleic Acid Related Substances.....	135
11) Amino Acids, Peptides and Proteins.....	136
12) The others .....	143
4. Reference list .....	150

# 1. COSMOSIL Chromatogram Index

More than 6,100 single compound elution profiles with full chromatographic condition description are available. They are not only an incredible help for chromatographers, but also can be used as references in choosing conditions for similar compounds.

These data are available at our web site: <http://www.nacalai.com>

Category:    
 Alcohols   
 Aldehydes & Ketones   
 Amines & Amides

Column name:    
 C18-AR-II   
 C18-MS-II   
 C18-PAQ

Sample Name:  begins with

Molecular Formula:  contains (Keyword search)

CAS number:  (ex: 498-02-2)

Result/Page:

## Data example

**COSMOSIL Chromatogram Index**

Sample: Tolnafate   
 CAS No.: [2398-96-1]   
 Molecular formula: C<sub>19</sub>H<sub>17</sub>NOS   
 Column: 5C<sub>18</sub>-MS-II   
 Column size: 4.6mm I.D.-150mm   
 Mobile phase: Methanol/ H<sub>2</sub>O=80/20   
 Flow rate: 1.0 ml/min   
 Temperature: 30°C   
 Detection: UV254 nm   
 Attenuation: 0.128 aufs   
 Sample conc.: 0.09mg/ml   
 Injection volume: 1.0µl   
 Retention time: 6.92min   
 Capacity factor: 3.20

NACALAI TESQUE, INC

# 2. Application data of substances in Japanese Pharmacopoeia, 15<sup>th</sup> version

We prepare data of drugs using three kinds of C<sub>18</sub> columns that are specified in HPLC analysis in Application Data of Substances in Japanese Pharmacopoeia, 15<sup>th</sup> version. The data are available at our web site.

<http://www.nacalai.co.jp/en/cosmosil/TheJP15.htm>, or type "Cosmosil Japanese Pharmacopoeia" at a search site.

## Interpretation of application data

THE JAPANESE PHARMACOPOEIA (Fifteenth Edition)

**1 Acetaminophen (Purity)**

Acetaminophen [103-90-2]   
 C<sub>9</sub>H<sub>9</sub>NO<sub>2</sub> 151.16   
*N*-(4-Hydroxyphenyl) acetamide

	5C <sub>18</sub> -MS-II	5C <sub>18</sub> -AR-II	5C <sub>18</sub> -PAQ
3 Suitability	10.7	10.5	11.2
5			recommendation

4 HPLC chromatograms for 5C<sub>18</sub>-MS-II, 5C<sub>18</sub>-AR-II, and 5C<sub>18</sub>-PAQ columns.

6 Condition and JP description

Column size: 4.6 mm I.D.-150 mm   
 Mobile phase: Methanol : 0.05 mol/L potassium dihydrogenphosphate (pH 4.7) = 20 : 80   
 Temperature: 40°C   
 Detection: UV 225 nm   
 Flow rate: 1.0 ml/min   
 Sample: 1: *p*-Aminophenol Hydrochloride (0.02 mg/ml) [Internal standard]   
 2: Acetaminophen (0.02 mg/ml) [Purity sample]   
 3: 4'-Acetoxycetamide (0.02 mg/ml) [Impurity]

Test solution: Mobile phase   
 Injection volume: 10 µL

Items described in The Japanese Pharmacopoeia   
 (Flow rate) Adjust the flow rate so that the retention time of Sample 2 is about 5 minutes.   
 (Selection of column) (\*) Use a column eluting off Sample 1 before Sample 2, and resolution (R<sub>s</sub>) between these 2 peaks should not be less than 7.   
 (Detection sensitivity) Adjust the detection sensitivity so that the peak height of Sample 2 obtained from 10 µL of the standard solution is about 15% of the full scale.   
 (Time span of measurement) About 5 times as long as the retention time of Sample 2 after the solvent peak.   
 (Point to notice) \*2 0.05 mol/L potassium dihydrogenphosphate (pH 4.7) Dissolve 0.50 g of potassium dihydrogenphosphate in 900 mL of water, adjust pH to exactly 4.7 with dilute sodium hydrochloride test solution, and add water to make 1000 mL solution.

NACALAI TESQUE, INC

- ① Substance name
- ② Substance information
- ③ Suitability (○ : suitable, × : unsuitable, \*\* : depend on condition)
- ④ HPLC chromatogram
- ⑤ Recommending column
- ⑥ Condition and JP description

# 3. COSMOSIL Application Data

1) Drugs	P80	(Application data of substances in Japanese Pharmacopoeia ...P80-109, P113-117)
2) Crude Drugs	P113	
3) Natural Compounds	P118	
4) Pesticides	P121	
5) Food Additives	P123	
6) Vitamins	P125	
7) Metabolites	P128	
8) Carbohydrates	P130	
9) Lipids	P133	
10) Nucleic Acid Related Substances	P135	
11) Amino Acids, Peptides and Proteins	P136	
12) The others	P143	


## 1) Drugs

### • Acrinol

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 4.6mmol/l Sodium *I*-Octanesulfonate, 65mmol/l NaH<sub>2</sub>PO<sub>4</sub> (pH2.8 with H<sub>3</sub>PO<sub>4</sub>) = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV268nm

Sample: 1; Acrinol (0.04μg)




AP-0457

### • Azathioprine

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 25mmol/l KH<sub>2</sub>PO<sub>4</sub> (pH2.5 with H<sub>3</sub>PO<sub>4</sub>) = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV296nm

Sample: 1; 6-Mercaptopurine Hydrate (0.16μg)  
 2; Azathioprine (0.16μg)  
 3; Benzoic Acid (9.6μg)




AP-0471

### • Aztreonam

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Methanol/ 5mmol/l Tetra-*n*-butylammonium Hydrogensulfate (pH3.0 with 0.5mol/l Na<sub>2</sub>HPO<sub>4</sub>) = 35/65  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV280nm

Sample: 1; 4-Aminobenzoic Acid (0.5μg)  
 2; Aztreonam (5.0μg)




AP-0474

### • Acetaminophen

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 50mmol/l KH<sub>2</sub>PO<sub>4</sub> (pH4.7 with NaOH) = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV225nm

Sample: 1; *p*-Aminophenol Hydrochloride (0.2μg)  
 2; Acetaminophen (0.2μg)  
 3; 4'-Acetoxyacetanilide (0.2μg)



AP-0452


## 1) Drugs

### • Atenolol

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/Tetrahydrofuran/ 4.6mmol/l Sodium *I*-Octanesulfonate, 1.2mmol/l Tetra-*n*-butylammonium Hydrogensulfate, 25mmol/l KH<sub>2</sub>PO<sub>4</sub> (pH3.0 with H<sub>3</sub>PO<sub>4</sub>) = 9/1/40  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV226nm

Sample: 1; Atenolol (0.1μg)




AP-0468

### • Amikacin Sulfate

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Methanol/ 20mmol/l KH<sub>2</sub>PO<sub>4</sub> (pH6.5 with KOH) = 72/28  
 Flow rate: 1.5 ml/min  
 Temperature: 35°C  
 Detection: UV340nm

Sample: 1; Amikacin Sulfate 2,4,6-trinitrobenzenesulfonic Acid Derivative (0.6μg)  
 2; Kanamycin Sulfate 2,4,6-trinitrobenzenesulfonic Acid Derivative (0.6μg)




AP-0460

### • Meglumine Sodium Amidotrizoate

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 6.25mmol/l Tetra-*n*-butylammonium Phosphate, 50mmol/l K<sub>2</sub>HPO<sub>4</sub> (pH7.0 with H<sub>3</sub>PO<sub>4</sub>) = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Amidotrizoic Acid (0.25μg)  
 2; Acetizoic Acid (0.30μg)




AP-0739

### • Amlodipine Besilate

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 30mmol/l KH<sub>2</sub>PO<sub>4</sub> = 65/35  
 Flow rate: 0.5 ml/min  
 Temperature: 30°C  
 Detection: UV237nm

Sample: 1; Amlodipine Besilate (0.56μg)  
 2; Isobutyl *p*-Hydroxybenzoate (0.60μg)




AP-1084

### • Amoxicillin

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 10mmol/l CH<sub>3</sub>COONa (pH4.5 with CH<sub>3</sub>COOH) = 5/95  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV230nm

Sample: 1; Amoxicillin (3.0μg)




AP-0463

### • Ampicillin

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: 10% Acetonitrile/ 90% 50mmol/l (NH<sub>4</sub>)<sub>2</sub>HPO<sub>4</sub> (pH5.0 with H<sub>3</sub>PO<sub>4</sub>)  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV230nm

Sample: 1; Ampicillin Sodium Salt (10μg)  
 2; Guaiacol Glycerol Ether (5μg)




AP-0464

1) Drugs

• Isoxsuprine Hydrochloride

**COSMOSIL Application Data**


Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 32.6mmol/l (NH<sub>4</sub>)<sub>2</sub>HPO<sub>4</sub>, 18.4mmol/l Sodium *L*-Pentanesulfonate (pH2.5 with H<sub>3</sub>PO<sub>4</sub>) = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV269nm  
 Sample: 1; Isoxsuprine Hydrochloride (2.0µg)



• Isoniazid

**COSMOSIL Application Data**


Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: 10mmol/l-Sodium *L*-Tridecanesulfonate-Methanol/ 50mmol/l Phosphate buffer (pH2.5) =60/40  
 Flow rate: 0.5 ml/min  
 Temperature: 40°C  
 Detection: UV265nm  
 Sample: 1; Isonicotinic Acid (0.4µg)  
 2; Isoniazid (0.5µg)



• Idoxuridine

**COSMOSIL Application Data**


Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 13/87  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm  
 Sample: 1; 2'-Deoxyuridine (0.3µg)  
 2; 5-Iodouracil (1.2µg)  
 3; Idoxuridine (3.0µg)  
 4; Sulfathiazole (0.5µg)



• Idoxuridine

**COSMOSIL Application Data**


Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 4/96  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm  
 Sample: 1; 2'-Deoxyuridine (0.04µg)  
 2; 5-Iodouracil (0.12µg)  
 3; Idoxuridine (0.40µg)



• Ipratropium Bromides

**COSMOSIL Application Data**


Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 0.5%H<sub>3</sub>PO<sub>4</sub>/ Methanesulfonic Acid = 120/1000/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm  
 Sample: 1; Ipratropium Bromide Derivative  
 2; Ipratropium Bromide



• Indomethacin

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 0.1%H<sub>3</sub>PO<sub>4</sub> = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm  
 Sample: 1; 4-Chlorobenzoic Acid (1.0µg)  
 2; Butyl *p*-Hydroxybenzoate (0.6µg)  
 3; Indometacin(Indomethacin) (1.0µg)




1) Drugs

• Estradiol Benzoate

**COSMOSIL Application Data**


Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 35°C  
 Detection: UV230nm  
 Sample: 1; Progesterone (0.15µg)  
 2; Estradiol Benzoate (0.65µg)



• Estriol

**COSMOSIL Application Data**


Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 50/50  
 Flow rate: 0.5 ml/min  
 Temperature: 25°C  
 Detection: UV280nm  
 Sample: 1; Estriol (0.5µg)  
 2; Methyl Benzoate (0.5µg)



• Epirubicin Hydrochloride

**COSMOSIL Application Data**


Column: 5TMS-MS  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: 6.9mmol/l Sodium Lauryl Sulfate-H<sub>2</sub>O/Acetonitrile/Methanol/H<sub>3</sub>PO<sub>4</sub> = 540/290/170/1  
 Flow rate: 2.0 ml/min  
 Temperature: 35°C  
 Detection: UV254nm  
 Sample: 1; 2-Naphthalenesulfonic Acid Sodium(5.0µg)  
 2; Doxorubicin Hydrochloride (5.0µg)  
 3; Epirubicin Hydrochloride (10µg)  
 4; Daunorubicin Hydrochloride (5.0µg)



• Oxytocin

**COSMOSIL Application Data**


Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A; 100mmol/l NaH<sub>2</sub>PO<sub>4</sub>  
 B; Acetonitrile/H<sub>2</sub>O = 50/50  
 B conc. 30→60%(0→30min)  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV220nm  
 Sample: 1; (Arg<sup>8</sup>)-Vasopressin (5.0µg)  
 2; Oxytocin Acetate Salt (5.0µg)



• Omeprazole

**COSMOSIL Application Data**


Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 7.9mmol/l Na<sub>2</sub>HPO<sub>4</sub>, 1.35mmol/l NaH<sub>2</sub>PO<sub>4</sub>(pH7.6 with 1% H<sub>3</sub>PO<sub>4</sub>) = 11/29  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm  
 Sample: 1; Omeprazole (0.10µg)  
 2; 1,2-Dinitrobenzene (0.25µg)



• Captopril

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 0.1%Acetic Acid = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV220nm  
 Sample: 1; Captopril (2.6µg)




## 1) Drugs

### • Gabexate Mesilate

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ (0.1% Sodium Lauryl Sulfate/ 0.5% Sodium *I*-Heptanesulfonate/ Acetic Acid = 200/20/1) = 71/29  
 Flow rate: 0.5 ml/min  
 Temperature: 25°C  
 Detection: UV245nm

Sample: 1; Ethyl *p*-Hydroxybenzoate (0.39µg)  
 2; Butyl *p*-Hydroxybenzoate (0.39µg)  
 3; Gabexate Mesilate (0.75µg)




AP-0673

### • Camostat Mesilate

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ (0.2% Sodium *I*-Heptane Sulfonate/ 0.1% Sodium Lauryl Sulfate/ Acetic Acid = 1000/500/10) = 55/45  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV265nm

Sample: 1; Camostat Mesilate (1.02µg)  
 2; Butyl *p*-Hydroxybenzoate (0.54µg)




AP-0528

### • Carbazochrome Sodium Sulfonate

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Ethanol/ 10.4mmol/l (NH<sub>4</sub>)<sub>2</sub>HPO<sub>4</sub> = 75/925 (pH7 with H<sub>3</sub>PO<sub>4</sub>)  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV360nm

Sample: 1; Carbazochrome Sodium Sulfonate (1.0µg)  
 2; Carbazochrome (1.0µg)




AP-0535

### • Carbidopa

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Ethanol/ 50mmol/l NaH<sub>2</sub>PO<sub>4</sub> = 5/95 (pH2.7 with H<sub>3</sub>PO<sub>4</sub>)  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm

Sample: 1; Methyldopa (10µg)  
 2; Carbidopa (10µg)




AP-0538

### • Quinidine Sulfate and Quinine Hydrochloride

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 0.088% Acetic Acid, 0.16% Methanesulfonic Acid, 0.22% Diethylamine = 10/90  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV235nm

Sample: 1; Cinchonidine (1.0µg)  
 2; Quinidine Sulfate (10µg)  
 3; Quinine Hydrochloride (10µg)  
 4; Hydroquinidine Hydrochloride (10µg)  
 5; Dihydroquinine




AP-0828

### • Potassium Clavulanate

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: A; 50mmol/l NaH<sub>2</sub>PO<sub>4</sub> (pH4.0 with H<sub>3</sub>PO<sub>4</sub>)  
 B; Methanol/ 50mmol/l NaH<sub>2</sub>PO<sub>4</sub> (pH4.0 with H<sub>3</sub>PO<sub>4</sub>) = 50/50  
 B conc. 0% (0-4 min), 0→100% (4→15 min), 100% (15-25 min)  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV230nm

Sample: 1; Potassium Clavulanate  
 2; Sodium Amoxicillin



AP-0810


## 1) Drugs

### • Clarithromycin

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 67mmol/l KH<sub>2</sub>PO<sub>4</sub> = 35/65  
 Flow rate: 1.0 ml/min  
 Temperature: 50°C  
 Detection: UV210nm

Sample: 1; Clarithromycin (2.8µg)




AP-0578

### • Griseofulvin

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Griseofulvin (2.2µg)  
 2; Butyl *p*-Hydroxybenzoate (1.6µg)



AP-0676

### • Clindamycin Phosphate

#### COSMOSIL Application Data

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 100mmol/l KH<sub>2</sub>PO<sub>4</sub> (pH2.5 with H<sub>3</sub>PO<sub>4</sub>) = 22/78  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV210nm

Sample: 1; Lincomycin Hydrochloride (4.0µg)  
 2; Clindamycin Phosphate (4.0µg)  
 3; Methyl *p*-Hydroxybenzoate (0.3µg)




AP-0584

### • Clindamycin Hydrochloride

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 50mmol/l KH<sub>2</sub>PO<sub>4</sub> (pH7.5 with KOH) = 45/55  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV210nm

Sample: 1; Clindamycin Hydrochloride (15µg)




AP-0581

### • Cloxacillin Sodium

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 50mmol/l (NH<sub>4</sub>)<sub>2</sub>HPO<sub>4</sub> = 25/75 (pH4.0 with H<sub>3</sub>PO<sub>4</sub>)  
 Flow rate: 1.5 ml/min  
 Temperature: 25°C  
 Detection: UV230nm

Sample: 1; Guaiacol Glycerol Ether (5.0µg)  
 2; Cloxacillin Sodium Monohydrate (10µg)




AP-0585

### • Clofibrate

#### COSMOSIL Application Data

Column: 5CN-MS  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Hexane/2-Propanol/Acetic Acid = 1970/30/1  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV275nm

Sample: 1; Clofibrate (2.0µg)  
 2; 4-Chlorophenol (0.12µg)  
 3; 4-Ethoxyphenol (0.12µg)



AP-0588


## 1) Drugs

### • Clofibrate

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 0.1%Acetic Acid = 60/40  
 Flow rate: 1.5 ml/min  
 Temperature: 25°C  
 Detection: UV275nm

Sample: 1; 4-Chlorophenol (4.0µg)  
 2; Ibuprofen (60µg)  
 3; Clofibrate (10µg)




NACALAI TESQUE, INC  
 AP-0590

### • Clobetasol Propionate

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ Methanol/ 50mmol/l NaH<sub>2</sub>PO<sub>4</sub> (pH2.5 with H<sub>3</sub>PO<sub>4</sub>) = 47.5/10/42.5  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV240nm

Sample: 1; Clobetasol Propionate (1.0µg)  
 2; Beclomethasone Dipropionate (2.0µg)




NACALAI TESQUE, INC  
 AP-1098

### • Cloperastine Hydrochloride

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 100mmol/l KH<sub>2</sub>PO<sub>4</sub>, 0.16%Perchloric acid = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV222nm

Sample: 1; Cloperastine Hydrochloride (0.2µg)  
 2; Benzophenone (0.4µg)  
 3; 4-Chlorobenzophenone (0.2µg)




NACALAI TESQUE, INC  
 AP-0593

### • Chlorpheniramine Maleate

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 74.5mmol/l (NH<sub>4</sub>)H<sub>2</sub>PO<sub>4</sub>, 14.7mmol/l H<sub>3</sub>PO<sub>4</sub> = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV225nm

Sample: 1; Maleic Acid (Impurity)  
 2; Chlorpheniramine Maleate (0.06µg)




NACALAI TESQUE, INC  
 AP-0568

### • Chlorpropamide

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 1%Acetic Acid = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV240nm

Sample: 1; Chlorpropamide (1.0µg)




NACALAI TESQUE, INC  
 AP-0570

### • Chlormadinone Acetate

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/H<sub>2</sub>O = 65/35  
 Flow rate: 0.5 ml/min  
 Temperature: 30°C  
 Detection: UV236nm

Sample: 1; Butyl *p*-Hydroxybenzoate (0.2µg)  
 2; 17 $\alpha$ -Acetoxyprogesterone (0.2µg)  
 3; Chlormadinone Acetate (0.8µg)



NACALAI TESQUE, INC  
 AP-0563


## 1) Drugs

### • Ketoprofen

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/H<sub>2</sub>O/ 0.5mol KH<sub>2</sub>PO<sub>4</sub> (pH3.5 with H<sub>3</sub>PO<sub>4</sub>) = 43/55/2  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV233nm

Sample: 1; Ketoprofen (0.4µg)



NACALAI TESQUE, INC  
 AP-0709


### • Cortisone

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A; Acetonitrile/ H<sub>2</sub>O = 30/70  
 B; Acetonitrile/ H<sub>2</sub>O = 70/30  
 B conc. 10% (0-5min)  
 10→90% (5→25min)  
 90% (25-30min)

Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Cortisone (3.6µg)  
 2; Cortisone Acetate (3.75µg)  
 3; Hydrocortisone Acetate (3.6µg)




NACALAI TESQUE, INC  
 AP-0595

### • Salicylic Acid

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 1.67% Acetic Acid = 40/60  
 Flow rate: 0.5 ml/min  
 Temperature: 35°C  
 Detection: UV270nm

Sample: 1; *p*-Hydroxybenzoic Acid (0.025µg)  
 2; Phenol (0.10µg)  
 3; Salicylic Acid (0.50µg)




NACALAI TESQUE, INC  
 AP-1123

### • Santonin

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/H<sub>2</sub>O = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Santonin (0.50µg)  
 2; Ethyl *p*-Hydroxybenzoate (0.24µg)




NACALAI TESQUE, INC  
 AP-0849

### • Cyanocobalamins

#### COSMOSIL Application Data

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 70mmol/l Na<sub>2</sub>HPO<sub>4</sub> (pH3.5 with H<sub>3</sub>PO<sub>4</sub>) = 53/147  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV361nm

Sample: 1; Cyanocobalamin  
 2; Cyanocobalamin Derivative




NACALAI TESQUE, INC  
 AP-1101

### • Digitoxin

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 75/25  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV230nm

Sample: 1; Digoxin (1.0µg)  
 2; Gitoxin (0.5µg)  
 3; Digitoxin (0.5µg)  
 4; Acenaphthene (0.03µg)



NACALAI TESQUE, INC  
 AP-0620


## 1) Drugs

### • Cyclosporin

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 0.22% H<sub>3</sub>PO<sub>4</sub> /  
*tert*-Butyl Methyl Ether = 50/45/5  
 Flow rate: 1.0 ml/min  
 Temperature: 50°C  
 Detection: UV210nm

Sample: 1; Cyclosporin(Cyclosporin) (24µg)




NACALAI TESQUE, INC  
 AP-0574

### • Diclofenac Sodium

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 0.12%Acetic Acid = 70/30  
 Flow rate: 0.5 ml/min  
 Temperature: 40°C  
 Detection: UV240nm

Sample:  
 1; Ethyl *p*-Hydroxybenzoate (0.14µg)  
 2; Propyl *p*-Hydroxybenzoate (0.2µg)  
 3; Diclofenac Sodium (0.4µg)




NACALAI TESQUE, INC  
 AP-0613

### • Diclofenamide

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 40mmol/l Na<sub>2</sub>HPO<sub>4</sub>,  
 52mmol/l NaH<sub>2</sub>PO<sub>4</sub> = 50/50  
 Flow rate: 0.5 ml/min  
 Temperature: 30°C  
 Detection: UV280nm

Sample:  
 1; Diclofenamide (10µg)  
 2; Butyl *p*-Hydroxybenzoate (1.2µg)




NACALAI TESQUE, INC  
 AP-0615

### • Disulfiram

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV210nm

Sample: 1; Benzophenone (0.05µg)  
 2; Disulfiram (0.05µg)




NACALAI TESQUE, INC  
 AP-0625

### • Zidovudine

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV265nm

Sample: 1; Thymine (0.08µg)  
 2; Zidovudine (2.0µg)




NACALAI TESQUE, INC  
 AP-1130

### • Dipyrindamole

#### COSMOSIL Application Data

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 7.5mmol/l KH<sub>2</sub>PO<sub>4</sub> = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV280nm

Sample: 1; Dipyrindamole (2.8µg)  
 2; *p*-Terphenyl (1.2µg)



NACALAI TESQUE, INC  
 AP-0624


## 1) Drugs

### • Josamycin

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 0.8mol Sodium Perchlorate  
 (pH2.5 with HCl) = 40/60  
 Flow rate: 2.0 ml/min  
 Temperature: 40°C  
 Detection: UV231nm

Sample: 1; Josamycin (10µg)




NACALAI TESQUE, INC  
 AP-0706

### • Cilastatin Sodium Salt

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: A; 0.1% H<sub>3</sub>PO<sub>4</sub>  
 B; Acetonitrile/ 0.1% H<sub>3</sub>PO<sub>4</sub> = 30/70  
 B conc. 15→30%(0→30min), 100%(30-40min)  
 Flow rate: 2.0 ml/min  
 Temperature: 50°C  
 Detection: UV210nm

Sample:  
 1; Cilastatin Sodium Salt (1.0µg)




NACALAI TESQUE, INC  
 AP-0576

### • Diltiazem Hydrochloride

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/Acetonitrile/  
 118mmol/l CH<sub>3</sub>COONa,  
 6.5mmol/l *d*-Camphorsulfonic Acid  
 = 25/25/50(pH6.68)  
 Flow rate: 1.0 ml/min  
 Temperature: 50°C  
 Detection: UV240nm

Sample: 1; Diltiazem Hydrochloride (3.0µg)  
 2; Phenyl Benzoate (2.0µg)




NACALAI TESQUE, INC  
 AP-0621

### • Sulbactam Sodium

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 5mmol/l Tetra *n*-Butylammonium  
 Hydroxide(pH4.0 with H<sub>3</sub>PO<sub>4</sub>) = 25/75  
 Flow rate: 1.0 ml/min  
 Temperature: 35°C  
 Detection: UV220nm

Sample: 1; Sulbactam Sodium (10µg)  
 2; Ethyl *p*-Hydroxybenzoate (7.0µg)




NACALAI TESQUE, INC  
 AP-0853

### • Cetirizine Hydrochloride

#### COSMOSIL Application Data

Column: 5SL-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 40mmol/l Sulfuric Acid = 94/6  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV230nm

Sample: 1; Cetirizine Dihydrochloride (0.5µg)  
 2; 4-Dimethylaminoantipyrine (0.9µg)




NACALAI TESQUE, INC  
 AP-1096

### • Cetirizine Hydrochloride

#### COSMOSIL Application Data

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 1.7mmol/l Sodium *I*-Heptanesulfonate  
 = 42/58 (pH3.0 with 0.5mol/l Sulfuric Acid)  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV230nm

Sample: 1; Cetirizine Dihydrochloride (2.0µg)  
 2; Propyl *p*-Hydroxybenzoate (4.0µg)



NACALAI TESQUE, INC  
 AP-1097


1) Drugs

● Cefaclor

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 50mmol/l KH<sub>2</sub>PO<sub>4</sub> (pH3.4 with H<sub>3</sub>PO<sub>4</sub>) = 6/94  
 Flow rate: 2.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Cefaclor (2.0µg)  
 2; 4-Aminoacetophenone (2.9µg)




● Cefaclor

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: A; 50mmol/l NaH<sub>2</sub>PO<sub>4</sub>(pH4.0 with H<sub>3</sub>PO<sub>4</sub>)  
 B; Acetonitrile/Buffer A = 45/55  
 B conc. 5→25%(0-30min), 25→100%(30-45min), 100%(45-55min)

Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV220nm

Sample: 1; Cefaclor (1.0µg)



● Cefadroxil

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 10mmol/l KH<sub>2</sub>PO<sub>4</sub> = 60/340  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV262nm

Sample: 1; Cefadroxil (1.3µg)




● Cephalothin Sodium Salt

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/Ethanol/ 158mmol/l CH<sub>3</sub>COONa, 0.076%Acetic Acid(pH5.9 with NaOH) = 15/7/78

Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV254nm

Sample: 1; Similar compound (0.25µg)  
 2; Cephalothin Sodium Salt (0.25µg)




● Cefsulodin Sodium Salt

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A; Acetonitrile/1%(NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> = 3/97  
 B; Acetonitrile/1%(NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> = 8/92  
 B conc. 0%(0-14min), 100%(14-30min)

Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Isonicotinamide (0.2µg)  
 2; Cefsulodin Sodium Salt Hydrate (0.2µg)




● Ceftazidime

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 50mmol/l (NH<sub>4</sub>)<sub>2</sub>PO<sub>4</sub> (pH3.5 with H<sub>3</sub>PO<sub>4</sub>) = 13/87

Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Ceftazidime (2.5µg)  
 2; Acetanilide (2.5µg)




1) Drugs

● Ceftazidime

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 50mmol/l (NH<sub>4</sub>)<sub>2</sub>PO<sub>4</sub> = 30/70(pH7.0 with NH<sub>3</sub>)  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV254nm

Sample: 1; Ceftazidime (0.5µg)  
 2; Pyridine (0.5µg)




● Ceftriaxone Sodium

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: 18.1mmol/l tetra-*n*-heptylammonium bromide-Acetonitrile / 4.1mmol/l Na<sub>2</sub>HPO<sub>4</sub>, 2.6mmol/l KH<sub>2</sub>PO<sub>4</sub>, 0.88mmol/l Citric Acid, 1.8mmol/l NaOH = 45/55

Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Ceftriaxone Sodium (0.5µg)  
 2; Diethyl Terephthalate (0.9µg)




● Daunorubicin Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 38/62 (pH2.2 with H<sub>3</sub>PO<sub>4</sub>)  
 Flow rate: 0.5 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Daunorubicin Hydrochloride (5.0µg)  
 2; 2-Naphthalenesulfonic Acid (10µg)




● Dantrolene Sodium

**COSMOSIL Application Data**

Column: 5SL-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Hexane/Acetic Acid/Ethanol = 90/10/9  
 Flow rate: 2.0 ml/min  
 Temperature: 30°C  
 Detection: UV300nm

Sample: 1; Theophylline (1.0µg)  
 2; Dantrolene Sodium (0.06µg)




● Thiamine Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/Acetonitrile/ 5mmol/l Sodium *I*-Octanesulfonate, 1% H<sub>3</sub>PO<sub>4</sub> = 24/16/60  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Thiamine Hydrochloride (4.0µg)  
 2; Methyl Benzoate (20µg)




● Thiopental Sodiums

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 7mmol/l KH<sub>2</sub>PO<sub>4</sub> (pH3.0 with H<sub>3</sub>PO<sub>4</sub>) = 30/70  
 Flow rate: 1.5 ml/min  
 Temperature: 40°C  
 Detection: UV254nm

Sample: 1; Isopropyl *p*-Hydroxybenzoate (1.0µg)  
 2; Propyl *p*-Hydroxybenzoate (1.0µg)  
 3; Thiopental Sodium Isomer (4.8µg)  
 4; Thiopental Sodium (4.8µg)




1) Drugs

• Ticarcillin Sodium

**COSMOSIL Application Data**


Column: 5TMS-MS  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 5mmol/l Tetra-*n*-Butylammonium Bromide, 25mmol/l NaH<sub>2</sub>PO<sub>4</sub>(pH3.0 with H<sub>3</sub>PO<sub>4</sub>) / Acetic Acid = 225/1000/2.5  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV230nm  
 Sample: 1; *o*-Toluic Acid (2.5µg)  
 2; Ticarcillin Sodium (7.6µg)



• Titepidine Hibenzate

**COSMOSIL Application Data**


Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Tetrahydrofuran/ 1%Ammonium Acetate = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 50°C  
 Detection: UV254nm  
 Sample: 1; Hibenzic Acid  
 2; Titepidine  
 3; Propyl *p*-Hydroxybenzoate (0.6µg)



• Titepidine Hibenzate

**COSMOSIL Application Data**


Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 0.2%Ammonium Acetate = 65/35  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV254nm  
 Sample: 1; Hibenzic Acid  
 2; Titepidine  
 3; Xanthene (0.8µg)



• Titepidine Hibenzate

**COSMOSIL Application Data**


Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: (Acetonitrile/ 2-Propanol = 2/1) / 0.1%Sodium Lauryl Sulfate, 0.2%H<sub>3</sub>PO<sub>4</sub> = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV254nm  
 Sample: 1; Hibenzic Acid  
 2; Titepidine  
 3; Dibucaine Hydrochloride (1.6µg)



• Timolol Maleate

**COSMOSIL Application Data**


Column: SPE-MS  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/Methanol/ 5mmol/l Sodium *l*-Hexanesulfonate, 0.3%Triethylamine, 0.4%Formic Acid(pH3.0 with Formic Acid) = 5/25/70  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV280nm  
 Sample: 1; Maleic Acid  
 2; Timolol Maleate (7.5µg)



• Dexamethasone

**COSMOSIL Application Data**

Column: SPE-MS  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 21 mmol/l Ammonium Formate (pH3.6 with Formic Acid) = 33/67  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm  
 Sample: 1; Dexamethasone (0.18µg)




1) Drugs

• Deferoxamine Mesilate

**COSMOSIL Application Data**


Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: 2-Propanol/ 10.5mmol/l (NH<sub>4</sub>)<sub>2</sub>HPO<sub>4</sub>, 1.05mmol/l EDTA, 5.6mmol/l Sodium *l*-Heptanesulfonate (pH2.8 with H<sub>3</sub>PO<sub>4</sub>) = 10/90  
 Flow rate: 0.5 ml/min  
 Temperature: 25°C  
 Detection: UV230nm  
 Sample: 1; Deferoxamine Mesilate (6.4µg)  
 2; Methyl *p*-Hydroxybenzoate (1.6µg)



• Doxycycline Hydrochloride

**COSMOSIL Application Data**


Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ *N,N*-Dimethyloctylamine/ 100mmol/l NaH<sub>2</sub>PO<sub>4</sub> = 550/3/450 (pH8.0 with NaOH)  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm  
 Sample: 1; Doxycycline Hydrochloride (10µg)



• Doxifluridine

**COSMOSIL Application Data**


Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 35/65  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm  
 Sample: 1; Doxifluridine (0.5µg)  
 2; Caffeine (1.0µg)



• Doxorubicin Hydrochloride

**COSMOSIL Application Data**


Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 10.4mmol/l Sodium Lauryl Sulfate, 0.14%H<sub>3</sub>PO<sub>4</sub>=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm  
 Sample: 1; Doxorubicin Hydrochloride (5.0µg)



• Tocopherols

**COSMOSIL Application Data**


Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 98/2  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV292nm  
 Sample: 1;  $\alpha$ -Tocopherol (20µg)  
 2; Tocopherol Acetate (20µg)



• Tocopherol Nicotinate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol  
 Flow rate: 1.0 ml/min  
 Temperature: 35°C  
 Detection: UV264nm  
 Sample: 1; Nicotinic Acid (0.75µg)  
 2; Tocopherol (12.5µg)  
 3; Tocopherol Nicotinate (2.5µg)




## 1) Drugs

### • Todralazine Hydrochloride

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 5.4mmol/l Sodium *I*-Heptanesulfonate = 40/60 (pH3.5 with Acetic Acid)  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV240nm  
 Sample: 1; Potassium Biphthalate (0.4µg)  
 2; Todralazine Hydrochloride (0.4µg)




NACALAI TESQUE, INC  
AP-0879

### • Dopamine Hydrochloride

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: 200mmol/l Na<sub>2</sub>HPO<sub>4</sub> (pH3.0 with Citric Acid)  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV280nm  
 Sample: 1; Uracil (0.5µg)  
 2; Dopamine Hydrochloride (1.0µg)




NACALAI TESQUE, INC  
AP-0631

## 1) Drugs

### • Triamcinolone Acetonide

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 35/65  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV240nm  
 Sample: 1; Prednisolone (0.4µg)  
 2; Triamcinolone Acetonide (0.8µg)




NACALAI TESQUE, INC  
AP-0892

### • Trimetazidine Hydrochloride

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A; Methanol/ 14.2mmol/l Sodium *I*-Heptanesulfonate (pH 3.0 with 10% H<sub>3</sub>PO<sub>4</sub>) = 2/3  
 B; Methanol  
 B conc. 5%→25% 50min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV240nm  
 Sample: 1; *I*-(2,3,4-Trimethoxybenzyl)piperazine Dihydrochloride (40µg)




NACALAI TESQUE, INC  
AP-1125

### • Dobutamine Hydrochloride

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 20mmol/l Tartrate buffer (pH3.0) = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm  
 Sample: 1; Dobutamine Hydrochloride (11µg)  
 2; Salicylamide (8.5µg)




NACALAI TESQUE, INC  
AP-0628

### • Tolazamide

#### COSMOSIL Application Data

Column: 5SL-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Hexane/ Water-saturated Hexane/ Tetrahydrofuran/ Ethanol/ Acetic Acid = 475/ 475/ 20/ 15/ 9  
 Flow rate: 2.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm  
 Sample: 1; Tolbutamide (7.5µg)  
 2; Tolazamide (30µg)




NACALAI TESQUE, INC  
AP-0882

### • Trimetazidine Hydrochloride

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 50mmol/l KH<sub>2</sub>PO<sub>4</sub> (pH3.0 with H<sub>3</sub>PO<sub>4</sub>) = 15/85  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV230nm  
 Sample: 1; Trimetazidine Hydrochloride (0.15µg)  
 2; *p*-Hydroxybenzoic Acid (0.18µg)




NACALAI TESQUE, INC  
AP-0895

### • Tolnaftate

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm  
 Sample: 1; Diphenyl Phthalate (18µg)  
 2; Tolnaftate (4.0µg)




NACALAI TESQUE, INC  
AP-0885

### • Tranexamic Acid

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 152mmol/l NaH<sub>2</sub>PO<sub>4</sub>, 8mmol/l Sodium Lauryl Sulfate, 0.83%Triethylamine(pH2.5 with H<sub>3</sub>PO<sub>4</sub>) = 40/60  
 Flow rate: 0.5 ml/min  
 Temperature: 35°C  
 Detection: UV220nm  
 Sample: 1; Tranexamic Acid (6.0µg)  
 2; 4-(Aminomethyl) benzoic Acid(0.06µg)




NACALAI TESQUE, INC  
AP-0888

### • Triamcinolone

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 25/75  
 Flow rate: 0.5 ml/min  
 Temperature: 30°C  
 Detection: UV254nm  
 Sample: 1; Triamcinolone (1.0µg)  
 2; Methyl *p*-Hydroxybenzoate (0.4µg)




NACALAI TESQUE, INC  
AP-0889

### • Domperidone

#### COSMOSIL Application Data

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 15.6mmol/l K<sub>2</sub>HPO<sub>4</sub> (pH3.5 with 20mmol/l H<sub>3</sub>PO<sub>4</sub>) = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 35°C  
 Detection: UV287nm  
 Sample: 1; Domperidone (0.1µg)  
 2; Ethyl *p*-hydroxybenzoate (0.2µg)




NACALAI TESQUE, INC  
AP-1102

### • Naphazoline and Chlorpheniramine

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 0.2% Sodium Lauryl Sulfate, 0.1% H<sub>3</sub>PO<sub>4</sub> = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm  
 Sample: 1; *o*-Ethoxybenzamide (5.1µg)  
 2; Naphazoline Nitrate (2.2µg)  
 3; Chlorpheniramine Maleate (4.1µg)



NACALAI TESQUE, INC  
AP-0954


## 1) Drugs

### • Nabumetone

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ Tetrahydrofuran/ 0.1% Acetic Acid = 28/12/60  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV254nm

Sample: 1; Nabumetone (10.0µg)




AP-1113

### • Nalidixic Acid

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 40mmol/l NaH<sub>2</sub>PO<sub>4</sub> (pH2.8 with H<sub>3</sub>PO<sub>4</sub>) = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV260nm

Sample: 1; Methyl *p*-Hydroxybenzoate (0.13µg)  
 2; Nalidixic Acid (0.005µg)




AP-0770

### • Nicardipine Hydrochloride

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 0.05%Perchloric Acid = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Nicardipine Hydrochloride (0.4µg)  
 2; Nifedipine (0.4µg)




AP-0775

### • Nicardipine Hydrochloride

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 10mmol/l KH<sub>2</sub>PO<sub>4</sub> = 68/32  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV254nm

Sample: 1; Nicardipine Hydrochloride (0.4µg)  
 2; Di-*n*-butyl Phthalate (1.3µg)




AP-0778

### • Nicotinic Acid

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: 5mmol/l Sodium *I*-Octanesulfonate-Methanol/ 50mmol/l NaH<sub>2</sub>PO<sub>4</sub>(pH3.0) = 20/80  
 Flow rate: 0.5 ml/min  
 Temperature: 35°C  
 Detection: UV260nm

Sample: 1; Nicotinic Acid (1.0µg)  
 2; Caffeine (1.0µg)



AP-0788

### • Nicotinamide

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 5mmol/l Sodium *I*-Heptanesulfonate = 30/70  
 Flow rate: 0.5 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Nicotinic Acid (4.0µg)  
 2; Nicotinamide (0.8µg)



AP-0784


## 1) Drugs

### • Nicergolin

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/Acetonitrile/ 50mmol/l KH<sub>2</sub>PO<sub>4</sub> (pH7.0 with Triethylamine) = 35/30/35  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV288nm

Sample: 1; Nicergoline (20µg)




AP-0781

### • Nitrendipine

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/Tetrahydrofuran/H<sub>2</sub>O = 20/24/56  
 Flow rate: 0.5 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Propyl *p*-Hydroxybenzoate (0.15µg)  
 2; Nitrendipine (0.50µg)




AP-0790

### • Neostigmine Methylsulfate

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 20mmol/l NaH<sub>2</sub>PO<sub>4</sub> (pH3.0 with H<sub>3</sub>PO<sub>4</sub>), 5mmol/l Sodium *I*-Pentanesulfonate = 11/89  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV259nm

Sample: 1; Dimethylaminophenol (0.8µg)  
 2; Neostigmine Methylsulfate (5.0µg)




AP-0772

### • Baclofen

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 0.11%Acetic acid = 60/40  
 Flow rate: 0.5 ml/min  
 Temperature: 25°C  
 Detection: UV268nm

Sample: 1; Baclofen (5.0µg)  
 2; Methyl *p*-Hydroxybenzoate (0.063µg)




AP-0477

### • Haloperidol

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: 3.5mmol/l Sodium Lauryl Sulfate-Methanol/ 10mmol/l *tri*-Sodium Citrate(pH3.5 with HCl) = 75/25  
 Flow rate: 0.5 ml/min  
 Temperature: 40°C  
 Detection: UV220nm

Sample: 1; *4*-(4-Chlorophenyl)-4-hydroxypiperidine (0.20µg)  
 2; Haloperidol (0.20µg)  
 3; Biphenyl (0.16µg)



AP-0680


### • Vancomycin Hydrochloride

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: A: Acetonitrile/Tetrahydrofuran/ 0.2%Triethylamine (pH3.2 with H<sub>3</sub>PO<sub>4</sub>) = 7/1/92  
 B: Acetonitrile/Tetrahydrofuran/ 0.2%Triethylamine (pH3.2 with H<sub>3</sub>PO<sub>4</sub>) = 29/1/70  
 B conc. 0%(0-12min), 0→100%(12→20min), 100%(20-22min)

Flow rate: 1.5 ml/min  
 Temperature: 25°C  
 Detection: UV280nm

Sample: 1; Similar compound 1  
 2; Vancomycin Hydrochloride (5.0µg)  
 3; Similar compound 2



AP-0901


1) Drugs

• Bisoprolol Fumarate

**COSMOSIL Application Data**

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 30mmol/l KH<sub>2</sub>PO<sub>4</sub> (pH2.5 with H<sub>3</sub>PO<sub>4</sub>) = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV225nm

Sample: 1; Fumaric Acid  
 2; Bisoprolol




NACALAI TESQUE, INC  
 AP-1090

• Hydrochlorothiazide

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 100mmol/l NaH<sub>2</sub>PO<sub>4</sub> (pH3.0 with H<sub>3</sub>PO<sub>4</sub>) = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; 4-Amino-6-chlorobenzene-1,3-disulfonamide (1.0µg)  
 2; Hydrochlorothiazide (3.0µg)  
 3; 4-Aminoacetophenone (3.6µg)



NACALAI TESQUE, INC  
 AP-0683


1) Drugs

• Piperacillin Sodium

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 50mmol/l Acetic Acid, 25mmol/l Triethylamine = 21/79  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Ampicillin (10µg)  
 2; Acetanilide (0.75µg)  
 3; Piperacillin Sodium (2.5µg)




NACALAI TESQUE, INC  
 AP-0799

• Piperacillin Sodium

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A: Acetonitrile/H<sub>2</sub>O/ 200mmol/l KH<sub>2</sub>PO<sub>4</sub> = 4/45/1  
 B: Acetonitrile/H<sub>2</sub>O/ 200mmol/l KH<sub>2</sub>PO<sub>4</sub> = 25/24/1  
 B conc. 0%(0-7 min), 0→17%(7→13 min), 17%(13-25 min), 17→80%(25→40min), 80%(40-44 min)  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm

Sample: 1; Piperacillin Sodium (0.4µg)




NACALAI TESQUE, INC  
 AP-0802

• Hydrocortisone

**COSMOSIL Application Data**

Column: 5SL-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Chloroform/Methanol/Acetic Acid = 1000/20/1  
 Flow rate: 2.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Prednisone (0.9µg)  
 2; Hydrocortisone (2.0µg)




NACALAI TESQUE, INC  
 AP-0685

• Hydrocortisone Succinate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 40mmol/l CH<sub>3</sub>COONa (pH4.0 with Acetic Acid) = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Hydrocortisone (0.4µg)  
 2; Hydrocortisone Succinate (1.0µg)  
 3; Butyl p-Hydroxybenzoate (0.3µg)




NACALAI TESQUE, INC  
 AP-0689

• Pimaricin

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: 0.1%Ammonium Acetate-Methanol/Tetrahydrofuran/H<sub>2</sub>O = 44/2/47  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV303nm

Sample: 1; Pimaricin (0.02µg)




NACALAI TESQUE, INC  
 AP-0796

• Vincristine Sulfate

**COSMOSIL Application Data**

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 1.67% Diethylamine (pH7.5 with H<sub>3</sub>PO<sub>4</sub>) = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV297nm

Sample: 1; Vincristine Sulfate (10µg)




NACALAI TESQUE, INC  
 AP-1128

• Hydrocortisone Acetate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 45/55  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Hydrocortisone (2.0µg)  
 2; Hydrocortisone Acetate (4.0µg)  
 3; Benzyl p-Hydroxybenzoate (2.0µg)




NACALAI TESQUE, INC  
 AP-0686

• Hypromellose Phthalate (impurity)

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 100mmol/l Cyanoacetic Acid = 15/85  
 Flow rate: 2.0 ml/min  
 Temperature: 30°C  
 Detection: UV235nm

Sample: 1; Phthalic Acid (0.5µg)




NACALAI TESQUE, INC  
 AP-0933

• Famotidine

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/Methanol/ 10mmol/l Sodium l-Heptanesulfonate (pH3.0 with Acetic Acid) = 19/3/78  
 Flow rate: 0.5 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Famotidine (0.5µg)  
 2; Methyl p-Hydroxybenzoate (0.8µg)




NACALAI TESQUE, INC  
 AP-0651

• Scopolamine Butyl Bromide

**COSMOSIL Application Data**

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 19mmol/l Sodium Lauryl Sulfate = 68/37(pH3.6 with H<sub>3</sub>PO<sub>4</sub>)  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm

Sample: 1; Scopolamine Hydrobromide (2.0µg)  
 2; Scopolamine n-Butyl Bromide(2.0µg)




NACALAI TESQUE, INC  
 AP-0852

1) Drugs

• Bufexamac

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/Acetonitrile/  
 13.6mmol/l Sodium *I*-Octane Sulfonate  
 0.94%Acetic Acid, 1.9mmol/l EDTA  
 = 24/24/52  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV275nm  
 Sample: 1; Bufexamac (2μg)  
 2; Diphenylimidazole (0.16μg)




AP-0522

• Pravastatin

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/H<sub>2</sub>O/Acetic Acid/Triethylamine  
 = 500/500/1/1  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV238nm  
 Sample: 1; Ethyl *p*-Hydroxybenzoate (0.7μg)  
 2; Pravastatin Sodium (1.6μg)




AP-0811

• Flavin Adenine Dinucleotide Sodium

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 0.2%KH<sub>2</sub>PO<sub>4</sub> = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 35°C  
 Detection: UV260nm  
 Sample: 1; Adenosine (1.0μg)  
 2; FAD (4.0μg)  
 3; FMN (4.0μg)  
 4; Riboflavin (1.0μg)




AP-0653

• Fluoxymesterone

**COSMOSIL Application Data**

Column: 5SL-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: *n*-Butyl Chloride/Water-saturated *n*-Butyl  
 Chloride/Tetrahydrofuran/Methanol/  
 Acetic Acid = 95/95/14/7/6  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm  
 Sample: 1; Fluoxymesterone (2.7μg)  
 2; Methyl Prednisolone (2.0μg)




AP-0669

• Fluocinonide

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV254nm  
 Sample: 1; Fluocinonide (4.0μg)  
 2; Propyl Benzoate (12.6μg)




AP-0663

• Fluocinolone Acetonide

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV254nm  
 Sample: 1; Ethyl *p*-Hydroxybenzoate (0.6μg)  
 2; Isopropyl *p*-Hydroxybenzoate (1.0μg)  
 3; Fluocinolone Acetonide (4.0μg)  
 4; Propyl *p*-Hydroxybenzoate (1.0μg)




AP-0660

1) Drugs

• Fluocinolone Acetonide

**COSMOSIL Application Data**

Column: 5SL-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Water-saturated Chloroform/Methanol/  
 Acetic Acid = 200/3/2  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm  
 Sample: 1; Triamcinolone Acetonide (3.0μg)  
 2; Fluocinolone Acetonide (3.0μg)




AP-0662

• Fluorometholone

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 35°C  
 Detection: UV254nm  
 Sample: 1; Fluorometholone (0.20μg)  
 2; Butyl *p*-Hydroxybenzoate (0.16μg)




AP-0666

• Flurbiprofen

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 7.7%H<sub>3</sub>PO<sub>4</sub> = 35/65  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm  
 Sample: 1; Flurbiprofen (0.2μg)




AP-0671

• Prednisolone Acetate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm  
 Sample: 1; Prednisolone (1.5μg)  
 2; Prednisolone Acetate (1.0μg)  
 3; Cortisone Acetate (1.5μg)  
 4; Butyl *p*-Hydroxybenzoate (0.5μg)




AP-0814

• Procaine Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 0.1% Sodium *I*-Pentanesulfonate,  
 50mmol/l KH<sub>2</sub>PO<sub>4</sub>(pH3.0 with H<sub>3</sub>PO<sub>4</sub>)  
 = 20/80  
 Flow rate: 0.5 ml/min  
 Temperature: 40°C  
 Detection: UV254nm  
 Sample: 1; Procaine Hydrochloride (1.25μg)  
 2; Caffeine (1.25μg)




AP-0817

• Procatamol Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 5mmol/l Sodium *I*-Pentanesulfonate/  
 Acetic Acid = 23/76/1  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV254nm  
 Sample: 1; Procatamol Hydrochloride (0.06μg)



AP-0821


1) Drugs

• Propranolol Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 12mmol/l Sodium Lauryl Sulfate, 0.7mmol/l Tetra-*n*-butylammonium Phosphate/ Sulfuric Acid = 550/450/1(pH3.3 with NaOH)  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV292nm

Sample:  
 1; Propranolol hydrochloride (0.08µg)




NACALAI TESQUE, INC  
 AP-0823

• Flopropione

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 1.2% H<sub>3</sub>PO<sub>4</sub> = 57/43  
 Flow rate: 1.0 ml/min  
 Temperature: 35°C  
 Detection: UV267nm

Sample: 1; Flopropione (0.40µg)  
 2; Ethyl 4-Hydroxybenzoate (0.26µg)




NACALAI TESQUE, INC  
 AP-0657

• Bromhexine Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 7.3mmol/l KH<sub>2</sub>PO<sub>4</sub> (pH7.0 with NaOH) = 80/20  
 Flow rate: 2.0 ml/min  
 Temperature: 40°C  
 Detection: UV245nm

Sample: 1; Bamethane Sulfate (25µg)  
 2; Bromhexine Hydrochloride (1.25µg)




NACALAI TESQUE, INC  
 AP-0518

• Beclometasone Dipropionate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Beclometasone Dipropionate (1.6µg)  
 2; Testosterone Propionate (0.8µg)




NACALAI TESQUE, INC  
 AP-0479

• Bezafibrate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 0.1% Acetic Acid = 9/4  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV230nm

Sample: 1; *p*-Chlorobenzoic Acid (0.55µg)  
 2; Bezafibrate (1.05µg)




NACALAI TESQUE, INC  
 AP-0510

• Bezafibrate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 0.1% Acetic Acid = 9/4  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV230nm

Sample: 1; *p*-Nitrophenol (0.32µg)  
 2; Bezafibrate (0.44µg)



NACALAI TESQUE, INC  
 AP-0513


1) Drugs

• Betahistine Mesilate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: 8mmol/l Sodium Lauryl Sulfate- Acetonitrile/ 2% Acetic Acid, 0.5% Diethylamine = 37/63  
 Flow rate: 1.0 ml/min  
 Temperature: 35°C  
 Detection: UV261nm

Sample:  
 1; 2-Vinylpyridine (0.16µg)  
 2; Betahistine Mesilate (0.16µg)




NACALAI TESQUE, INC  
 AP-0488

• Betamethasone

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV241nm

Sample: 1; Betamethasone (0.06µg)




NACALAI TESQUE, INC  
 AP-0497

• Betamethasone Sodium Phosphate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 5mmol/l Tetra-*n*-butylammonium Bromide 8.9mmol/l Na<sub>2</sub>HPO<sub>4</sub> 50.7mmol/l KH<sub>2</sub>PO<sub>4</sub> = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Betamethasone Sodium Phosphate (1.0µg)  
 2; Butyl *p*-Hydroxybenzoate (0.2µg)




NACALAI TESQUE, INC  
 AP-0494

• Betamethasone Valerate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 65/35  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Betamethasone Valerate (1.5µg)  
 2; Beclometasone Dipropionate (1.5µg)




NACALAI TESQUE, INC  
 AP-0503

• Berberine Chloride Hydrate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 6mmol/l Sodium Lauryl Sulfate 25mmol/l KH<sub>2</sub>PO<sub>4</sub> = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV345nm

Sample: 1; Palmatine Chloride (1.0µg)  
 2; Berberine Chloride (1.0µg)




NACALAI TESQUE, INC  
 AP-0487

• Benzylpenicillin Potassium

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 50mmol/l (NH<sub>4</sub>)<sub>2</sub>HPO<sub>4</sub> = 24/76 (pH8 with H<sub>3</sub>PO<sub>4</sub>)  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample:  
 1; Benzylpenicillin Potassium (4.0µg)  
 2; Methyl *p*-Hydroxybenzoate (0.06µg)




NACALAI TESQUE, INC  
 AP-0482

## 1) Drugs

### • Calcium Folate

#### **COSMOSIL Application Data**


Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-250mm  
Mobile phase: [H<sub>2</sub>O/Acetonitrile/Methanol/  
40%Tetra-*n*-butylammonium Hydroxide  
= 760/200/8.6/9.4 (pH7.5 with NaH<sub>2</sub>PO<sub>4</sub>)  
→1000(with H<sub>2</sub>O)]  
Flow rate: 0.5 ml/min  
Temperature: 30°C  
Detection: UV254nm  
Sample: 1; Calcium Folate (2.8μg)  
2; Folic Acid (0.8μg)



### • Calcium Folate

#### **COSMOSIL Application Data**


Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Methanol/ 8mmol/l Na<sub>2</sub>HPO<sub>4</sub>/ Tetrabutylammonium  
Hydroxide = 110/385/4 (pH 7.5 with H<sub>3</sub>PO<sub>4</sub>)  
Flow rate: 1.0 ml/min  
Temperature: 45°C  
Detection: UV254nm  
Sample: 1; Calcium Folate (2.0μg)  
2; Folic Acid (2.0μg)



### • Mitomycin C

#### **COSMOSIL Application Data**


Column: 5PE-MS  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Methanol/ 20mmol/l Ammonium Acetate,  
0.025%Acetic Acid = 40/60  
Flow rate: 1.0 ml/min  
Temperature: 25°C  
Detection: UV365nm  
Sample: 1; Mitomycin C (5.0μg)  
2; Ethyl Vanillin (75μg)



### • Mizoribine

#### **COSMOSIL Application Data**


Column: 5C<sub>18</sub>-PAQ  
Column size: 4.6mmI.D.-250mm  
Mobile phase: 0.067% H<sub>3</sub>PO<sub>4</sub> aq.  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV279nm  
Sample: 1; Mizoribine (1.0μg)



### • Minocycline Hydrochloride

#### **COSMOSIL Application Data**


Column: 5C<sub>8</sub>-MS  
Column size: 4.6mmI.D.-150mm  
Mobile phase: 2.8%Ammonium Oxalate/*N,N*-dimethylformamide/  
100mmol/l EDTA = 11/5/4  
(pH6.2 with Tetrabutylammonium Hydroxide)  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV280nm  
Sample: 1; 4-*epi*-Minocycline  
2; Minocycline Hydrochloride



### • Mexiletine Hydrochloride

#### **COSMOSIL Application Data**

Column: 5C<sub>8</sub>-MS  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ 14.4mmol/l Sodium Lauryl Sulfate,  
25mmol/l NaH<sub>2</sub>PO<sub>4</sub> = 21/30  
Flow rate: 2.0 ml/min  
Temperature: 30°C  
Detection: UV210nm  
Sample: 1; Phenetylamine Hydrochloride (0.6μg)  
2; Mexiletine Hydrochloride (1.0μg)




## 1) Drugs

### • Mecobalamin

#### **COSMOSIL Application Data**


Column: 5C<sub>18</sub>-PAQ  
Column size: 4.6mmI.D.-250mm  
Mobile phase: 100mmol/l Sodium *L*-Hexanesulfonate,  
Acetonitrile/ 20mmol/l Phosphate  
buffer(pH3.5) = 20/80  
Flow rate: 0.5 ml/min  
Temperature: 40°C  
Detection: UV266nm  
Sample: 1; Cyanocobalamin (0.5μg)  
2; Hydroxocobalamin Acetate (0.5μg)  
3; Mecobalamin (1.0μg)



### • Meticrane

#### **COSMOSIL Application Data**


Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ H<sub>2</sub>O = 15/85  
Flow rate: 1.0 ml/min  
Temperature: 40°C  
Detection: UV230nm  
Sample: 1; Caffeine (0.2μg)  
2; Meticrane (0.2μg)



### • Meticrane

#### **COSMOSIL Application Data**


Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ H<sub>2</sub>O = 50/50  
Flow rate: 1.0 ml/min  
Temperature: 40°C  
Detection: UV230nm  
Sample: 1; Meticrane (0.4μg)  
2; Methyl *p*-Hydroxybenzoate (0.4μg)



### • Methotrexate

#### **COSMOSIL Application Data**


Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 200mmol/l Na<sub>2</sub>HPO<sub>4</sub>  
(pH6.0 with 100mmol/l Citric Acid) = 10/90  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV302nm  
Sample: 1; Folic Acid (1.0μg)  
2; Methotrexate (1.0μg)



### • Metformin Hydrochloride

#### **COSMOSIL Application Data**


Column: 5C<sub>18</sub>-AR-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ 4.5mmol/l Sodium Lauryl Sulfate,  
0.04%H<sub>2</sub>PO<sub>4</sub> = 38/62  
Flow rate: 1.0 ml/min  
Temperature: 40°C  
Detection: UV235nm  
Sample: 1; Metformin Hydrochloride (0.50μg)  
2; Isobutyl *p*-Hydroxybenzoate (0.95μg)



### • Metronidazole

#### **COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Methanol/ H<sub>2</sub>O = 20/80  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV320nm  
Sample: 1; Metronidazole (0.24μg)




1) Drugs

• Menatetrenone

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Methanol  
Flow rate: 1.0 ml/min  
Temperature: 40°C  
Detection: UV270nm

Sample: 1; Menatetrenone (0.6µg)  
2; Phytanadione (0.6µg)




NACALAI TESQUE, INC  
AP-0742

• Mepivacaine Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: 10mmol/l Sodium Lauryl Sulfate-Acetonitrile/  
20mmol/l Phosphate buffer(pH3.0) = 45/55  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254nm

Sample: 1; Mepivacaine Hydrochloride (2.0µg)  
2; Benzophenone (0.04µg)



NACALAI TESQUE, INC  
AP-0746


1) Drugs

• Latamoxef Sodiums

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Methanol/ 100mmol/l Ammonium  
Acetate = 5/95  
Flow rate: 1.0 ml/min  
Temperature: 25°C  
Detection: UV254nm

Sample: 1; Latamoxef Sodium (isomer 1)  
2; Latamoxef Sodium (isomer 2)




NACALAI TESQUE, INC  
AP-0716

• Liothyronine Sodium

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Methanol/ H<sub>2</sub>O = 60/40  
Flow rate: 1.0 ml/min  
Temperature: 25°C  
Detection: UV220nm

Sample: 1; Propyl *p*-Hydroxybenzoate (0.14µg)  
2; Liothyronine Sodium (0.08µg)




NACALAI TESQUE, INC  
AP-0728

• Meropenem

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Methanol/ 0.1%Triethylamine  
(pH5.0 with H<sub>3</sub>PO<sub>4</sub>) = 10/90  
Flow rate: 1.0 ml/min  
Temperature: 25°C  
Detection: UV220nm

Sample: 1; Meropenem Trihydrate (2.5µg)  
2; Benzyl Alcohol (1.5µg)




NACALAI TESQUE, INC  
AP-0749

• Ubidecarenone

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Methanol/ Ethanol = 50/50  
Flow rate: 0.5 ml/min  
Temperature: 35°C  
Detection: UV275nm

Sample: 1; Ubidecarenone (2.5µg)




NACALAI TESQUE, INC  
AP-0900

• Lisinopril

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
Column size: 4.6mmI.D.-150mm  
Mobile phase: A; 25mmol/l NaH<sub>2</sub>PO<sub>4</sub>  
B; Acetonitrile/ 25mmol/l NaH<sub>2</sub>PO<sub>4</sub> = 40/60  
B conc. 10→50%(0→10min), 50%(10-25min)  
Flow rate: 1.5 ml/min  
Temperature: 60°C  
Detection: UV215nm

Sample: 1; Lisinopril (0.45µg)  
2; Caffeine (0.015µg)




NACALAI TESQUE, INC  
AP-0732

• Lisinopril

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ 25mmol/l NaH<sub>2</sub>PO<sub>4</sub>  
= 5/95  
Flow rate: 1.0 ml/min  
Temperature: 60°C  
Detection: UV215nm

Sample: 1; Lisinopril (2.0µg)  
2; Caffeine (0.5µg)




NACALAI TESQUE, INC  
AP-0733

• Iodine, Salicylic Acid and Phenol

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Methanol/ 100mmol/l Phosphate  
buffer(pH7.0) = 25/75  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV270nm

Sample: 1; Benzoic Acid (2.0µg)  
2; Theophylline (0.5µg)  
3; Salicylic Acid (2.0µg)  
4; Phenol (2.0µg)




NACALAI TESQUE, INC  
AP-0935

• Latamoxef Sodiums

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Methanol/ 51mmol/l KH<sub>2</sub>PO<sub>4</sub>,  
9mmol/l Na<sub>2</sub>HPO<sub>4</sub>,  
5mmol/l Tetra-*n*-Butylammonium Bromide  
= 25/75  
Flow rate: 1.0 ml/min  
Temperature: 25°C  
Detection: UV254nm

Sample: 1; 5-Mercapto-*l*-methyltetrazole (1.65µg)  
2; Latamoxef Sodium (isomer 1)  
3; Latamoxef Sodium (isomer 2)  
4; *m*-Cresol (7.5µg)




NACALAI TESQUE, INC  
AP-0711

• Lidocaine

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: 10mmol/l Sodium Lauryl Sulfate-Acetonitrile/  
20mmol/l Phosphate buffer(pH3.0) =45/55  
Flow rate: 1.5 ml/min  
Temperature: 25°C  
Detection: UV254nm

Sample: 1; Lidocaine (8.5µg)  
2; Benzophenone (0.25µg)




NACALAI TESQUE, INC  
AP-0724

• Ritodrine Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Methanol/ 71mmol/l (NH<sub>4</sub>)<sub>2</sub>HPO<sub>4</sub>,  
7.8mmol/l Sodium *l*-Heptanesulfonate  
= 30/70(pH3.0 with H<sub>3</sub>PO<sub>4</sub>)  
Flow rate: 1.0 ml/min  
Temperature: 25°C  
Detection: UV220nm

Sample: 1; Ritodrine Hydrochloride  
2; *threo*-Ritodrine Hydrochloride  
3; by-product



NACALAI TESQUE, INC  
AP-0835


## 1) Drugs

### • Ritodrine Hydrochloride

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 71mmol/l (NH<sub>4</sub>)<sub>2</sub>HPO<sub>4</sub>,  
 7.8mmol/l Sodium *J*-Heptanesulfonate  
 = 30/70 (pH 3.0 with H<sub>3</sub>PO<sub>4</sub>)  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV274nm

Sample:  
 1; Ritodrine Hydrochloride (0.30µg)  
 2; Methyl *p*-Hydroxybenzoate (0.05µg)




NACALAI TESQUE, INC  
 AP-0838

### • Rifampicin

#### COSMOSIL Application Data

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: 20mmol/l Citric Acid, 11mmol/l Sodium Perchlorate-  
 Acetonitrile/H<sub>2</sub>O/ 1mol/l KH<sub>2</sub>PO<sub>4</sub>,  
 55mmol/l H<sub>3</sub>PO<sub>4</sub> (pH 3.1) = 7/11/2  
 Flow rate: 2.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample:  
 1; Butyl *p*-Hydroxybenzoate (0.15µg)  
 2; Rifampicin (0.80µg)




NACALAI TESQUE, INC  
 AP-0832

### • Rifampicin

#### COSMOSIL Application Data

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 15mmol/l Sodium Perchlorate,  
 28mmol/l Citric Acid, 17mmol/l NaH<sub>2</sub>PO<sub>4</sub> = 45/55  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Rifampicin (0.08µg)




NACALAI TESQUE, INC  
 AP-0833

### • Lincomycin Hydrochloride

#### COSMOSIL Application Data

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/Methanol/ 1.35% H<sub>3</sub>PO<sub>4</sub>  
 (pH 6.0 with Ammonia) = 15/15/78  
 Flow rate: 1.0 ml/min  
 Temperature: 46°C  
 Detection: UV210nm

Sample: 1; Lincomycin Hydrochloride (20µg)




NACALAI TESQUE, INC  
 AP-0729

### • Reserpine

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 50mmol/l KH<sub>2</sub>PO<sub>4</sub>  
 (pH 3.0 with H<sub>3</sub>PO<sub>4</sub>) = 45/55  
 Flow rate: 0.5 ml/min  
 Temperature: 40°C  
 Detection: UV268nm

Sample:  
 1; Reserpine (0.20µg)  
 2; Butyl *p*-Hydroxybenzoate (0.08µg)




NACALAI TESQUE, INC  
 AP-0830

### • Levallorphan Tartrate

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 7mmol/l Sodium Lauryl Sulfate,  
 0.1% H<sub>3</sub>PO<sub>4</sub> (pH 3.0 with NaOH) = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV280nm

Sample:  
 1; *iso*-Butyl *p*-Hydroxybenzoate (0.3µg)  
 2; Levallorphan Tartrate (1.7µg)



NACALAI TESQUE, INC  
 AP-0718


## 1) Drugs

### • Levothyroxine

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Methanol/ 0.15% H<sub>3</sub>PO<sub>4</sub> = 67/33  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV230nm

Sample:  
 1; Levothyroxine Sodium (0.1µg)  
 2; Ethinylestradiol (0.2µg)




NACALAI TESQUE, INC  
 AP-0721

### • Roxatidine Acetate Hydrochloride

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O/ Triethylamine/  
 Acetic Acid = 60/340/2/1  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV274nm

Sample: 1; Roxatidine Acetate Hydrochloride (4.2µg)




NACALAI TESQUE, INC  
 AP-0840

### • Roxithromycin

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 0.2mmol/l (NH<sub>4</sub>)<sub>2</sub>HPO<sub>4</sub>  
 (pH 5.3 with NaOH) = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV205nm

Sample: 1; Roxithromycin (1.0µg)




NACALAI TESQUE, INC  
 AP-0845

### • Tricyclic Drugs

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate  
 buffer (pH 7) = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.2AUFS

Sample:  
 1; Carbamazepine  
 2; Desipramine Hydrochloride  
 3; Imipramine Hydrochloride




NACALAI TESQUE, INC  
 AP-0066

### • Bronchodilators

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 1.0AUFS

Sample:  
 1; Uracil (0.5µg)  
 2; Theobromine (1.5µg)  
 3; Theophylline (2.0µg)  
 4; Caffeine (2.0µg)  
 5; Phenol (0.8µg)




NACALAI TESQUE, INC  
 AP-0076

### • Antiarrhythmic Drugs

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate  
 buffer (pH 7) = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.05AUFS

Sample:  
 1; Phenytoin  
 2; Ketamine Hydrochloride  
 3; Quinidine  
 4; Lidocaine



NACALAI TESQUE, INC  
 AP-0067


## 1) Drugs

### ● Antiarrhythmic Drugs

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH2) = 10/90  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.2AUFS

Sample: 1; Procainamide (0.5μg)  
 2; N-Acetylprocainamide (0.5μg)




NACALAI TESQUE, INC  
 AP-0069

### ● Antiepileptics

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH7) = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.5AUFS

Sample: 1; Barbital (7.59μg)  
 2; Phenytoin (10.38μg)  
 3; Carbamazepine (1.02μg)



NACALAI TESQUE, INC  
 AP-0077


## 1) Drugs

### ● Antihyperlipidemic Drugs

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphoric Acid = 10/90  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.5AUFS

Sample: 1; Hydralazine Hydrochloride (1.0μg)  
 2; Todalazine Hydrochloride (1.0μg)




NACALAI TESQUE, INC  
 AP-0080

### ● Profens

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 20mmol/l Acetic Acid = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 1.0AUFS

Sample: 1; Ketoprofen (1.60μg)  
 2; Ibuprofen (1.69μg)  
 3; Flurbiprofen (1.57μg)




NACALAI TESQUE, INC  
 AP-0081

### ● Analgesics

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphoric Acid = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.5AUFS

Sample: 1; Aminoantipyrine (0.74μg)  
 2; Acetaminophen (2.97μg)  
 3; Antipyrine (0.85μg)  
 4; Acetylsalicylic Acid (Aspirine) (0.54μg)  
 5; Phenacetin (1.0μg)  
 6; Salicylic acid (1.0μg)




NACALAI TESQUE, INC  
 AP-0073

### ● Analgesics

#### COSMOSIL Application Data

Column: 5SL-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Ethyl Acetate/Hexane = 1/1  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.2AUFS

Sample: 1; Acetanilide (1.0μg)  
 2; Phenacetin (1.0μg)  
 3; Acetaminophen (1.0μg)




NACALAI TESQUE, INC  
 AP-0074

### ● Cardiac Glycosides

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV230nm, 0.32AUFS

Sample: 1; Digitoxigenin (2.5μg)  
 2; Digitoxin (5.0μg)




NACALAI TESQUE, INC  
 AP-0083

### ● Anticancer Drugs

#### COSMOSIL Application Data

Column: Sugar-D  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/H<sub>2</sub>O = 80/20  
 Flow rate: 0.3 ml/min  
 Temperature: Room temperature  
 Detection: UV226nm

Sample: 1; cis-Platin (CDDP) (1.46μg)  
 2; Guanosine (0.50μg)




NACALAI TESQUE, INC  
 AP-0380

### ● Histamine H1-Receptor Blockers

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 10mmol/l SDS, 0.1% Phosphoric Acid = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Uracil  
 2; Diphenhydramine Hydrochloride




NACALAI TESQUE, INC  
 AP-0078

### ● Histamine H2-Receptor Blockers

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH7) = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.2AUFS

Sample: 1; Famotidine (1.3μg)  
 2; Cimetidine (33.5μg)  
 3; Ranitidine (1.6μg)




NACALAI TESQUE, INC  
 AP-0079

### ● Tetracyclines Antibiotics

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH3) = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 20°C  
 Detection: UV254nm

Sample: 1; Oxytetracycline  
 2; Tetracycline  
 3; Chlortetracycline




NACALAI TESQUE, INC  
 AP-0085

### ● Penicillin Antibiotics

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH7) = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV235nm, 0.2AUFS

Sample: 1; Carbenicillin (1.5μg)  
 2; Ampicillin (3.0μg)  
 3; Methicillin Sodium Salt (1.5μg)  
 4; Penicillin G Potassium Salt (3.0μg)



NACALAI TESQUE, INC  
 AP-0086


## 1) Drugs

### • Quinolone Antimicrobials

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH3) = 55/45  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV270nm, 0.16AUFS

Sample: 1; Ofloxacin (0.39µg)  
 2; Oxolinic Acid (0.08µg)  
 3; Flumequine (1.08µg)  
 4; Piromidic Acid (0.08µg)




NACALAI TESQUE, INC  
 AP-0087

### • Nitrofurantoin Antimicrobials

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol / 20mmol/l Phosphoric Acid = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV260nm, 0.5AUFS

Sample: 1; Nitrofurantoin (1.0µg)  
 2; Nitrofurazone (1.0µg)




NACALAI TESQUE, INC  
 AP-0089

### • Streptomycin Sulfate

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 5mmol/l Sodium l-Hexanesulfonate, 20mmol/l KH<sub>2</sub>PO<sub>4</sub> = 10/90  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV205nm, 0.2AUFS

Sample: Streptomycin Sulfate (5.0µg)




NACALAI TESQUE, INC  
 AP-0091

### • Sulfa Drugs

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol / 20mmol/l Phosphoric Acid = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV230nm, 0.5AUFS

Sample: 1; Sulfaisomidin (0.24µg)  
 2; Sulfathiazole (0.24µg)  
 3; Sulfamethazine (0.24µg)  
 4; Sulfamethoxy pyridazine (0.24µg)  
 5; Sulfamethoxazole (0.24µg)  
 6; Sulfachloropyridazine (0.24µg)  
 7; Sulfamonomethoxine (0.24µg)  
 8; Sulfadimethoxine (0.24µg)




NACALAI TESQUE, INC  
 AP-0090

### • Allantoin

#### COSMOSIL Application Data

Column: HILIC  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7.0) = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm

Sample: 1; Allantoin  
 2; Allantoic Acid



NACALAI TESQUE, INC  
 AP-1060


## 2) Crude Drugs

### • Bearberry Leaf

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 1mmol/l HCl = 5/95  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm

Sample: 1; Arbutin (5.0µg)  
 2; Hydroquinone (5.0µg)  
 3; Gallic Acid (5.0µg)




NACALAI TESQUE, INC  
 AP-0906

### • Scutellaria Root

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 0.68% Phosphoric Acid = 28/72  
 Flow rate: 1.0 ml/min  
 Temperature: 50°C  
 Detection: UV277nm

Sample: 1; Baicalin (0.1µg)  
 2; Methyl p-Hydroxybenzoate (0.2µg)




NACALAI TESQUE, INC  
 AP-1134

### • Phellodendron Bark

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 5.9mmol/l Sodium Lauryl Sulfate, 25mmol/l KH<sub>2</sub>PO<sub>4</sub> = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV345nm

Sample: 1; Palmatine Chloride Hydrate (1.0µg)  
 2; Berberine Chloride (1.0µg)




NACALAI TESQUE, INC  
 AP-0964

### • Pueraria Root

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 50mmol/l NaH<sub>2</sub>PO<sub>4</sub> = 10/90  
 Flow rate: 0.5 ml/min  
 Temperature: 40°C  
 Detection: UV250nm

Sample: 1; Puerarin (1.3µg)




NACALAI TESQUE, INC  
 AP-0969

### • Kamishoyosan Extract

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 0.11% H<sub>3</sub>PO<sub>4</sub> = 10/90  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV240nm

Sample: 1; Geniposide (1.2µg)




NACALAI TESQUE, INC  
 AP-0942

### • Glycyrrhiza

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 2% H<sub>3</sub>PO<sub>4</sub> = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Glycyrrhizic Acid (5.0µg)  
 2; Propyl p-Hydroxybenzoate (1.0µg)



NACALAI TESQUE, INC  
 AP-0927


## 2) Crude Drugs

### • Dried Yeast

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ 9.2mmol/l Sodium *l*-Octanesulfonate, 20mmol/l KH<sub>2</sub>PO<sub>4</sub>(pH3.5 with H<sub>3</sub>PO<sub>4</sub>) = 20/80  
Flow rate: 1.0 ml/min  
Temperature: 40°C  
Detection: UV254nm

Sample: 1; Thiamine Hydrochloride (0.10µg)  
2; Phenacetin (0.06µg)




NACALAI TESQUE, INC  
AP-0919

### • Keishibukuryougan Extract

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Methanol/ 50mmol/l NaH<sub>2</sub>PO<sub>4</sub> = 1/5  
Flow rate: 0.8 ml/min  
Temperature: 45°C  
Detection: UV210nm

Sample: 1; *D*-(-)-Amygdalin (2.0µg)



NACALAI TESQUE, INC  
AP-0944


## 2) Crude Drugs

### • Saireito Extract

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ 50mmol/l NaH<sub>2</sub>PO<sub>4</sub> = 3/5  
Flow rate: 1.0 ml/min  
Temperature: 40°C  
Detection: UV254nm

Sample: 1; Saikosaponin b2 (1.0µg)



NACALAI TESQUE, INC  
AP-0978

### • Saireito Extract

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ 0.5% H<sub>3</sub>PO<sub>4</sub> = 24/76  
Flow rate: 1.0 ml/min  
Temperature: 40°C  
Detection: UV277nm

Sample: 1; Baicalin (0.5µg)




NACALAI TESQUE, INC  
AP-0981

### • Red Ginseng and Ginseng

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ H<sub>2</sub>O = 20/80  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV203nm

Sample: 1; Ginsenoside Rg1 (2.5µg)  
2; Ginsenoside Re (2.5µg)




NACALAI TESQUE, INC  
AP-0922

### • Red Ginseng and Ginseng

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ H<sub>2</sub>O = 30/70  
Flow rate: 1.0 ml/min  
Temperature: 40°C  
Detection: UV203nm

Sample: 1; Ginsenoside Rb1 (2.5µg)  
2; Ginsenoside Rc (2.5µg)




NACALAI TESQUE, INC  
AP-0923

### • Powdered Gardenia Fruit

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ H<sub>2</sub>O = 60/40  
Flow rate: 0.5 ml/min  
Temperature: 30°C  
Detection: UV240nm

Sample: 1; Caffeine (0.175µg)  
2; Geniposide (0.175µg)



NACALAI TESQUE, INC  
AP-0966

### • Peony Root

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ 0.12% H<sub>3</sub>PO<sub>4</sub> = 15/85  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV232nm

Sample: 1; Albiflorin (1.0µg)  
2; Paeoniflorin (1.0µg)




NACALAI TESQUE, INC  
AP-0960

### • Magnolia Bark

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ 2% Acetic Acid = 50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV289nm

Sample: 1; Honokiol (1.0µg)  
2; Magnolol (1.0µg)




NACALAI TESQUE, INC  
AP-0948

### • Bupleurum Root

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ H<sub>2</sub>O = 40/60  
Flow rate: 1.0 ml/min  
Temperature: 50°C  
Detection: UV206nm

Sample: 1; Saikosaponin a (1.0µg)  
2; Saikosaponin d (1.0µg)




NACALAI TESQUE, INC  
AP-0910

### • Toad Venom

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ 0.1% H<sub>3</sub>PO<sub>4</sub> = 45/55  
Flow rate: 1.0 ml/min  
Temperature: 40°C  
Detection: UV300nm

Sample: 1; Bufalin (0.4µg)  
2; Cinobufagin (0.4µg)  
3; Resibufogenin (0.6µg)  
4; Indometacin(Indomethacin) (0.6µg)




NACALAI TESQUE, INC  
AP-1000

### • Senna Leaf

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-250mm  
Mobile phase: 5mmol/l Tetra-*n*-heptylammonium Bromide-Acetonitrile/ 0.1mol/l Acetic Acid, 1% Sodium Acetate(pH5.0) = 32/68  
Flow rate: 1.0 ml/min  
Temperature: 50°C  
Detection: UV340nm

Sample: 1; Sennoside B (1.0µg)  
2; Sennoside A (0.5µg)



NACALAI TESQUE, INC  
AP-0993


## 2) Crude Drugs

### • Swertia Herb

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ H<sub>2</sub>O = 10/90  
Flow rate: 0.5 ml/min  
Temperature: 50°C  
Detection: UV238nm

Sample: 1; Theophylline (1.0µg)  
2; Swertiamarin (1.0µg)




AP-0996

### • Rhubarb

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ 1.25%Acetic Acid = 20/80  
Flow rate: 0.5 ml/min  
Temperature: 40°C  
Detection: UV340nm

Sample: 1; Senoside A (2.0µg)  
2; Naringin (2.0µg)



AP-0972


## 2) Crude Drugs

### • Moutan Bark

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 3.1%Acetic Acid = 35/65  
Flow rate: 1.0 ml/min  
Temperature: 40°C  
Detection: UV274nm

Sample: 1; Paeonol (0.5µg)  
2; Butyl *p*-Hydroxybenzoate (2.2µg)




AP-0950

### • Hochuekkito Extract

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ 1.2%*H*<sub>3</sub>*P*O<sub>4</sub> = 18/82  
Flow rate: 1.0 ml/min  
Temperature: 40°C  
Detection: UV285nm

Sample: 1; Naringin (0.1µg)  
2; Hesperidin (0.1µg)




AP-0930

### • Daiokanzoto Extract

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ 0.04%*H*<sub>3</sub>*P*O<sub>4</sub> = 540/2460  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV340nm

Sample: 1; Senoside A (0.26µg)




AP-0916

### • Capsicum

#### COSMOSIL Application Data

Column: 5PE-MS  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 0.1%*H*<sub>3</sub>*P*O<sub>4</sub> = 40/60  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV281nm

Sample: 1; Vanillynonamide (0.4µg)  
2; Capsaicin (0.4µg)




AP-0913

### • Povidone

#### COSMOSIL Application Data

Column: 5C<sub>8</sub>-MS  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Methanol/ H<sub>2</sub>O = 20/80  
Flow rate: 1.0 ml/min  
Temperature: 40°C  
Detection: UV254nm

Sample: 1; *l*-Vinyl-2-pyrrolidone (0.1µg)  
2; Vinyl Acetate (5.0µg)




AP-0965

### • Nux Vomica

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ 50mmol/l KH<sub>2</sub>PO<sub>4</sub>,  
1.11%Triethylamine = 10/90(pH3.0 with H<sub>3</sub>PO<sub>4</sub>)  
Flow rate: 0.5 ml/min  
Temperature: 30°C  
Detection: UV210nm

Sample: 1; Strychnine (0.8µg)  
2; Barbitol Sodium (1.1µg)




AP-0958

### • Ipecac

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Methanol/ 20mmol/l Sodium *l*-Heptanesulfonate  
(pH4.0 with Acetic Acid) = 50/50  
Flow rate: 1.0 ml/min  
Temperature: 50°C  
Detection: UV283nm

Sample: 1; Cephaeline Hydrobromide (2.0µg)  
2; Emetine Hydrochloride (0.5µg)




AP-0938

### • Belladonna Root

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ 50mmol/l KH<sub>2</sub>PO<sub>4</sub>,  
1%Triethylamine(pH3.5 with H<sub>3</sub>PO<sub>4</sub>)  
= 10/90  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210nm

Sample: 1; Atropine Sulfate (2.0µg)  
2; Brucine Dihydrate (0.5µg)




AP-0909

### • Ryokeijutukanto Extract

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ 0.13%*H*<sub>3</sub>*P*O<sub>4</sub> = 25/75  
Flow rate: 1.0 ml/min  
Temperature: 40°C  
Detection: UV273nm

Sample: 1; *trans*-Cinnamic Acid (0.1µg)




AP-0976

### • Scopolia Rhizome and Scopolia Extract

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ 50mmol/l KH<sub>2</sub>PO<sub>4</sub>,  
0.1%Triethylamine(pH3.5 with H<sub>3</sub>PO<sub>4</sub>)  
= 10/90  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210nm

Sample: 1; Scopolamine Hydrobromide (0.6µg)  
2; Atropine Sulfate (2.0µg)  
3; Brucine Dihydrate (0.4µg)



AP-0989


### 3) Natural Compounds

#### • Carotenes

##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Tetrahydrofuran/Methanol = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV470nm

Sample: 1;  $\alpha$ -Carotene  
 2;  $\beta$ -Carotene




#### • Flavanones

##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 20mmol/l Phosphate buffer(pH2.5) = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm

Sample: 1; Naringin (0.4 $\mu$ g)  
 2; Naringenin (0.2 $\mu$ g)  
 3; Apigenin (0.2 $\mu$ g)




#### • Saikosaponins

##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 45/55  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD, Gain=6

Sample: 1; Saikosaponin c (1.5 $\mu$ g)  
 2; Saikosaponin a (1.5 $\mu$ g)  
 3; Saikosaponin b<sub>2</sub> (1.5 $\mu$ g)  
 4; Saikosaponin b<sub>1</sub> (1.5 $\mu$ g)  
 5; Saikosaponin d (1.5 $\mu$ g)




### 3) Natural Compounds

#### • Hydroxyflavones

##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile / 20mmol/l Phosphate buffer(pH2.5) = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm

Sample: 1; 7-Hydroxyflavone (0.2 $\mu$ g)  
 2; 6-Hydroxyflavone (0.1 $\mu$ g)  
 3; 3-Hydroxyflavone (0.5 $\mu$ g)  
 4; 5-Hydroxyflavone (0.1 $\mu$ g)




#### • Coumarins

##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol / H<sub>2</sub>O = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.16AUFS

Sample: 1; Coumarin (0.84 $\mu$ g)  
 2; Dihydrocoumarin (3.86 $\mu$ g)




#### • Anthraquinone dyes

##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH3) = 75/25  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.16AUFS

Sample: 1; Alizalin (0.05 $\mu$ g)  
 2; Chrysazin (0.1 $\mu$ g)  
 3; Anthrarufin (0.3 $\mu$ g)  
 6; Amylbenzene (10 $\mu$ g)




#### • Chlorophyll

##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: A: Methanol/ 20mmol/l Phosphate buffer(pH7.0) = 60/40  
 B: Methanol/ 20mmol/l Phosphate buffer(pH7.0) = 95/5  
 B conc. 0 $\rightarrow$ 100% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV405nm

Sample: Sodium Copper Chlorophyll (2.0 $\mu$ g)




#### • Hinokitiol

##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 1mmol/l EDTA, 20mmol/l Phosphoric Acid = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.5AUFS

Sample: Hinokitiol (1.0 $\mu$ g)




#### • Anthraquinone dyes

##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH3) = 75/25  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.16AUFS

Sample: 4; Purpurin (0.2 $\mu$ g)  
 5; Quinizarin (0.1 $\mu$ g)  
 6; Amylbenzene (10 $\mu$ g)




### 3) Natural Compounds

#### • Capsaicin

##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm

Sample: Capsaicin (2.0µg)




NACALAI TESQUE, INC  
AP-1034

#### • Alkaloids

##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH3) = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm

Sample: 1; Papaverine Hydrochloride (0.17µg)  
 2; Aconitine (8.1µg)  
 3; Reserpine (1.7µg)



NACALAI TESQUE, INC


AP-0133

#### • Ginkgo Biloba

##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD

Sample: 1; Ginkgolide C (1.3µg)  
 2; Ginkgolide A (1.3µg)  
 3; Ginkgolide B (1.3µg)



NACALAI TESQUE, INC


AP-1063

#### • Zingiberis Rhizoma

##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV225nm

Sample: 1; 6-Gingerol (0.31µg)  
 2; 6-Shogaol (0.13µg)



NACALAI TESQUE, INC


AP-1016

#### • Watermelon

##### COSMOSIL Application Data

Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 20mmol/l Phosphate buffer(pH7.0) = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm

Sample: 1; L-Citrulline (7.5µg)  
 2; Malic Acid (3.8µg)



NACALAI TESQUE, INC

AP-1062


### 4) Pesticides

#### • Asulam

##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 50mmol/l Phosphate buffer(pH3) = 15/85  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV270nm

Sample: Asulam



NACALAI TESQUE, INC


AP-0144

#### • Chlorophenoxyacetic Herbicides

##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 0.1% Acetic Acid = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV275nm, 0.02AUFS

Sample: 1; MCP (0.22µg)  
 2; MCPP (0.20µg)  
 3; MCPB (0.24µg)



NACALAI TESQUE, INC


AP-0147

#### • Pesticides used at Golf Course

##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 50mmol/l Phosphate buffer(pH3) = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV230nm

Sample: 1; Thiram  
 2; Iprodione  
 3; Bensulide



NACALAI TESQUE, INC


AP-0145

#### • Diphenyl Ether Herbicides

##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.1AUFS

Sample: 1; Fluorodifen (0.6µg)  
 2; Chlormethoxynil (1.0µg)  
 3; Nitrofen (1.0µg)  
 4; Oxyfluorfen (1.2µg)  
 5; CNP (1.0µg)



NACALAI TESQUE, INC


AP-0146

#### • Aniline Herbicides

##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.1AUFS

Sample: 1; DCMU (0.11µg)  
 2; DCPA (0.08µg)  
 3; Linuron (0.11µg)  
 4; MCC (0.18µg)



NACALAI TESQUE, INC


AP-0148

#### • Dithiocarbamate Herbicides

##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 20mmol/l KH<sub>2</sub>PO<sub>4</sub> = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample: 1; Maneb (60µg)  
 2; Thiram (10µg)



NACALAI TESQUE, INC

AP-0149


## 4) Pesticides

### • Carbamate Herbicides

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.05AUFS

Sample:  
 1: Methomyl (0.5µg)  
 2: Pirimicarb (0.5µg)  
 3: MTMC (2.0µg)  
 4: PHC (2.0µg)  
 5: MPMC (2.0µg)  
 6: NAC (1.0µg)  
 7: XMC (1.0µg)  
 8: Ethiofencarb (2.0µg)  
 9: Isoprocarb (2.0µg)  
 10: BPMC (2.0µg)




NACALAI TESQUE, INC  
 AP-0150

### • Carbamate Herbicides

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.05AUFS

Sample: Carbosulfan (6.0µg)




NACALAI TESQUE, INC  
 AP-0151

### • Triazine and Urea Herbicides

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 20mmol/l KH<sub>2</sub>PO<sub>4</sub> = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.2AUFS

Sample:  
 1: PAC (0.5µg)  
 2: Bromacil (0.5µg)  
 3: CAT (0.5µg)  
 4: Methabenzthiazuron (0.5µg)  
 5: Chlorotoluron (0.5µg)  
 6: Isoproturon (0.5µg)  
 7: Monolinuron (0.5µg)  
 8: Metabromuron (0.5µg)  
 9: Dimetufun (0.5µg)  
 10: Propazine (0.5µg)  
 11: Terbutylazine (0.5µg)  
 12: Linuron (0.5µg)  
 13: Chloroxuron (0.5µg)  
 14: Prometryn (0.5µg)  
 15: Terbutryn (0.5µg)  
 16: Chloro IPC (2.5µg)  
 17: Ethofumesate (5.0µg)




NACALAI TESQUE, INC  
 AP-0152

### • Triazine and Urea Herbicides

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 20mmol/l KH<sub>2</sub>PO<sub>4</sub> = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.2AUFS

Sample:  
 1: Ethidimuron (1.0µg)  
 2: Metoxuron (1.0µg)  
 3: Cyanazine (1.0µg)  
 4: Methabenzthiazuron (1.0µg)  
 5: Chlorotoluron (1.0µg)  
 6: Atrazine (1.0µg)  
 7: Isoproturon (1.0µg)  
 8: Metabromuron (1.0µg)  
 9: Metazachlor (5.0µg)  
 10: Propazine (1.0µg)  
 11: Terbutylazine (1.0µg)  
 12: Linuron (1.0µg)



NACALAI TESQUE, INC  
 AP-0153


## 5) Food Additives

### • Natural Colorants (Chlorophyll)

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol / 20mmol/l Phosphoric Acid = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.02AUFS

Sample:  
 1: Sodium Copper Chlorophyllin (30µg)  
 2: Sodium Iron Chlorophyllin (30µg)




NACALAI TESQUE, INC  
 AP-0095

### • Natural Colorants (Carotenoid)

#### COSMOSIL Application Data

Column: 5SL-II  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Chloroform/Hexane = 1/9  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.16AUFS

Sample:  
 1: β-Carotene (0.5µg)  
 2: all-trans-Retinol Acetate (1.5µg)




NACALAI TESQUE, INC  
 AP-0094

### • Synthetic Sweeteners

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile / 20mmol/l Phosphoric Acid = 10/90  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm, 1.0AUFS

Sample:  
 1: Acesulfame K (1.0µg)  
 2: Saccharin (1.0µg)  
 3: Diketopiperazine (1.0µg)  
 4: Aspartame (1.0µg)




NACALAI TESQUE, INC  
 AP-0107

### • Synthetic Colorants

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: A: Methanol/ 20mmol Phosphate buffer(pH7) = 10/90  
 B: Methanol/ 20mmol Phosphate buffer(pH7) = 80/20  
 B conc. 15→100% Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.64AUFS

Sample:  
 1: Tartrazine(黄色4号) (1.0µg)  
 2: Amaranth(赤色2号) (1.0µg)  
 3: Indigo Carmine(青色2号) (1.0µg)  
 4: New Coccine(赤色102号) (1.0µg)  
 5: Sunset Yellow FCF(黄色5号) (1.0µg)  
 6: Fast Green FCF(绿色3号) (3.0µg)  
 7: Brilliant Blue FCF(青色1号) (3.0µg)  
 8: Acid Red(赤色106号) (1.0µg)  
 9: Erythrosine(赤色3号) (1.0µg)  
 10: Phloxine(赤色104号) (1.0µg)  
 11: Rose Bengale(赤色105号) (1.0µg)




NACALAI TESQUE, INC  
 AP-0092

### • Synthetic Colorants

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: A: Methanol/ 20mmol Phosphate buffer(pH7) = 10/90  
 B: Methanol/ 20mmol Phosphate buffer(pH7) = 80/20  
 B conc. 15→100% Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.64AUFS

Sample:  
 1: Tartrazine(黄色4号) (1.0µg)  
 2: Amaranth(赤色2号) (1.0µg)  
 3: Indigo Carmine(青色2号) (1.0µg)  
 4: New Coccine(赤色102号) (1.0µg)  
 5: Sunset Yellow FCF(黄色5号) (1.0µg)  
 6: Fast Green FCF(绿色3号) (3.0µg)  
 7: Brilliant Blue FCF(青色1号) (3.0µg)  
 8: Acid Red(赤色106号) (1.0µg)  
 9: Erythrosine(赤色3号) (1.0µg)  
 10: Phloxine(赤色104号) (1.0µg)  
 11: Rose Bengale(赤色105号) (1.0µg)




NACALAI TESQUE, INC  
 AP-0093

### • Food Preservatives

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/Methanol/ 20mmol/l Phosphate buffer(pH4.0) = 20/10/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV230nm

Sample:  
 1: Benzoic Acid (0.075µg)  
 2: Sorbic Acid (0.1µg)  
 3: Dehydroacetic Acid (0.175µg)



NACALAI TESQUE, INC  
 AP-0378


## 5) Food Additives

### • Preservatives (Fungicides)

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV230nm, 0.5AUFS

Sample: 1; Thiabendazole (0.3µg)  
 2; *o*-Phenylphenol (0.3µg)  
 3; Imazalil (2.1µg)  
 4; Diphenyl (0.3µg)




NACALAI TESQUE, INC  
 AP-0097

### • Glycyrrhizic Acid

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile / (Acetic Acid/H<sub>2</sub>O=1/15) = 2/3  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.64AUFS

Sample: 1; Glycyrrhizic Acid (5.0µg)  
 2; Propyl *p*-Hydroxybenzoate (1.0µg)




NACALAI TESQUE, INC  
 AP-0109

### • Antioxidants

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol / H<sub>2</sub>O = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm, 0.16AUFS

Sample: 1; 2-*tert*-Butylhydroquinone (0.56µg)  
 2; 3-*tert*-Butyl-4-hydroxyanisole (0.58µg)  
 3; 2,6-Di-*tert*-butyl-4-hydroxytoluene (1.85µg)




NACALAI TESQUE, INC  
 AP-0100

### • Antioxidants

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol / H<sub>2</sub>O = 90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.32AUFS

Sample: 1; Propyl Gallate (1.0µg)  
 2; Rthoxyquin (1.0µg)



NACALAI TESQUE, INC  
 AP-0101

### • Repellents

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol / H<sub>2</sub>O = 90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.32AUFS

Sample: Piperonyl Butoxide (10µg)




NACALAI TESQUE, INC  
 AP-0104

### • Umami Seasonings (Nucleic Acids)

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 5mmol/l Tributylammonium bromide, 20mmol/l Phosphate buffer(pH7) = 5/95  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV270nm, 0.32AUFS

Sample: 1; Sodium 5'-Cytidylate (1.0µg)  
 2; Sodium 5'-Uridylate (1.0µg)  
 3; Sodium 5'-Guanylate (1.0µg)  
 4; Sodium 5'-Inosinate (1.0µg)



NACALAI TESQUE, INC  
 AP-0106


## 6) Vitamins

### • Hydrosoluble Vitamins

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile / 5mmol/l Sodium *I*-Hexanesulfonate, 20mmol/l Phosphoric Acid = 10/90  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm, 0.16AUFS

Sample: 1; Vitamin C (0.23µg)  
 2; Niacin (0.057µg)  
 3; Nicotinamide (0.042µg)  
 4; Vitamin B<sub>6</sub> (0.040µg)  
 5; Riboflavin Phosphate Sodium Salt (0.19µg)  
 6; Vitamin B<sub>1</sub> Hydrochloride (0.19µg)  
 7; Folic Acid (0.084µg)  
 8; Vitamin B<sub>2</sub> (0.57µg)




NACALAI TESQUE, INC  
 AP-0059

### • Hydrosoluble Vitamins

#### COSMOSIL Application Data

Column: HILIC  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 100mmol/l Ammonium Acetate = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm

Sample: 1; Nicotinamide (0.125µg)  
 2; Pyridoxine(Vitamin B<sub>6</sub>) (0.25µg)  
 3; Riboflavin (Vitamin B<sub>2</sub>) (0.25µg)  
 4; Nicotinic Acid (0.125µg)  
 5; *D*-Pantothenic Acid (3.125µg)  
 6; *L*(+)-Ascorbic Acid (0.875µg)




NACALAI TESQUE, INC  
 AP-0315

### • Water-soluble Vitamins

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A; 20mmol/l Phosphate buffer(pH2.5)  
 B; Methanol/ 20mmol/l Phosphate buffer(pH2.5) = 60/40  
 B conc. 0→80% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm

Sample: 1; Vitamin B<sub>1</sub> Hydrochloride (0.17µg)  
 2; *L*(+)-Ascorbic Acid (0.33µg)  
 3; Nicotinic Acid (0.05µg)  
 4; Nicotinamide (0.05µg)  
 5; Pyridoxine Hydrochloride (0.27µg)  
 6; *D*-Pantothenic Acid Sodium Salt (2.01µg)  
 7; Vitamin B<sub>2</sub> (0.20µg)  
 8; Folic Acid (0.26µg)  
 9; *D*-Biotin [Vitamin H] (2.02µg)  
 10; Flavin Mononucleotide Sodium Salt (0.26µg)  
 11; Vitamin B<sub>2</sub> (Riboflavin) (0.13µg)




NACALAI TESQUE, INC  
 AP-1055

### • Ascorbic Acids

#### COSMOSIL Application Data

Column: HILIC  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 100mmol/l Ammonium Acetate = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Sorbic Acid (0.05µg)  
 2; Isoascorbic Acid (0.3µg)  
 3; Ascorbic Acid (0.3µg)




NACALAI TESQUE, INC  
 AP-0318

### • Energy Drink

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 5mmol/l Sodium *I*-Hexanesulfonate, 20mmol/l Phosphoric Acid = 15/85  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm

Sample: 1; Citric Acid (10mg/ml)  
 2; Carnitine (20mg/ml)  
 3; Nicotinamide (0.2mg/ml)  
 4; Vitamin B<sub>6</sub> (0.2mg/ml)  
 5; Vitamin B<sub>1</sub> (0.2mg/ml)  
 6; Caffeine (0.2mg/ml)  
 7; Riboflavin Phosphate (0.2mg/ml)



NACALAI TESQUE, INC  
 AP-1048


## 6) Vitamins

### • Fruit Juice

#### COSMOSIL Application Data

Column: HILIC  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7.0) = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm

Sample: 1; Ascorbic Acid (1.5µg)  
 2; Malic Acid (3.0µg)




NACALAI TESQUE, INC  
 AP-0313

### • Fat-soluble Vitamins

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/Methanol = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm, 0.64AUFS

Sample: 1; Vitamin K<sub>3</sub> (1.71µg)  
 2; Vitamin A Acetate (0.44µg)  
 3; Vitamin D<sub>3</sub> (0.57µg)  
 4; Vitamin E (4.97µg)  
 5; Vitamin E Acetate (4.79µg)  
 6; Vitamin K<sub>1</sub> (1.68µg)




NACALAI TESQUE, INC  
 AP-0060

### • Vitamin D

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV265nm

Sample: 1; Vitamin D<sub>2</sub> (0.3µg)  
 2; Vitamin D<sub>3</sub> (0.1µg)




NACALAI TESQUE, INC  
 AP-1035

### • Vitamin E

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: 5C<sub>18</sub>-MS-II Methanol  
 πNAP Methanol/ H<sub>2</sub>O = 90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV295nm

Sample: 1; δ-Tocotrienol  
 2; γ-Tocotrienol  
 3; β-Tocotrienol  
 4; α-Tocotrienol  
 5; δ-Tocopherol  
 6; γ-Tocopherol  
 7; β-Tocopherol  
 8; α-Tocopherol



NACALAI TESQUE, INC  
 AP-1071


## 6) Vitamins

### • Tocopherols

#### COSMOSIL Application Data

Column: SSL-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: n-Hexane/2-Propanol = 99.5/0.5  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm, 0.16AUFS

Sample: 1; α-Tocopherol (3.0µg)  
 2; β-Tocopherol (7.5µg)  
 3; γ-Tocopherol (7.5µg)  
 4; δ-Tocopherol (7.5µg)




NACALAI TESQUE, INC  
 AP-0061

### • Vitamin A

#### COSMOSIL Application Data

Column: SSL-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: n-Hexane/Chloroform = 9/1  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.16AUFS

Sample: 1; β-Carotene (0.5µg)  
 2; Vitamin A (1.5µg)




NACALAI TESQUE, INC  
 AP-0062

### • Vitamin A Acid

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol / 20mmol/l Phosphate buffer(pH2.5) = 90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV350nm

Sample: 1; 13-cis-Retinoic Acid (0.04µg)  
 2; all-trans-Retinoic Acid (0.04µg)



NACALAI TESQUE, INC  
 AP-1036


## 7) Metabolites

### • Androgens

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Methanol / H<sub>2</sub>O = 70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV235nm, 0.5AUFS

Sample: 1; 4-Androstene-3,17-dione (1.0µg)  
2; Testosterone (1.0µg)  
3; 17-Methyltestosterone (1.0µg)  
4; Testosterone Propionate (1.0µg)




NACALAI TESQUE, INC  
AP-0111

### • Catecholamines

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: 20mmol/l NaH<sub>2</sub>PO<sub>4</sub>  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV270nm, 0.1AUFS

Sample: 1; Norepinephrine (1.0µg)  
2; Epinephrine (1.2µg)  
3; L-DOPA (1.2µg)  
4; Dopamine (1.0µg)  
5; 3,4-Dihydroxyphenylacetic Acid (1.0µg)  
6; 3-Methoxytyramine (1.0µg)




NACALAI TESQUE, INC  
AP-0120

### • Estrogens

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ H<sub>2</sub>O = 35/65  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV280nm, 0.2AUFS

Sample: 1; Estriol (1.5µg)  
2; Estradiol (1.5µg)  
3; Estrone (1.5µg)




NACALAI TESQUE, INC  
AP-0114

### • Estradiols

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ H<sub>2</sub>O = 40/60  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV280nm

Sample: 1; 17-β-Estradiol (2.0µg)  
2; 17-α-Estradiol (2.0µg)




NACALAI TESQUE, INC  
AP-0243

### • Hippuric Acids

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH3) = 15/85  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV225nm, 0.32AUFS

Sample: 1; Creatinine (0.14µg)  
2; Mandelic Acid (1.89µg)  
3; Hippuric Acid (0.61µg)  
4; o-Methyl Hippuric Acid(0.52µg)




NACALAI TESQUE, INC  
AP-0112

### • Methylhippuric Acids

#### COSMOSIL Application Data

Column: SPYE  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH2.5) = 40/60  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254nm

Sample: N-(o-Toluoyl) glycine (3.06µg)  
N-(m-Toluoyl) glycine (0.65µg)  
N-(p-Toluoyl) glycine (0.33µg)



NACALAI TESQUE, INC  
AP-0289


## 7) Metabolites

### • Urate Metabolites

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: 20mmol/l Phosphoric Acid  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV280nm, 0.1AUFS

Sample: 1; Hypoxanthine (0.9µg)  
2; Uric Acid (18.3µg)  
3; Xanthine (9.0µg)  
4; Allopurinol (1.8µg)




NACALAI TESQUE, INC  
AP-0117

### • Prostaglandins

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: 0.05%TFA-40%Acetonitrile  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: ELSD

Sample: 1; Prostaglandin I<sub>2</sub> (2.0µg)  
2; Prostaglandin F<sub>2</sub> α (2.0µg)  
3; Prostaglandin E<sub>2</sub> (2.0µg)  
4; Prostaglandin D<sub>2</sub> (2.0µg)




NACALAI TESQUE, INC  
AP-0247

### • Adrenal Cortical Hormones

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: 5C<sub>18</sub>-MS-II Methanol/ H<sub>2</sub>O = 50/50  
πNAP Methanol/ H<sub>2</sub>O = 60/40  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254nm

Sample: 1; Triamcinolone  
2; Cortisone  
3; Prednisolone  
4; Hydrocortisone  
5; Betamethasone  
6; Dexamethasone  
7; Triamcinolone Acetonide  
8; Cortisone Acetate  
9; Prednisolone Acetate  
10; Fluocinolone Acetonide  
11; Fluorometholone  
(each 0.23µg)




NACALAI TESQUE, INC  
AP-1037

### • Adrenal Cortical Hormones

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
Column size: 4.6mmI.D.-150mm  
Mobile phase: Acetonitrile/ H<sub>2</sub>O = 30/70  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV240nm

Sample: 1; Prednisolone (0.1µg)  
2; Hydrocortisone (0.1µg)  
3; Cortisone (0.1µg)  
4; Betamethasone (0.1µg)  
5; Dexamethasone (0.1µg)



NACALAI TESQUE, INC  
AP-0233


## 8) Carbohydrates

### • Oligosaccharides

#### COSMOSIL Application Data

Column: Sugar-D  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile / H<sub>2</sub>O = 75/25  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: RI

Sample: 1; Sucrose (10µg)  
2; Maltose (10µg)  
3; Trehalose (10µg)  
4; Raffinose (10µg)




NACALAI TESQUE, INC  
AP-0335

### • Cyclodextrins

#### COSMOSIL Application Data

Column: Sugar-D  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile / H<sub>2</sub>O = 65/35  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: RI

Sample: 1;  $\alpha$ -Cyclodextrin (10µg)  
2;  $\beta$ -Cyclodextrin (10µg)  
3;  $\gamma$ -Cyclodextrin (10µg)



NACALAI TESQUE, INC  
AP-0336

## 8) Carbohydrates


### • Energy Drink

#### COSMOSIL Application Data

Column: Sugar-D  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile / H<sub>2</sub>O = 85/15  
Flow rate: 1.0 ml/min  
Temperature: 40°C  
Detection: ELSD

Sample: 1; Fructose (1.0mg/ml)  
2; Glucose (1.0mg/ml)  
3; Glucitol(Sorbitol) (1.0mg/ml)  
4; Sucrose (1.0mg/ml)  
5; Inositol (1.0mg/ml)

Injection Vol. 5.0µl




NACALAI TESQUE, INC  
AP-1049

### • Xylooligosaccharides

#### COSMOSIL Application Data

Column: Sugar-D  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile / H<sub>2</sub>O = 75/25  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: RI

Sample: Xylooligosaccharides (50µg)  
(Xylobiose, Xylotriose etc.)




NACALAI TESQUE, INC  
AP-0321

### • Isomaltooligosaccharides

#### COSMOSIL Application Data

Column: Sugar-D  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile / H<sub>2</sub>O = 75/25  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: RI

Sample: Isomaltooligosaccharides (50µg)  
(Isomaltose, Isomaltotriose, Panose etc.)




NACALAI TESQUE, INC  
AP-0322

### • Chocolate Components

#### COSMOSIL Application Data

Column: Sugar-D  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile / H<sub>2</sub>O = 75/25  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: RI

Sample: 1; Sucrose (10µg)  
2; Lactose (10µg)




NACALAI TESQUE, INC  
AP-0324

### • Gum Components

#### COSMOSIL Application Data

Column: Sugar-D  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile / H<sub>2</sub>O = 75/25  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: RI

Sample: 1; Xylitol (10µg)  
2; Mannitol (10µg)  
3; Maltitol (10µg)  
4; Palatinin (10µg)  
5; Maltotriitol (10µg)  
6; Maltotetraitol (10µg)




NACALAI TESQUE, INC  
AP-0325

### • Anticarious Foods Components

#### COSMOSIL Application Data

Column: Sugar-D  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile / H<sub>2</sub>O = 75/25  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: RI

Sample: 1; *meso*-Erythritol (10µg)  
2; Xylitol (10µg)  
3; Palatinose (10µg)  
4; Maltitol (10µg)




NACALAI TESQUE, INC  
AP-0323

### • Infusion Solution Components

#### COSMOSIL Application Data

Column: Sugar-D  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile / H<sub>2</sub>O = 85/15  
Flow rate: 1.0 ml/min  
Temperature: 50°C  
Detection: RI

Sample: 1; Xylitol (10µg)  
2; Fructose (10µg)  
3; Glucose (10µg)




NACALAI TESQUE, INC  
AP-0328

### • Cold Beverage Components

#### COSMOSIL Application Data

Column: Sugar-D  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile / H<sub>2</sub>O = 75/25  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: RI

Sample: 1; Fructose (10µg)  
2; Glucose (10µg)  
3; Sucrose (10µg)  
4; Maltose (10µg)




NACALAI TESQUE, INC  
AP-0326

### • Sports Drink Components

#### COSMOSIL Application Data

Column: Sugar-D  
Column size: (4.6mmI.D.-250mm) × 2  
Mobile phase: Acetonitrile / H<sub>2</sub>O = 75/25  
Flow rate: 1.0 ml/min  
Temperature: 50°C  
Detection: RI

Sample: 1; *meso*-Erythritol (10µg)  
2; Fructose (10µg)  
3; Glucitol (10µg)  
4; Glucose (10µg)  
5; Sucrose (10µg)  
6; Maltose (10µg)  
7; Maltotriitol (10µg)  
8; Maltotriose (10µg)



NACALAI TESQUE, INC  
AP-0327


## 8) Carbohydrates

### • PA-Glucose Oligomer

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A: 20mmol/l Acetate buffer(pH3.3)  
 B: 20mmol/l Acetate buffer(pH3.3)+0.5%Butanol  
 B(0→100%) 45min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: FLS at Ex.320nm  
 Em. 400nm

Sample:  
 PA-Glucose Oligomer(DP=3-22)



NACALAI TESQUE, INC


AP-0379

### • Phosphorylated Sugars

#### COSMOSIL Application Data

Column: HILIC  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 20mmol/l Phosphate  
 buffer(pH7.0) = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample:  
 1; D-Fructose-6-phosphate (50µg)  
 2; D-Glucose-6-phosphate (50µg)



NACALAI TESQUE, INC


AP-0317

### • Phosphorylated Sugars

#### COSMOSIL Application Data

Column: HILIC  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 20mmol/l Phosphate  
 buffer(pH7.0) = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample:  
 1; Glucose (10µg)  
 2; α-D-Glucose-1-phosphate (50µg)  
 3; D-Glucose-6-phosphate (50µg)



NACALAI TESQUE, INC


AP-0314

### • Hyaluronic Acid

#### COSMOSIL Application Data

Column: CNT-1000  
 Column size: 7.5mmI.D.-300mm  
 Mobile phase: 20mmol/l Phosphate buffer(pH7),100mmol/l Na<sub>2</sub>SO<sub>4</sub>  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample: Hyaluronic Acid (100µg)



NACALAI TESQUE, INC


AP-1081

### • Enzyme digests of Chondroitin Sulfate A

#### COSMOSIL Application Data

Column: Cholesterol  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A: 1mmol/l Tetrabutylammonium bisulfate/H<sub>2</sub>O  
 B: 1mmol/l Tetrabutylammonium bisulfate  
 -Acetonitrile/H<sub>2</sub>O=67/33  
 B conc. 20%→65%(7min)→65%(12min)  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV240nm

Sample: Chondroitin Sulfate A  
 Chondroitinase AC-II digested



NACALAI TESQUE, INC


AP-1082

### • Enzyme digests of Chondroitin Sulfate A

#### COSMOSIL Application Data

Column: Cholesterol  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A: 1mmol/l Tetrabutylammonium bisulfate/H<sub>2</sub>O  
 B: 1mmol/l Tetrabutylammonium bisulfate  
 -Acetonitrile/H<sub>2</sub>O=67/33  
 B conc. 20%→65%(7min)→65%(12min)  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV240nm

Sample: Chondroitin Sulfate C  
 Chondroitinase AC-II digested



NACALAI TESQUE, INC

AP-1083


## 9) Lipids

### • Fatty Acids

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: 0.05%TFA-90%Methanol  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD

Sample: 1; Oleic Acid (3.0µg)  
 2; Elaidic Acid (3.0µg)  
 3; Stearic Acid (3.0µg)



NACALAI TESQUE, INC


AP-1038

### • Fatty Acids

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: 0.05%TFA-90%Methanol  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD

Sample: 1; cis-Vaccenic Acid (3.0µg)  
 2; trans-Vaccenic Acid (3.0µg)



NACALAI TESQUE, INC


AP-1039

### • Methylated Fatty Acids

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/H<sub>2</sub>O = 95/5  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample: Methyl Linolenate  
 Methyl Linoleate  
 Methyl Oleate  
 Methyl Margarate  
 Methyl Stearate  
 Methyl n-Nonadecanoate  
 (each 10µg)




NACALAI TESQUE, INC


AP-1040

## 9) Lipids


### • Fatty Acid Derivatives




### • Phosphatides




### • Phosphatides




### • Phosphatides



### • Phosphatides




### • Phosphatides




## 10) Nucleic Acid Related Substances


### • Nucleobases and Nucleosides




### • Nucleobases, Nucleosides and Nucleotides




### • Nucleotides




### • Nucleotides



### • DNA Oligomer



### • DNA Ladder




• Dansyl Amino Acids

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A; Acetonitrile/ 20mmol Phosphate buffer(pH7.0) =10/90  
 B; Acetonitrile/ 20mmol Phosphate buffer(pH7.0) =40/60  
 B conc. 0→100% 30min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: Ex.365nm Em.530nm RANGE 2

Sample: Dansyl Amino Acids (1.0µg each)



NACALAI TESQUE, INC  
AP-0003



• PTH-Amino Acids

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A; Acetonitrile/ 20mmol/l Acetate buffer(pH4.8) =10/90  
 B; Acetonitrile/ 20mmol/l Acetate buffer(pH4.8) =50/50  
 Step wise gradient

Flow rate: 1.0 ml/min  
 Temperature: 60°C  
 Detection: UV270nm, 0.16AUFS  
 Sample: PTH-Asp (0.4mg/ml)  
 PTH-His (0.3mg/ml)  
 PTH-Arg (0.5mg/ml)  
 PTH-Ser (0.3mg/ml)  
 PTH-Trp (0.25mg/ml)  
 PTH-Lys (0.35mg/ml)  
 Others (0.2mg/ml)

Injection 2.0µl




NACALAI TESQUE, INC  
AP-0004

• Amino-Acid Derivatives (Diastereomer)

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Boc-(L)Phe-(L)Phe-OMe (7.5µg)  
 2; Boc-(D)Phe-(L)Phe-OMe (7.5µg)




NACALAI TESQUE, INC  
AP-0241

• The umami of Vesitables

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: 5mmol/l Sodium L-Hexanesulfonate, 20mmol/l Phosphoric Acid  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm

Sample: 1; Aspartic Acid (10µg)  
 2; Glutamic Acid (10µg)  
 3; 4-Amino-n-butyric Acid [GABA] (20µg)



NACALAI TESQUE, INC  
AP-1015


• Energy Drink

**COSMOSIL Application Data**

Column: HILIC  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium Acetate = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD

Sample: 1; Taurine (10mg/ml)

Standard (10mg/ml)




NACALAI TESQUE, INC  
AP-1050

• Peptides

**COSMOSIL Application Data**

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A; 0.05%TFA-H<sub>2</sub>O  
 B; 0.05%TFA-Acetonitrile  
 B conc. 10→40% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm

Sample: 1; Oxytocin (0.8µg)  
 2; Angiotensin II(Human) (0.8µg)  
 3; Angiotensin I(Human) (0.8µg)  
 4; Substance P (0.8µg)




NACALAI TESQUE, INC  
AP-0350

• Peptides

**COSMOSIL Application Data**

Column: Cholesterol  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A; 0.05%TFA-H<sub>2</sub>O  
 B; 0.05%TFA-Acetonitrile  
 B conc. 10→40% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm

Sample: 1; Oxytocin (0.8µg)  
 2; Angiotensin II(Human) (0.8µg)  
 3; Angiotensin I(Human) (0.8µg)  
 4; Substance P (0.8µg)



NACALAI TESQUE, INC  
AP-0238


## 11) Amino Acids, Peptides and Proteins

### ● Peptides

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: A; 0.05%TFA-10%Acetonitrile  
 B; 0.05%TFA-30%Acetonitrile  
 B conc. 0→100% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm

Sample:  
 1; Angiotensin II, [Sar<sup>1</sup>,Thr<sup>8</sup>] (0.49µg)  
 2; Angiotensin II, [Sar<sup>1</sup>,Ala<sup>8</sup>] (0.49µg)  
 3; Angiotensin II, Des-Asp<sup>1</sup>-[Ile<sup>6</sup>] (0.49µg)  
 4; Angiotensin II, [Sar<sup>1</sup>,Ile<sup>8</sup>] (0.49µg)  
 5; Angiotensin II, [Asn<sup>1</sup>,Val<sup>7</sup>] (0.49µg)  
 6; Angiotensin II, [Val<sup>7</sup>] (0.49µg)  
 7; Angiotensin II (Human) (0.49µg)




NACALAI TESQUE, INC  
AP-1005

### ● Peptides

#### COSMOSIL Application Data

Column: HILIC  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7.0) = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm

Sample:  
 3; Angiotensin II, Des-Asp<sup>1</sup>-[Ile<sup>8</sup>] (0.5µg)  
 4; Angiotensin II, [Sar<sup>1</sup>,Ile<sup>8</sup>] (0.5µg)  
 5; Angiotensin II, [Asn<sup>1</sup>,Val<sup>7</sup>] (0.5µg)  
 6; Angiotensin II, [Val<sup>7</sup>] (0.5µg)




NACALAI TESQUE, INC  
AP-0316

### ● Peptides

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-300  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: 0.05%TFA-22%Acetonitrile  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.16AUFS

Sample:  
 1; Bradykinin (1.0µg)  
 2; Angiotensin II (1.0µg)  
 3; Neurotensin (1.0µg)  
 4; Bombesin (1.0µg)  
 5; Substance P (1.0µg)




NACALAI TESQUE, INC  
AP-0009

### ● Peptides

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-300  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: 0.05%TFA-22%Acetonitrile  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.16AUFS

Sample:  
 1; Met-Enkephalin (0.5µg)  
 2; [Ala<sup>2</sup>]-Met-Enkephalin (0.5µg)  
 3; Leu-Enkephalin (0.5µg)  
 4; [Ala<sup>2</sup>]-Leu-Enkephalin (0.5µg)




NACALAI TESQUE, INC  
AP-0005

### ● Peptide Mappings

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A; 5mmol/l TFA-H<sub>2</sub>O  
 B; 5mmol/l TFA-60%Acetonitrile  
 B conc. 0→100% 30min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm, 0.32AUFS

Sample: Cytochrome C  
 Lysyl Endopeptidase digested  
 (2.5µl)




NACALAI TESQUE, INC  
AP-0002

### ● Semi-purified Myosin

#### COSMOSIL Application Data

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 10min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.256AUFS

Sample: Myosin (20µg)



NACALAI TESQUE, INC  
AP-0346

## 11) Amino Acids, Peptides and Proteins

### ● Milk Protein


#### COSMOSIL Application Data

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A; 0.05%TFA-H<sub>2</sub>O  
 B; 0.05%TFA-Acetonitrile  
 B conc. 20→80% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm

Sample: Milk

Sample Preparation:  
 •Ultracentrifuged at 90,000 g for 1 hr.  
 •Clear supernatant solution was injected.

Injection vol.: 1.0µl



NACALAI TESQUE, INC  
AP-1079

### ● Soymilk Protein


#### COSMOSIL Application Data

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A; 0.05%TFA-H<sub>2</sub>O  
 B; 0.05%TFA-Acetonitrile  
 B conc. 20→80% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm

Sample: Soybean Milk

Sample Preparation:  
 •Ultracentrifuged at 90,000 g for 1 hr.  
 •Clear supernatant solution was injected.

Injection vol.: 1.0µl




NACALAI TESQUE, INC  
AP-1080

### ● Bacteria-derived Proteins

#### COSMOSIL Application Data

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 10min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.256AUFS

Sample:  
 1; Choline Oxidase (6.0µg)  
 2; α-Amylase (3.0µg)  
 3; Glucose Oxidase (6.0µg)  
 4; Thermolysin (9.0µg)




NACALAI TESQUE, INC  
AP-0337

### ● Bacteria-derived Proteins

#### COSMOSIL Application Data

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 10min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.256AUFS

Sample:  
 1; Actinase E (13.4µg)  
 2; Alcohol Dehydrogenase (6.6µg)




NACALAI TESQUE, INC  
AP-0338

### ● Human-derived Proteins

#### COSMOSIL Application Data

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 10min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.256AUFS

Sample:  
 1; Transferrin (1.5µg)  
 2; Albumin(Human) (2.0µg)  
 3; Carbonic Anhydrase (1.0µg)




NACALAI TESQUE, INC  
AP-0339

### ● Bovine-derived Proteins

#### COSMOSIL Application Data

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 10min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.256AUFS

Sample:  
 1; Albumin(Bovine) (2.0µg)  
 2; L-Glutamic Dehydrogenase (1.5µg)  
 3; Carbonic Anhydrase (1.5µg)



NACALAI TESQUE, INC  
AP-0340


## 11) Amino Acids, Peptides and Proteins

### • Bovine-derived Proteins

#### COSMOSIL Application Data

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 10min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.256AUFS

Sample: 1; Fibrinogen (4.0µg)  
 2; Catalase (2.0µg)  
 3; Thyroglobulin (6.0µg)




NACALAI TESQUE, INC  
 AP-0341

### • Bovine Spleen-derived Proteins

#### COSMOSIL Application Data

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 10min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.256AUFS

Sample: 1; Ribonuclease A (1.0µg)  
 2; α-Chymotrypsinogen A (1.0µg)  
 3; Deoxyribonuclease I (3.0µg)




NACALAI TESQUE, INC  
 AP-0342

### • Egg-derived Proteins

#### COSMOSIL Application Data

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 10min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.256AUFS

Sample: 1; Lysozyme (1.0µg)  
 2; Conalbumin (2.0µg)  
 3; Albumin(Ovalbumin)(1.8µg)




NACALAI TESQUE, INC  
 AP-0343

### • Other Proteins

#### COSMOSIL Application Data

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 10min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.256AUFS

Sample: 1; Cytochrome C (1.0µg)  
 2; Albumin (Goat) (2.0µg)  
 3; Myoglobin (1.0µg)  
 4; Concanavalin A (0.5µg)




NACALAI TESQUE, INC  
 AP-0344

### • Other Proteins

#### COSMOSIL Application Data

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 10min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.256AUFS

Sample: 1; Albumin(Rabbit) (3.2µg)  
 2; Peroxidase (3.2µg)




NACALAI TESQUE, INC  
 AP-0345

### • Pyruvate Kinase

#### COSMOSIL Application Data

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 10min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.256AUFS

Sample: Pyruvate Kinase (4.0µg)



NACALAI TESQUE, INC  
 AP-0348


## 11) Amino Acids, Peptides and Proteins

### • Semi-purified Diaphorase

#### COSMOSIL Application Data

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 10min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.256AUFS

Sample: Diaphorase (6.0µg)




NACALAI TESQUE, INC  
 AP-0347

### • Glycoproteins

#### COSMOSIL Application Data

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm

Sample: Concanavalin A (1.0µg)




NACALAI TESQUE, INC  
 AP-0352

### • Proteins

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-300  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.64AUFS

Sample: 1; Insulin (1.5µg)  
 2; Trypsinogen (6.0µg)  
 3; Transferrin (4.0µg)  
 4; Trypsin Inhibitor, Soybean (5.0µg)  
 5; α-Chymotrypsinogen A (4.0µg)  
 6; Carbonic Anhydrase (3.0µg)




NACALAI TESQUE, INC  
 AP-0017

### • Proteins

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-300  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.64AUFS

Sample: 1; Insulin (1.5µg)  
 2; Trypsinogen (6.0µg)  
 3; Transferrin (4.0µg)  
 4; Trypsin Inhibitor, Soybean (5.0µg)  
 5; α-Chymotrypsinogen A (4.0µg)  
 6; Carbonic Anhydrase (3.0µg)




NACALAI TESQUE, INC  
 AP-0018

### • Proteins

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-300  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.64AUFS

Sample: 1; Insulin (1.5µg)  
 2; Trypsinogen (6.0µg)  
 3; Transferrin (4.0µg)  
 4; Trypsin Inhibitor, Soybean (5.0µg)  
 5; α-Chymotrypsinogen A (4.0µg)  
 6; Carbonic Anhydrase (3.0µg)




NACALAI TESQUE, INC  
 AP-0019

### • Proteins

#### COSMOSIL Application Data

Column: 5Ph-AR-300  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.64AUFS

Sample: 1; Insulin (1.5µg)  
 2; Trypsinogen (6.0µg)  
 3; Transferrin (4.0µg)  
 4; Trypsin Inhibitor, Soybean (5.0µg)  
 5; α-Chymotrypsinogen A (4.0µg)  
 6; Carbonic Anhydrase (3.0µg)



NACALAI TESQUE, INC  
 AP-0020


## 11) Amino Acids, Peptides and Proteins

### • Human Serum

#### COSMOSIL Application Data

Column: 5Diol-120-II  
 Column size: 7.5mmI.D.-600mm  
 Mobile phase: 20mmol/l Phosphate buffer(pH 7.0)  
 +100mmol/l Na<sub>2</sub>SO<sub>4</sub>  
 Flow rate: 1.0 ml/min  
 Temperature: Room temperature  
 Detection: UV280nm

Sample: Human Serum




NACALAI TESQUE, INC  
 AP-0382

### • Human Serum

#### COSMOSIL Application Data

Column: 5Diol-300-II  
 Column size: 7.5mmI.D.-600mm  
 Mobile phase: 20mmol/l Phosphate buffer(pH 7.0)  
 +100mmol/l Na<sub>2</sub>SO<sub>4</sub>  
 Flow rate: 1.0 ml/min  
 Temperature: Room temperature  
 Detection: UV280nm

Sample: Human Serum




NACALAI TESQUE, INC  
 AP-0383

### • Bovine Serum

#### COSMOSIL Application Data

Column: 5Diol-120-II  
 Column size: 7.5mmI.D.-600mm  
 Mobile phase: 20mmol/l Phosphate buffer(pH 7.0)  
 +100mmol/l Na<sub>2</sub>SO<sub>4</sub>  
 Flow rate: 1.0 ml/min  
 Temperature: Room temperature  
 Detection: UV280nm

Sample: Bovine Serum




NACALAI TESQUE, INC  
 AP-0386

### • Bovine Serum

#### COSMOSIL Application Data

Column: 5Diol-300-II  
 Column size: 7.5mmI.D.-600mm  
 Mobile phase: 20mmol/l Phosphate buffer(pH 7.0)  
 +100mmol/l Na<sub>2</sub>SO<sub>4</sub>  
 Flow rate: 1.0 ml/min  
 Temperature: Room temperature  
 Detection: UV280nm

Sample: Bovine Serum




NACALAI TESQUE, INC  
 AP-0387

### • Egg White

#### COSMOSIL Application Data

Column: 5Diol-120-II  
 Column size: 7.5mmI.D.-600mm  
 Mobile phase: 20mmol/l Phosphate buffer(pH 7.0)  
 +100mmol/l Na<sub>2</sub>SO<sub>4</sub>  
 Flow rate: 1.0 ml/min  
 Temperature: Room temperature  
 Detection: UV280nm

Sample: Egg White




NACALAI TESQUE, INC  
 AP-0384

### • Egg White

#### COSMOSIL Application Data

Column: 5Diol-300-II  
 Column size: 7.5mmI.D.-600mm  
 Mobile phase: 20mmol/l Phosphate buffer(pH 7.0)  
 +100mmol/l Na<sub>2</sub>SO<sub>4</sub>  
 Flow rate: 1.0 ml/min  
 Temperature: Room temperature  
 Detection: UV280nm

Sample: Egg White



NACALAI TESQUE, INC  
 AP-0385


## 12) The others

### • Methoxyphenols

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: *o*-Methoxyphenol (3.3µg)  
*m*-Methoxyphenol (3.3µg)  
*p*-Methoxyphenol (3.3µg)




NACALAI TESQUE, INC  
 AP-1041

### • Terphenyls

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: *o*-Terphenyl (0.15µg)  
*m*-Terphenyl (0.05µg)  
*p*-Terphenyl (0.075µg)




NACALAI TESQUE, INC  
 AP-1042

### • Chlorophenols

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm


Sample: *o*-Chlorophenol (2.0µg)  
*m*-Chlorophenol (2.0µg)  
*p*-Chlorophenol (4.0µg)




NACALAI TESQUE, INC  
 AP-1043

## 12) The others


### • Stilbenes



### • Halogenated Benzenes




### • Fluorinated Benzenes




## 12) The others


### • Pyridine and Phenol




### • Pyridine and Phenol




### • Uracil and Uridine




### • Glycerol



### • Ethylene Glycol



### • Diethylene Glycol




## 12) The others

### • Oxalic Acid

#### COSMOSIL Application Data

Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7.0) = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm

Sample: 1; Oxamic Acid (0.2µg)  
 2; Oxalic Acid (1.0µg)




NACALAI TESQUE, INC  
 AP-0307

### • Hydrophilic Compounds (Ionicity)

#### COSMOSIL Application Data

Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 50mmol/l Ammonium Acetate = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD

Sample: 1; *meso*-Erythritol (3.0µg)  
 2; Tris(hydroxymethyl)aminomethane (3.0µg)  
 3; Glyceric Acid (6.0µg)




NACALAI TESQUE, INC  
 AP-0305

### • Fluorine Compounds

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: 0.1%TFA-95%Methanol  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm

Sample: Perfluorotetradecanoic Acid (50µg)




NACALAI TESQUE, INC  
 AP-0381

### • Urea

#### COSMOSIL Application Data

Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/H<sub>2</sub>O = 90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm

Sample: Urea (20µg)




NACALAI TESQUE, INC  
 AP-1004

### • Acid Compounds

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 0.1% Phosphoric Acid = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.32AUFS

Sample: 1; Uracil (0.05µg)  
 2; Benzoic Acid (1.0µg)  
 3; *o*-Toluic Acid (1.0µg)  
 4; *p*-Ethylbenzoic Acid (0.2µg)  
 5; Benzene




NACALAI TESQUE, INC  
 AP-0164

### • Basic Compounds

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH3) = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.2AUFS

Sample: 1; 2-Ethylpyridine (0.4µg)  
 2; *N,N*-Dimethylaniline (0.6µg)



NACALAI TESQUE, INC  
 AP-0166


## 12) The others

### • Monosubstituted Benzenes (20 samples)

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol / H<sub>2</sub>O = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.16AUFS

Sample: 1; Benzamide (0.49µg) 14; Chlorobenzene (7.05µg)  
 2; Aniline (0.4µg) 15; Toluene (5.84µg)  
 3; Phenol (0.67µg) 16; Bromobenzene (15.37µg)  
 4; Benzonitrile (0.83µg) 17; Iodobenzene (3.66µg)  
 5; Acetophenone (0.04µg) 18; Ethylbenzene (6.87µg)  
 6; Styrene oxide (1.1µg) 19; *n*-Propylbenzene (14.1µg)  
 7; Nitrobenzene (0.06µg) 20; *n*-Butylbenzene (15.93µg)  
 8; Methylbenzoate (0.62µg)  
 9; Anisole (0.79µg)  
 10; Fluorobenzene (0.65µg)  
 11; Benzene (0.79µg)  
 12; *N,N*-Dimethylamine (0.15µg)  
 13; Thioanisole (0.12µg)




NACALAI TESQUE, INC  
 AP-0154

### • Furans

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.2AUFS

Sample: 1; Furfuryl alcohol (0.13µg)  
 2; Furfural (0.25µg)  
 3; Furan (0.23µg)




NACALAI TESQUE, INC  
 AP-0155

### • Phenolphthalein

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol / H<sub>2</sub>O = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.1AUFS

Sample: Phenolphthalein (0.6µg)




NACALAI TESQUE, INC  
 AP-0156

### • Anilines

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile / 0.1% Phosphoric Acid = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Aniline  
 2; *p*-Nitroaniline  
 3; 2,4-Dinitroaniline




NACALAI TESQUE, INC  
 AP-0160

### • Esters and others

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.32AUFS

Sample: 1; Uracil  
 2; Pyridine  
 3; Phenol  
 4; Ethyl *p*-Hydroxybenzoate  
 5; Propyl *p*-Hydroxybenzoate  
 6; Methyl Salicylate  
 7; Toluene  
 8; Ethyl Salicylate




NACALAI TESQUE, INC  
 AP-0163

### • Aminophenols

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: 20mmol/l Phosphate buffer(pH7.0)  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: *p*-Aminophenol (0.33µg)  
*m*-Aminophenol (1.00µg)  
*o*-Aminophenol (1.67µg)



NACALAI TESQUE, INC  
 AP-1059


## 12) The others

### • Polyaromatic Compounds

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A: Methanol / H<sub>2</sub>O = 70/30  
 B: Methanol  
 B conc. 0→100%  
 3→15min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.64AUFS

Sample:  
 1; Naphthalene  
 2; Acenaphthylene  
 3; Acenaphthene  
 4; Fluorene  
 5; Phenanthrene  
 6; Anthracene  
 7; Fluoranthene  
 8; Pyrene  
 9; Benz[*a*]anthracene  
 10; Chrysene  
 11; Benz[*b*]fluoranthene  
 12; Benz[*a*]pyrene  
 13; Benz[*a*]pyrene  
 14; Dibenz[*a,h*]anthracene  
 15; Benz[*g,h,i*]perylene  
 16; Indeno[1,2,3-*c,d*]pyrene




NACALAI TESQUE, INC  
 AP-0161

### • Benzylpyridines

#### COSMOSIL Application Data

Column: 5PYE  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l KH<sub>2</sub>PO<sub>4</sub> = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample:  
 1; 2-Benzylpyridine  
 2; 3-Benzylpyridine  
 3; 4-Benzylpyridine




NACALAI TESQUE, INC  
 AP-0182

### • Carvylacetate

#### COSMOSIL Application Data

Column: 5PYE  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol / H<sub>2</sub>O = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV215nm

Sample: Carvylacetate (*cis*, *trans*)





NACALAI TESQUE, INC  
 AP-0184

### • Diastereomers

#### COSMOSIL Application Data

Column: 5PYE  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol / H<sub>2</sub>O = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample:  
 1; *Threo* form  
 2; *Erythro* form




NACALAI TESQUE, INC  
 AP-0186

### • Surfactant

#### COSMOSIL Application Data

Column: 5PBB-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Tetrahydrofuran / Methanol = 10/90  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample: Polyoxyethylene hexadecyl ether (500µg)




NACALAI TESQUE, INC  
 AP-0354

### • Surfactant

#### COSMOSIL Application Data

Column: 5PBB-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Tetrahydrofuran / Methanol = 10/90  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample: Polyoxyethylene Lauryl Ether (500µg)



NACALAI TESQUE, INC  
 AP-0364


## 12) The others

### • Aromatic Compounds

#### COSMOSIL Application Data

Column: 5PBB-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/H<sub>2</sub>O = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample:  
 1; Diphenylmethane (5.80µg)  
 2; Fluorene (0.13µg)




NACALAI TESQUE, INC  
 AP-0358

### • Aromatic Compounds

#### COSMOSIL Application Data

Column: 5PBB-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/H<sub>2</sub>O = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample:  
 1; Dibenzosuberone (0.25µg)  
 2; Dibenzosuberone (0.10µg)




NACALAI TESQUE, INC  
 AP-0360

### • Aromatic Compounds

#### COSMOSIL Application Data

Column: 5PBB-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/H<sub>2</sub>O = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample:  
 1; Propylbenzene (27.9µg)  
 2; Allylbenzene (0.45µg)




NACALAI TESQUE, INC  
 AP-0362

### • Aromatic Compounds

#### COSMOSIL Application Data

Column: 5PBB-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/H<sub>2</sub>O = 90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample:  
 1; 1,1'-Binaphthyl (0.60µg)  
 2; Perylene (0.72µg)




NACALAI TESQUE, INC  
 AP-0356

### • Coloring agent (CBB)

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH2.5) = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: Coomassie Brilliant Blue G-250 (3.0µg)



NACALAI TESQUE, INC  
 AP-1061


## 4. Reference list

No.	Column	TITLE	AUTHOR	JOURNAL	YEAR	VOL. (ISSUE)	PAGE
1	Cholester	Improved separation of furocoumarins of essential oils by supercritical fluid chromatography	C. Desmorteux, M. Rothaupt, C. West and E. Lesellier	Journal of Chromatography A	2009	1216(42)	7088-7095
2	Cholester	Preparative Separation of cis- and trans-Isomers of Unsaturated Fatty Acid Methyl Esters Contained in Edible Oils by Reversed-Phase High-Performance Liquid Chromatography.	Tsuzuki W, Ushida K.	Lipids.	2009	44(4)	373-379
3	Cholester 5C <sub>18</sub> -AR-II	Isolation and Structure of a Galactocerebroside from the Sea Cucumber <i>Bohadschia argus</i>	Yuriko Ikeda, Masanori Inagaki, Koji Yamada, Xiao Wen Zhang, Bo Zhang, Tomofumi Miyamoto and Ryuichi Higuchi	CHEMICAL & PHARMACEUTICAL BULLETIN	2009	57(3)	315-317
4	Cholester	Influence of Cold Hardening on Chlorophyll and Carotenoid in <i>Chlorella vulgaris</i>	Watanabe Yuta, Yamada Naotaka, Machida Takeshi, Honjoh Ken-ichi, Kuwano Eiichi	Journal of the Faculty of Agriculture, Kyushu University	2009	54(1)	195-200
5	Cholester	Pyripyropenes, Fungal Sesquiterpenes Conjugated with $\alpha$ -Pyrone and Pyridine Moieties, Exhibits Anti-angiogenic Activity against Human Umbilical Vein Endothelial Cells	Asami Hayashi, Masayoshi Arai, Mayumi Fujita and Motomasa Kobayashi	Biological & Pharmaceutical Bulletin	2009	32(7)	1261-1265
6	Cholester	Practical Electrochemical Iodination of Aromatic Compounds	Kazuhide Kataoka, Yuji Hagiwara, Koji Midorikawa, Seiji Suga and Jun-ichi Yoshida	Org. Process Res. Dev.	2008	12(6)	1130-1136
7	Cholester	De novo synthesis of (Z)- and (E)-7-hexadecenylnitric acids by a selective lignin-degrading fungus, <i>Ceriporiopsis subvermispora</i>	Hiroshi Nishimura, Saeko Tsuda, Hito Shimizu, Yasunori Ohashi, Takahito Watanabe, Yoichi Honda and Takashi Watanabe	Phytochemistry	2008	69(14)	2593-2602
8	Cholester	Kadsuracoccinic Acids A-C, Ring-A seco-Lanostane Triterpenes from <i>Kadsura coccinea</i> and Their Effects on Embryonic Cell Division of <i>Xenopus laevis</i>	Heran Li, Liyan Wang, Syohei Miyata and Susumu Kitanaka	J. Nat. Prod.	2008	71(4)	739-741
9	Cholester	Peptide Thioester Synthesis via an Auxiliary-Mediated N-S Acyl Shift Reaction in Solution	Ken'ichiro Nakamura, Hiroaki Mori, Toru Kawakami, Hironobu Hojo, Yoshiaki Nakahara and Saburo Aimoto	International Journal of Peptide Research and Therapeutics	2007	13(1-2)	191-202
10	Cholester	Metachromins L-Q, new sesquiterpenoid quinones with an amino acid residue from sponge <i>Spongia</i> sp.	Yohei Takahashi, Takaaki Kubota, Jane Fromont and Jun'ichi Kobayashi	Tetrahedron	2007	63(36)	8770-8773
11	Cholester	Metachromins R-T, New Sesquiterpenoids from Marine Sponge <i>Spongia</i> sp.	Yohei Takahashi, Mika Yamada, Takaaki Kubota, Jane Fromont and Jun'ichi Kobayashi	CHEMICAL & PHARMACEUTICAL BULLETIN	2007	55(No.12)	1731-1733
12	$\pi$ NAP PYE PBB	Possibility of predicting separations in supercritical fluid chromatography with the solvation parameter model	C. West, J. Ogden and E. Lesellier	Journal of Chromatography A	2009	1216 (29)	5600-5607
13	PYE	Possibility of predicting separations in supercritical fluid chromatography with the solvation parameter model	C. West, J. Ogden and E. Lesellier	Journal of Chromatography A	2009	1216 (29)	5600-5607
14	PYE	Calculations and assignments of endohedral helium-3 chemical shifts of open-cage fullerenes and higher fullerenes	Guan-Wu Wang and Ping Wu	Theoretical Chemistry Accounts	2009	123 (5-6)	375-381
15	PYE	Effects of $\pi$ - $\pi$ Interactions on the Separation of PAHs on Phenyl-Type Stationary Phases	Paul G. Stevenson; Sindy Kayillo; Gary R. Dennis; R. Andrew Shalliker	Journal of Liquid Chromatography & Related Technologies	2008	31(3)	324-347
16	PYE	Organochlorine contaminants in endangered Steller sea lion pups ( <i>Eumetopias jubatus</i> ) from western Alaska and the Russian Far East	Matthew J. Myers, Gina M. Ylitalo, Margaret M. Krahn, Daryle Boyd, Don Calkins, Vladimir Burkanov and Shannon Atkinson	Science of The Total Environment	2008	396 (1)	60-69
17	PYE	Spatial extent, magnitude, and patterns of persistent organochlorine pollutants in Pacific herring ( <i>Clupea pallasii</i> ) populations in the Puget Sound (USA) and Strait of Georgia (Canada)	James E. West, Sandra M. O'Neill, and Gina M. Ylitalo	Science of The Total Environment	2008	394 (2-3)	369-378
18	PYE	Polychlorinated dibenzo-p-dioxins and dibenzofurans (PCDD/Fs) in marine mussels from French coasts: Levels, patterns and temporal trends from 1981 to 2005	C. Munsch, N. Guiot, K. Héas-Moisan, C. Tixier and J. Tronczynski	Chemosphere	2008	73 (6)	945-953
19	PYE	Analysis of chlorothalonil and degradation products in soil and water by GC/MS and LC/MS	Alicia Chaves, Damian Shea, David Danehower	Chemosphere	2008	71(4)	629-638
20	PYE	Polychlorinated Naphthalenes and Other Dioxin-Like Compounds in Elbe River Sediments	W Brack, L Bláha, J.P Giesy, M Grote, M Moeder, S Schrader and M Hecker	Environmental Toxicology and Chemistry	2008	27 (3)	519-528

No.	Column	TITLE	AUTHOR	JOURNAL	YEAR	VOL. (ISSUE)	PAGE
21	PYE	Extraction and high-performance liquid chromatographic analysis of C <sub>60</sub> , C <sub>70</sub> , and [6,6]-phenyl C <sub>61</sub> -butyric acid methyl ester in synthetic and natural waters	Dermont Bouchard, and Xin Ma	Journal of Chromatography A	2008	1203 (2)	153-159
22	PYE	Orthogonal screening system of columns for supercritical fluid chromatography	C. West, and E. Lesellier	Journal of Chromatography A	2008	1203 (1)	105-113
23	PYE	A unified classification of stationary phases for packed column supercritical fluid chromatography	C. West, E. Lesellier	Journal of Chromatography A	2008	1191(1-2)	21-39
24	PYE	Automated fractionation procedure for polycyclic aromatic compounds in sediment extracts on three coupled normal-phase high-performance liquid chromatography columns	Urte Lübcke-von Varel, Georg Streck, Werner Brack	Journal of Chromatography A	2008	1185 (1)	31-42
25	PYE	Characterisation of stationary phases in supercritical fluid chromatography with the solvation parameter model: V. Elaboration of a reduced set of test solutes for rapid evaluation	C. West, E. Lesellier	Journal of Chromatography A	2007	1169 (1-2)	205-219
26	PYE	Fractionation of chlorinated and brominated persistent organic pollutants in several food samples by pyrenyl-silica liquid chromatography prior to GC-MS determination	Belén Gómara, Carmen García-Ruiz, María José González, María Luisa Marina	Analytica Chimica Acta	2006	565 (2)	208-213
27	PYE	Rapid Method for Determination of Dioxin-Like Polychlorinated Biphenyls and Other Congeners in Marine Sediments Using Sonic Extraction and Photodiode Array Detection	J. Buzitis, G. M. Ylitalo and M. M. Krahn	Archives of Environmental Contamination and Toxicology	2006	51 (3)	337-346
28	PYE	Occurrence of Polybrominated Biphenyls, Polybrominated Dibenzo-p-dioxins, and Polybrominated Dibenzofurans as Impurities in Commercial Polybrominated Diphenyl Ether Mixtures	Nobuyasu Hanari, Kurunthachalam Kannan, Yuichi Miyake, Tsuyoshi Okazawa, Prasada Rao S. Kodavanti, Kenneth M. Aldous, and Nobuyoshi Yamashita	Environ. Sci. Technol.	2006	40 (14)	4400-4405
29	NPE 5C <sub>18</sub> -MS	Synthesis of cyclic bis(3'-5')diguanilyc acid (c-di-GMP) analogs	Mamoru Hyodo, Yumi Sato and Yoshihiro Hayakawa	Tetrahedron	2006	62(13)	3089-3094
30	HILIC	Development and validation of a reversed-phase high-performance liquid chromatographic method for quantification of peptide dendrimers in human skin permeation experiments	S. Mutalik, A.K. Hewavitharana, P.N. Shaw, Y.G. Anissimov, M.S. Roberts and H.S. Parekh	Journal of Chromatography B	2009	877(29)	3556-3562
31	HILIC	Direct Evidence for Efficient Transport and Minimal Metabolism of L-Cephalaxin by Oligopeptide Transporter 1 in Budded Baculovirus Fraction	Keisuke Mitsuoka, Ikumi Tamai, Yasushi Morohashi, Yoshiyuki Kubo, Ryoichi Saitoh, Akira Tsuji and Yukio Kato	Biological & Pharmaceutical Bulletin	2009	32(8)	1459-1461
32	HILIC 5C <sub>18</sub> -MS-II	Structures of Acetylated Oleanane-type Triterpene Saponins, Rarasaponins IV, V, and VI, and Anti-hyperlipidemic Constituents from the Pericarps of <i>Sapindus rarak</i>	Yasunobu ASAO, Toshio MORIKAWA, Yuanyuan XIE, Masaki OKAMOTO, Makoto HAMAOKA, Hisashi MATSUDA, Osamu MURAOKA, Dan YUAN, and Masayuki YOSHIKAWA	Chem. Pharm. Bull.	2009	57(2)	198-203
33	HILIC	Simultaneous measurement of diazolidinyl urea, urea, and allantoin in cosmetic samples by hydrophilic interaction chromatography	Takahiro Doi, Keiji Kajimura, Satoshi Takatori, Naoki Fukui, Shuzo Taguchi and Shozo Iwagami	Journal of Chromatography B	2009	877(10)	1005-1010
34	HILIC	Convergent synthesis of oligomers of triazole-linked DNA analogue (TLDNA) in solution phase	Tomoko Fujino, Naomi Yamazaki and Hiroyuki Isobe	Tetrahedron Letters	2009	50(28)	4101-4103
35	HILIC	Determination of para-aminohippuric acid (PAH) in human plasma and urine by liquid chromatography-tandem mass spectrometry	Phey Yen Han, P. Nicholas Shaw and Carl M.J. Kirkpatrick	Journal of Chromatography B	2009	877(27)	3215-3220
36	HILIC 5C <sub>18</sub> -MS-II	Medicinal Flowers. XXI. Structures of Perennisaponins A, B, C, D, E, and F, Acylated Oleanane-Type Triterpene Oligoglycosides, from the Flowers of <i>Bellis perennis</i>	Masayuki Yoshikawa, Xuezheng Li, Eriko Nishida, Seikou Nakamura, Hisashi Matsuda, Osamu Muraoka and Toshio Morikawa	CHEMICAL & PHARMACEUTICAL BULLETIN	2008	56(No.4)	559-568
37	HILIC	Separation efficiencies in hydrophilic interaction chromatography	Tohru Ikegami, Kouki Tomomatsu, Hirota Takubo, Kanta Horie and Nobuo Tanaka	Journal of Chromatography A	2008	1184(1-2)	474-503
38	HILIC	Triazole-Linked Analogue of Deoxyribonucleic Acid (TLDNA): Design, Synthesis, and Double-Strand Formation with Natural DNA	Hiroyuki Isobe, Tomoko Fujino, Naomi Yamazaki, Marine Guillot-Nieckowski and Eiichi Nakamura	ORGANIC LETTERS	2008	10(17)	3729-3732
39	HILIC	Tetradotoxin poisoning evidenced by solid-phase extraction combining with liquid chromatography-tandem mass spectrometry	Hsiao-Chin Jen, Shin-Jung Lin, Yung-Hsiang Tsai, Chun-Hsiang Chen, Zu-Chun Lin and Deng-Fwu Hwang	Journal of Chromatography B	2008	871(1)	95-100
40	CNT-300 CNT-1000 CNT-2000	Chromatographic Length-Separation and Photoluminescence Study on DNA-Wrapped Single-Wall and Double-Wall Carbon Nanotubes	Yuki Asada, Toshiaki Sugai, Ryo Kitaura and Hisanori Shinohara	Journal of nanomaterials		in press	

# IV TECHNICAL NOTE

1. Selectivity of packing materials in reversed phase liquid chromatography ...	154
2. Preparation of mobile phase for HPLC .....	161
3. Sample pretreatment for HPLC .....	164
4. Baseline noise in gradient elution .....	168
5. Troubleshooting for increased pressure .....	169
6. Effect of guard column .....	172
7. Troubleshooting for normal phase chromatography .....	173
8. Inner diameter of column (scale down and scale up) .....	175








# 1. Selectivity of packing materials in reversed phase liquid chromatography

Reversed phase chromatography is the most commonly used method of HPLC, because of the high theoretical plate number, excellent separation characteristics, reproducibility, and ease of use. Columns packed with octadecyl group bonded type silica gel (C<sub>18</sub>, ODS) are the most widely used reversed phase chromatography. However, C<sub>18</sub> columns provide insufficient separation for compounds similar in hydrophobicity because the main separation mechanism of C<sub>18</sub> column is based on hydrophobic interaction. It may improve separation of compounds with similar hydrophobicity by using longer columns, changing mobile phases or changing temperature. However, in many cases, it is probably most effective to use different packing materials which retain compounds based on a secondary interaction in addition to hydrophobic interaction.

At Nacalai Tesque, we offer a variety of COSMOSIL reversed phase packing materials. Summary of these packing materials and their respective retention mechanism are in Table 1. Retention of compounds in each stationary phase depends on summation of the interactions. Therefore, comprehension of each interaction leads to selection of an appropriate column.


Table 1. Stationary phase and interaction of packing materials

	C <sub>18</sub> -MS-II	C <sub>18</sub> -AR-II	C <sub>8</sub> -MS	PE-MS	πNAP	PYE	NPE	PBB-R	Cholester
Silica gel	High purity porous spherical silica								
Average particle size	3 · 5 · 15 μm			5 μm					
Average pore size	approx. 120 Å								
Specific surface area	approx. 300 m <sup>2</sup> /g								
Stationary phase									
Types	Monomeric	Polymeric	Monomeric	Monomeric	Monomeric	Monomeric	Monomeric	Monomeric	Monomeric
Interactions	Hydrophobic interaction	Hydrophobic interaction	Hydrophobic interaction	Hydrophobic interaction π-π interaction	Hydrophobic interaction π-π interaction	Hydrophobic interaction π-π interaction Dispersion interaction	Hydrophobic interaction π-π interaction Dipole-dipole interaction	Hydrophobic interaction Dispersion interaction	Hydrophobic interaction Shape selectivity
End capping	Near-perfect treatment								
Carbon content	approx. 16%	approx. 17%	approx. 10%	approx. 10%	approx. 11%	approx. 18%	approx. 9%	approx. 8%	approx. 20%

## 1) Selectivity for polar functional group

### Selectivity

Selectivity for polar functional group is evaluated based on the separation of benzene, nitrobenzene, which has a nitro group, and anisole, which has a methoxy group. The chromatograms below show separation of the three compounds on four COSMOSIL columns: C<sub>18</sub>-MS-II, PE-MS, πNAP and PYE. Elution order on the C<sub>18</sub> column is as following: nitrobenzene, anisole and benzene. Elution orders on the aromatic columns are reversed. Separation on the C<sub>18</sub> column is based on hydrophobic interaction only. On the other hand, the packing materials on the other three columns have aromatic rings and reverse the elution order by π-π interaction.





The graph of selectivity for polar functional group is shown below. Among nine COSMOSIL columns, PYE and NPE columns have the highest selectivity factors for polar groups. As to mobile phases, methanol is more effective than acetonitrile for separation using π-π interaction.

### Selectivity for polar functional group

Column: 5C<sub>18</sub>-MS-II, 5PE-MS, πNAP, 5PYE  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol / H<sub>2</sub>O = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Nitrobenzene (0.13 μg)  
 2; Anisole (1.5 μg)  
 3; Benzene (4.0 μg)



### Application data



#### Separation of toluenitrile position isomers

Toluenitriles have three position isomers. It is difficult to separate ortho and para isomers by C<sub>18</sub> or PE-MS column because of lack of poor π-π interaction. On the other hand, the isomers are well separated on PYE or NPE column which has strong π-π interaction.

### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II, 5PE-MS, 5NPE  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol / H<sub>2</sub>O = 40/60  
 5PYE Methanol / H<sub>2</sub>O = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: o-Toluenitrile (2.0 μg)  
 m-Toluenitrile (2.0 μg)  
 p-Toluenitrile (1.0 μg)



## 2) Selectivity for dipole

### Selectivity


Selectivity for dipole is evaluated based on the separation of 1,5-dinitronaphthalene and 1,8-dinitronaphthalene.

Dinitronaphthalenes (peak 1 and 2) were strongly retained on PYE and NPE because of  $\pi$ - $\pi$  interaction compared with dimethylnaphthalenes. However, there is a slight difference between these two columns. While 1,5-dinitronaphthalene (peak 2) was preferentially retained on PYE,

1,8-dinitronaphthalene (peak 1) was retained longer on NPE.



The results with NPE indicate the presence of strong dipole-dipole interaction. The two nitro group dipoles in

1,8-dinitronaphthalene are aligned for a much greater dipolar coupling with the bonded nitrophenyl group in NPE than 1,5-dinitronaphthalene.



### Selectivity for dipole

Column size: 4.6mm I.D.-150mm  
 Mobile phase: C<sub>18</sub>-MS-II Methanol / H<sub>2</sub>O = 80/20  
 NPE Methanol / H<sub>2</sub>O = 70/30  
 PYE Methanol / H<sub>2</sub>O = 90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm  
 Sample: 1; 1,8-Dinitronaphthalene (1,8-DNN) (0.21  $\mu$ g)  
 2; 1,5-Dinitronaphthalene (1,5-DNN) (0.11  $\mu$ g)  
 3; Naphthalene (0.25  $\mu$ g)  
 4; 1-Methylnaphthalene (0.35  $\mu$ g)  
 5; 1,5-Dimethylnaphthalene (0.42  $\mu$ g)



### Application data



#### ● Separation of phthalonitrile position isomers

Phthalonitriles have three position isomers. NPE or PYE completely separates these compounds due to  $\pi$ - $\pi$  interaction.

Furthermore, NPE strongly retains o-phthalonitrile due to dipole-dipole interaction.

### COSMOSIL Application Data

Column size: 4.6mm I.D.-150mm  
 Mobile phase: 5C<sub>18</sub>-MS-II Methanol / H<sub>2</sub>O = 30/70  
 5NPE Methanol / H<sub>2</sub>O = 40/60  
 5PYE Methanol / H<sub>2</sub>O = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm  
 Sample: o-; Phthalonitrile (0.3  $\mu$ g)  
 m-; Isophthalonitrile (3.0  $\mu$ g)  
 p-; Terephthalonitrile (0.15  $\mu$ g)




## 3) Selectivity for polyaromatic compounds

### Selectivity

Selectivity for polyaromatic compounds is evaluated based on the separation of benzene, naphthalene and anthracene.



The elution orders in all columns are the same : benzene, naphthalene and anthracene. Retention increases in all columns with increasing number of aromatic rings. In addition, highly dispersive packing materials such as PBB and PYE show much stronger retention for polyaromatic compounds due to dispersion interaction.



### Selectivity for polyaromatic compounds

Column size: 4.6mm I.D.-150mm  
 Mobile phase: 5C<sub>18</sub>-MS-II Methanol/H<sub>2</sub>O = 90/10  
 5PBB-R Methanol/H<sub>2</sub>O = 90/10  
 5PYE Methanol/H<sub>2</sub>O = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Benzene (1.67  $\mu$ g)  
 2; Naphthalene (0.11  $\mu$ g)  
 3; Anthracene (0.0063  $\mu$ g)



### Application data



#### ● Separation of dibenzosuberone and dibenzosuberone

C<sub>18</sub> retains dibenzosuberone (peak 1) longer than dibenzosuberone (peak 2). On the other hand, PBB-R and PYE retain dibenzosuberone (peak 2), which has a  $\pi$ -electron conjugated system, longer than dibenzosuberone (peak 1).

### COSMOSIL Application Data

Column size: 4.6mm I.D.-150mm  
 Mobile phase: 5C<sub>18</sub>-MS-II Methanol/H<sub>2</sub>O = 80/20  
 5PBB-R Methanol  
 5PYE Methanol/H<sub>2</sub>O = 90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm


Sample: 1; Dibenzosuberone (0.1  $\mu$ g)  
 2; Dibenzosuberone (0.025  $\mu$ g)



#### 4) Selectivity for molecular shape

##### Selectivity



Selectivity for molecular shape is evaluated based on the separation of chrysene and benz [a] anthracene. The isomers of two polyaromatic hydrocarbons, which consist of four benzene rings, are difficult to separate because of the similar hydrophobicity or aromaticity. However, PYE and Cholester columns, which recognize molecular shape, enable them to separate chrysene and benz [a] anthracene.



##### Selectivity for molecular shape

Column: 4.6mm I.D.-150mm  
 Mobile phase: 5C<sub>18</sub>-MS-II, 5PYE Methanol / H<sub>2</sub>O = 90/10  
 Cholester Methanol  
 5PE-MS Methanol / H<sub>2</sub>O = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Benz[a]anthracene (0.04μg)  
 2; Chrysene (0.04μg)



##### Application data



##### • Separation of diastereomers (threo- and erythro-)

C<sub>18</sub> cannot separate the threo and erythro forms. On the other hand, PYE retains the planar erythro form longer than the threo form.

##### COSMOSIL Application Data

Column: 4.6mm I.D.-150mm  
 Mobile phase: 5C<sub>18</sub>-MS-II, 5PYE Methanol / H<sub>2</sub>O = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm


Sample: 1; Threo form  
 2; Erythro form



#### 5) Selectivity for halides

##### Selectivity



Selectivity for halide is evaluated based on the separation of chlorobenzene and bromobenzene. PBB-R shows the highest selectivity factor due to dispersion interaction of the five bromine atoms.



##### Selectivity for halide

Column: 4.6mm I.D.-150mm  
 Mobile phase: Methanol / H<sub>2</sub>O = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Toluene (3.3μg)  
 2; Chlorobenzene (3.3μg)  
 3; Bromobenzene (3.3μg)



##### Application data


##### • Separation of halogen exchange reaction products

PYE and PBB-R retain dispersed iodine atom longer than bromine atom. As a result, PYE and PBB-R can separate the complicated bromine and iodine compounds that C<sub>18</sub> cannot separate.

##### COSMOSIL Application Data

Column: 4.6mm I.D.-150mm  
 Mobile phase: 5C<sub>18</sub>-MS-II Methanol / H<sub>2</sub>O = 90/10  
 5PBB-R, 5PYE Methanol  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: R=I or Br




Sample courtesy of Dr. H. Yamamoto, RIKEN, Condensed Molecular Materials Laboratory

## 6) Selectivity for hydrophobicity

### Selectivity


Selectivity for hydrophobicity is evaluated based on the separation of alkylbenzenes. Two C<sub>18</sub> and Cholester show similar high selectivity for hydrophobicity. Other columns show less hydrophobic selectivity than C<sub>18</sub>.



### Selectivity for hydrophobicity

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/H<sub>2</sub>O = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1: Benzene (1.67μg)  
 2: Toluene (1.67μg)  
 3: Ethylbenzene (1.67μg)  
 4: Propylbenzene (1.67μg)  
 5: Butylbenzene (1.67μg)  
 6: Amylbenzene (1.67μg)




Lower concentration of organic solvent in mobile phase leads to much retention in reversed phase chromatography. In case of NPE, when methanol concentration is reduced to 60%, the retention times increase to those similar to C<sub>18</sub> with 80% methanol.

### Control of retention time

Column: 5NPE  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase:  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1: Benzene (1.67μg)  
 2: Toluene (1.67μg)  
 3: Ethylbenzene (1.67μg)  
 4: Propylbenzene (1.67μg)  
 5: Butylbenzene (1.67μg)  
 6: Amylbenzene (1.67μg)



## 2. Preparation of mobile phase for HPLC

### 1) Organic solvent / aqueous mixed mobile phase

1)-1. Preparation of methanol : water = 70 : 30 (v/v) 1L

- ① Measure 700 ml of methanol in a measuring cylinder.
- ② Measure 300 ml of distilled water in a measuring cylinder.
- ③ Mix ① and ② thoroughly and degas.


The better approach is to prepare the mobile phase gravimetrically rather than volumetrically. Following is example of preparation.

Composition table for mobile phase 1L using methanol and water			Composition table for mobile phase 1L using acetonitrile and water		
Methanol / Water	Methanol (g)	Distilled water (g)	Acetonitrile / Water	Acetonitrile (g)	Distilled water (g)
90/10 (v/v)	711.9	99.8	90/10 (v/v)	707.4	99.8
80/20 (v/v)	632.8	199.6	80/20 (v/v)	628.8	199.6
70/30 (v/v)	553.7	299.5	70/30 (v/v)	550.2	299.5
60/40 (v/v)	474.6	399.3	60/40 (v/v)	471.6	399.3
50/50 (v/v)	395.5	499.1	50/50 (v/v)	393.0	499.1
40/60 (v/v)	316.4	598.9	40/60 (v/v)	314.4	598.9
30/70 (v/v)	237.3	698.7	30/70 (v/v)	235.8	698.7
20/80 (v/v)	158.2	798.6	20/80 (v/v)	157.2	798.6
10/90 (v/v)	79.1	898.4	10/90 (v/v)	78.6	898.4

Caution : Methanol and acetonitrile are hazardous substances, do not use for medical purpose. Always process in a laboratory hood and wear an eye protection and a mask.

Influence of organic solvent composition in mobile phase on the retention time.


1% difference in the composition significantly changes the retention.



Special attention should be paid to measure correct amount of organic solvent as retention time is significantly changed by 1% different composition.


Column: 5C<sub>18</sub>-MS-II 4.6mm I.D. x 150mm  
 Flow rate: 1.0ml/min  
 Detection: 254nm 0.16AUFS  
 Temperature: 30°C  
 Sample: 1. Uracil  
 2. Naphthalene

### Differences between acetonitrile and methanol in reversed phase liquid chromatography



Column: COSMOSIL 5C<sub>18</sub>-MS-II  
 Column size: 4.6 mm I.D. x 150 mm  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C

Back pressure differs depending on species of organic solvents and mixing ratio. Back pressure of acetonitrile is lower than that of methanol at the same concentration.



## 2) Organic solvent/ aqueous mixed mobile phase

### 2)-1. Preparation of methanol : 20 mmol/l phosphate buffer (pH2.5) = 80 : 20 (v/v) 1L

- Preparation of 20 mmol/l phosphate buffer (pH2.5)

- Prepare 20 mmol/l sodium dihydrogenphosphate aqueous solution.
- Prepare 20 mmol/l phosphoric acid aqueous solution.
- Adjust the pH to 2.5 by mixing ① with ②.
- Filter ③ under reduced pressure to remove insoluble substance that may deteriorate pump-seal and clog columns (0.45 μm or smaller pore size is recommended).

#### (Easy method)

- Dissolve 1.31 g of sodium dihydrogenphosphate and 1.05 g of phosphoric acid in distilled water to make 1 L solution.

- Filter the solution under reduced pressure to remove insoluble substance (0.45 μm or smaller pore size is recommended).

- Confirm that the solution is pH2.5.

- Preparation of methanol : 20 mmol/l phosphate buffer (pH2.5) = 80 : 20 1L

- Measure 800 ml of methanol in a measuring cylinder.
- Measure 200 ml of 20 mmol/l phosphate buffer (pH2.5) in a measuring cylinder.
- Mix ① and ② thoroughly and degas.

### 2)-2. Preparation of methanol : 20 mmol/l phosphate buffer (pH7.0) = 80 : 20 (v/v) 1L

- Preparation of 20 mmol/l phosphate buffer (pH7.0)

- Prepare 20 mmol/l sodium dihydrogenphosphate aqueous solution.
- Prepare 20 mmol/l di-sodium hydrogenphosphate aqueous solution.
- Adjust the pH7.0 by mixing ① with ②.
- Filter ③ under reduced pressure to remove insoluble substance that may deteriorate pump-seal and clog columns (0.45 μm or smaller pore size is recommended).

#### (Easy method)

- Dissolve 1.14 g of sodium dihydrogenphosphate and 1.49 g of di-sodium hydrogenphosphate in distilled water to make 1L solution.
- Filter the solution under reduced pressure to remove insoluble substance (0.45 μm or smaller pore size is recommended).
- Confirm that the solution is pH7.0.

- Preparation of Methanol : 20 mmol/l phosphate buffer (pH7.0) = 80 : 20 1 L

- Measure 800 ml of methanol in a measuring cylinder.
- Measure 200 ml of 20 mmol/l phosphate buffer (pH7.0) in a measuring cylinder.
- Mix ① and ② thoroughly and degas.

The better approach is to prepare the mobile phase gravimetrically rather than volumetrically. Following is example of preparation.

Methanol : 20 mmol/l Phosphate buffer	Methanol (g)	20 mmol/l Phosphate buffer (pH2.5) (g)	20 mmol/l Phosphate buffer (pH7.0) (g)
90 / 10 (v/v)	711.9	99.8	99.9
80 / 20 (v/v)	632.8	199.6	199.8
70 / 30 (v/v)	553.7	299.4	299.7
60 / 40 (v/v)	474.6	399.2	399.6
50 / 50 (v/v)	395.5	499.0	499.5
40 / 60 (v/v)	316.4	598.8	599.4
30 / 70 (v/v)	237.3	698.6	699.3
20 / 80 (v/v)	158.2	798.4	799.2
10 / 90 (v/v)	79.1	898.2	899.1

Caution : Methanol and acetonitrile are hazardous substances, do not use for medical purpose. Always process in a laboratory hood and wear an eye protection and a mask.

## 3) Preparation of ion pair reagent containing mobile phase

### 3)-1. Preparation of 5 mmol/l sodium 1-butanefulfonate containing 20 mmol/l phosphate buffer (pH2.5)

- Prepare 5 mmol/l sodium 1-butanefulfonate containing 20 mmol/l sodium dihydrogenphosphate aqueous solution.
- Prepare 5 mmol/l sodium 1-butanefulfonate containing 20 mmol/l phosphoric acid aqueous solution
- Adjust the pH to 2.5 by mixing ① with ②.
- Filter ③ under reduced pressure to remove insoluble substance that may deteriorate pump-seal and clog columns (0.45 μm or smaller pore size is recommended).

#### (Easy method)

- Dissolve 1.31 g of sodium dihydrogenphosphate, 1.05 g of phosphoric acid and 0.80 g of sodium 1-butanefulfonate in distilled water to make 1L solution.

- Filter the solution under reduced pressure to remove insoluble substance (0.45 μm or smaller pore size is recommended).

- Confirm that the solution is pH2.5.


※ 0.5 M Sodium 1-butanefulfonate aqueous solution is also available from Nacalai Tesque.

### 3. Sample pretreatment for HPLC

Pretreatment before HPLC analysis is often required for samples of low concentration or samples containing analytical contaminants. It improves reproducibility and sensitivity in analysis, and protects HPLC columns. The preparation methods are different according to the each sample. The followings are examples of different pretreatments.

#### 1) Filtration


Filtration is a common method used for separating solids from liquids. It extends a column's life by minimizing column damages from solid contaminants such as particles, sediments and colloid substances. It also improves reproducibility of analytical data. We offer both syringe-type and spin-type filters for sample filtration.

	■ Syringe filter	■ Centrifugal filter
Product	Cosmonice filter	Cosmospin filter
Configuration		
Usage	Easy to use. Just attach a filter on top of a syringe.	Easy to use by centrifugation.
Type	<ul style="list-style-type: none"> <li>• W (aqueous system)</li> <li>• S (solvent system)</li> </ul>	<ul style="list-style-type: none"> <li>• Pore diameter : 0.2 μm</li> <li>• Pore diameter : 0.45 μm</li> </ul>
Required equipment	Syringe · Sample bottle	Centrifuge
Page	refer to page 74	refer to page 74

#### Cosmonice filter

How to use :

- ① Fill a syringe with the sample you want to filter.
- ② Attach a Cosmonice filter to the syringe.
- ③ Push the syringe plunger to filter the sample.
- ④ Analyze the filtered sample by HPLC.




#### Cosmospin filter

Components : · Sample reservoir  
· Sample collection tube


How to use :

- ① Insert a Cosmospin sample reservoir into a Cosmospin sample collection tube.
- ② Add a sample into the Cosmospin sample reservoir.
- ③ Close the sample collection tube cap and centrifuge.
- ④ Remove the sample reservoir and collect the filtered sample in the sample collection tube.
- ⑤ Analyze the filtered sample by HPLC.



### 2) Protein precipitation

Protein precipitation is commonly used to remove proteins in samples for downstream analysis. For example, when analyzing drug concentration in blood samples, proteins have to be removed first. Otherwise, proteins may be adsorbed in columns and interfere with the analysis. Common methods for protein precipitation include salting out, isoelectric point precipitation and precipitation with organic solvents. The following shows a general procedure for protein precipitation with organic solvents.




Procedure for protein precipitation :

### 3) Ultrafiltration


Ultrafiltration is a method to concentrate proteins or other macromolecules through a semipermeable membrane with defined pores. Ultrafiltration is applicable for sample desalting, concentrating proteins from dilute solution such as urine samples, or deproteinizing samples with high protein concentration (i.e. blood serum or plasma). Following is a general procedure for ultrafiltration.

Procedure for ultrafiltration :




#### 4) Solvent extraction method

Solvent extraction is a method to separate compounds due to their unequal solubility in two immiscible liquid phases, usually water and an organic solvent. The method is used to concentrate highly hydrophobic compounds, and consequently increase analytical sensitivity. A buffer solution is added to sample to optimize the pH and target substance is then extracted by an organic solvent such as ether and chloroform. However, when target substance is combined with proteins, solvent extraction may not work well.




#### 5) Ion exchange

Pretreatment by ion-exchange resin may be effective for samples that the solvent extraction method cannot be adapted due to its emulsification. A preliminary experiment may be required for the selection of resin and experimental conditions. For example, a negatively charged compound is strongly adsorbed on an anion-exchange resin such as DEAE cellulose resin. Therefore, the target compound is collected by increasing salt concentration of buffer solution or adjusting pH of elution buffer after washing off other weakly adsorbed undesired substances.



## 4. Baseline noise in gradient elution

In gradient analysis, incomplete mixing of mobile phases or impurities in water of mobile phase can cause baseline noise. In the former case, it can be improved by using a proper mixer before injector (Baseline 1→2). In the latter case, it can be improved by using a pre-column. Impurities in water are adsorbed on the pre-column (Baseline 2→3). COSMOSIL 5C<sub>18</sub>-AR-II 4.6 mm I.D. x 10 mm or 10 mm I.D. x 20 mm as a pre-column.




Column COSMOSIL 5C<sub>18</sub>-AR-300 4.6 mm I.D. x 150mm  
 Precolumn COSMOSIL 5C<sub>18</sub>-AR-II 4.6 mm I.D. x 10mm  
 Mobile phase A: 0.1% TFA containing water  
 B: 0.1% TFA containing 95% acetonitrile  
 B: 0% → 100%/30 min liner gradient  
 Flow rate 1.0 ml/min  
 Temperature 30°C  
 Detection UV 220nm

## 5. Troubleshooting for increased pressure

Repeated analysis may increase back pressure. Continuous use of HPLC columns under high pressure can cause deterioration and overload of the equipment. Therefore, it is important to monitor column back pressure regularly and solve the problem timely.

The back pressure increase can be due to clogging of a column or clogging of the equipment. First of all, identify the clogging site.

Schematic diagram of HPLC



Remove analytical column first and connect the plumbing from the back of HPLC system directly to the detector. Measure the pressure of flowing mobile phase without an HPLC column. Generally equipment should hardly generate any pressure. If there is significant flow pressure, disconnect the system components one by one to identify the clogged component(s). Possible causes and solutions of clogged equipment are discussed in section II below.

If the flow pressure without a column is normal, then pressure increase is due to clogging of a column. In this case, one needs to determine the causes and whether it is time to replace the column. Possible causes and solutions of clogged column are discussed in greater details in sections I and III.

Symptom	Possible Cause
Pressure increase rapidly in short-term use	→ Flow pressure without a column is 0-0.3 MPa → Clogging of column Refer to section I
	Yes Flow pressure with a column is 0-0.3 MPa or higher → Clogging of equipment Refer to section II
Pressure increase gradually in long-term use	→ Deterioration of column due to long-term use refer to section III

### I. Solution in case an HPLC column is clogged in short-term use.

Select the possible cause of clogging according to the following flow chart.

STEP 1	• Salt deposition • Use mobile phase of high concentration organic solvent right after using buffer	→ YES	Cause 1
↓ NO			
STEP 2	• Forget to filter mobile phase • Sample is not dissolved enough	→ YES	Cause 2
↓ NO			
STEP 3	• Analyzing samples which tend to adsorb to a column (i.e. protein samples) • Sample deposition in column	→ YES	Cause 3

**Cause 1** Salt is deposited on a column.

**Solution :** Wash columns for 30 minutes at half of the analytical flow rate with 10% organic solvent (methanol or acetonitrile) in water to dissolve deposited salt. If the situation is not improved, wash with 100% water under the same condition.

**Prevention :** To switch to high concentration organic solvent after using a buffer, first wash a column with a mobile phase not containing salt (with the same concentration of organic solvent as the buffer), then switch to the mobile phase of higher organic solvent concentration. Example : To change mobile phase from 10/90 (v/v) acetonitrile/20mmol/l phosphate buffer (pH2.5) to 90/10 (v/v) acetonitrile/water, first wash for 15 minutes with 10/90 (v/v) acetonitrile/water, and then switch to 90/10 (v/v) acetonitrile/water.

**Cause 2** Column filter is clogged by sample or impurities.

**Solution :** Connect the column in reverse direction, and then wash the column for 30 minutes at half of the usual analytical flow rate with the mobile phase used for analysis. If the situation is not improved, change the end fitting in the front of column. (We can replace end fittings with a paid service fee.)

**Prevention :** We recommend filtering sample and/or mobile phase. For more information, please refer to page 164, TECHNICAL NOTE 3. Sample pretreatment for HPLC 1) filtration.

**Cause 3** Sample may be adsorbed to packing material or deposited in a column.

**Solution :** Wash for 30 minutes at half of analytical flow rates with a solvent which adsorbed substances are dissolved in. The followings are how to wash each type of columns.

[Reversed phase columns]

a) When an adsorbed substance is not protein, wash with methanol or tetrahydrofuran.

b) When an adsorbed substance is protein, wash with 50-70% of acetonitrile/water (containing 0.1% of trifluoroacetic acid). However proteins may be deposited in high concentration of organic solvent depending on varieties.

[COSMOSIL Sugar-D/NH<sub>2</sub>/HILIC columns]

Wash with 50/50 (v/v) acetonitrile/water for NH<sub>2</sub>-MS and 100% water for Sugar-D and HILIC columns.

[COSMOSIL SL-II]

Wash with methanol, tetrahydrofuran or ethanol.

**Prevention :** Choose appropriate pretreatment for each sample. For more information, please refer to page 164, TECHNICAL NOTE 3. Sample pretreatment for HPLC 1) filtration. We also recommend using guard column. For more information for guard columns, please refer to page 172, TECHNICAL NOTE 6. Effect of guard columns.

**Caution :**

- When wash columns, do not connect column exit and let the solution through.
- Long term of washing may deteriorate the performance of columns.
- Do not use strongly alkaline solution (more than pH 7.5) or strongly acidic solution (less than pH 1.5) for silica gel base packing material.
- Store columns with manufacturer recommended storage solvent after washing, When the situation is not improved, replace the column.

**II. Solutions in case pressure is too high because of clogged equipment.**

First, identify the specific clogging site by disconnecting the components in the system one by one and checking the flow pressure. The followings are possible common causes.

**Cause 1** Salt is deposited in plumbing.

**Solution :** Flow water to the plumbing without connecting a column and any other equipment. Washing out the plumbing in a reversing connection is also an effective way. If the situation is not improved, replace it with a new one.

**Cause 2** Check-valve of pump is clogged by stain

**Solution :** Wash the check-valve with a stain dissolving solvent. Take apart the washable part, soak it in the solvent, then clean in an ultrasonic cleaner.

**Cause 3** Manual injector is clogged with stain

**Solution :** Wash with a stain dissolving solvent. Soak rotor seal and line filter in water and clean them in an ultrasonic cleaner. If the situation is not improved, replace the injector.

**Prevention :** It extends the life time of an HPLC system to maintain regular wash of the system. Wash the system the same as wash an HPLC column. When the mobile phase contains salt, wash for 10-15 minutes with a mobile phase which has the same composition but not containing salt. For example, when using 50/50 (v/v) methanol/20mmol/l phosphate buffer, wash with 50/50 (v/v) methanol/water. When the mobile phase contains halogen, acid and/or base, wash for 10-15 minutes with mobile phase which has the same composition but not containing halogen, acid and/or base.

**III. Solutions in case a column is damaged from long term use**

Every column will have to be replaced eventually. Performance of a column is expected to deteriorate slowly after long term use. One has to decide whether it is time to replace the column.

**Cause 1** Column deterioration result from long term use.

**Solution :** Wash according to Section I, Cause 3.

**Prevention :** Same as prevention in Section I, Cause 3. When the column condition is not improved, you could continue to use the column if peak shapes do not change and the maximum pressure is less than 20 MPa. However, we recommend replacing the column because it place extra burden on the equipment.

**Cause 2**

Silica gel in the column may be cracked because of long term use.


**Solution :** Replace the columns.

## 6. Effect of guard column

The use of guard columns to protect both analytical and preparative columns is highly recommended. COSMOSIL guard columns are packed with packing materials identical to that used in analytical and preparative columns. As a result, the COSMOSIL guard columns do not contribute to any decrease in the performance of the main column.


### Example of using guard columns

The following chromatograms show analysis examples using a COSMOSIL 5C<sub>18</sub>-MS-II analytical column (4.6 mm I.D. × 150 mm) and the same column connected with its guard column (4.6 mm I.D. × 10 mm). There is no change in separation characteristics since the packing material of the guard column is identical to that used in the main column.



### Inspection of columns

We recommend that the performance of a column be examined regularly. Deteriorated columns shall be replaced timely. If a deteriorated guard column continues to be used, the packed column will also deteriorate.



### Ordering information

Please refer to respective pages of each column.

Please refer to page 76 for the parts necessary to use guard columns or guard column cartridges.

## 7. Troubleshooting for normal phase chromatography

**Q 1** : How can I convert from reversed phase mode to normal phase mode or vice versus using the same HPLC equipment?

**A 1** : To convert from reversed phase mode to normal phase mode, or vice versus, flush the equipment with a solvent that is miscible with both the current mobile phase and the intended mobile phase. Connect the HPLC pump directly with the detector, and replace the solvents according to following instructions.

Solvent conversion from reversed phase to normal phase

① To convert a mobile phase without buffer solution in reversed phase to normal phase, replace solvents according to the following steps :

- 1) Flush the equipment with a solvent for reversed phase. For example, methanol/H<sub>2</sub>O (v/v=50/50).
- 2) Flush the equipment with a solvent miscible to both mobile phases. For example, tetrahydrofuran, ethanol.
- 3) Flush the equipment with a solvent for normal phase.

② To convert a mobile phase with buffer solution in reversed phase to normal phase, replace solvents according to the following steps :

- 1) Flush the equipment with a solvent with buffer solution for reversed phase.
- 2) Flush the equipment with a solvent with the composition same as 1) and without salt.
- 3) Flush the equipment with a solvent miscible to both mobile phases. For example, tetrahydrofuran, ethanol.
- 4) Flush the equipment with a solvent for normal phase.

Solvent conversion from normal phase to reversed phase

① To convert from normal phase to a mobile phase without buffer solution in reversed phase, replace solvents according to following steps :

- 1) Flush the equipment with a solvent for normal phase.
- 2) Flush the equipment with a solvent miscible to both mobile phases. For example tetrahydrofuran, ethanol.
- 3) Flush the equipment with a solvent for reversed phase.

② To convert from normal phase to a mobile phase with buffer solution in reversed phase, replace solvents according to following steps :

- 1) Flush the equipment with a solvent for normal phase. For example, hexane/ethyl acetate.
- 2) Flush the equipment with a solvent miscible to both mobile phases. For example tetrahydrofuran, ethanol.
- 3) Flush the equipment with a solvent with the composition same as 4) and without salt.
- 4) Flush the equipment with a solvent for reversed phase.

**Q 2** : My flow rate is not stable. How can I troubleshoot?

**A 2** : Possible causes for unstable flow rate can be a malfunctioning check valve or air in a mobile phase, wash the check valve thoroughly by ultrasonic cleaner. Solvents with a low boiling point such as n-Hexan and n-Heptane generate air easily. To prevent air generation, degas the mobile phase sufficiently.

**Q 3** : In spite of using the same condition, the retention time is different. How can I solve the problem?

**A 3** : One possible cause is unstable flow rate. Please refer to Q2 section.

Another possible cause is variation of polar component in mobile phases. In normal phase chromatography, the retention time depends on the concentration of small amounts of very polar constituents in the mobile phase. This is especially true for water content in a mobile phase. In this case, always use fresh solvents in a mobile phase. If sample solvent includes water, change to a solvent without water or decrease injection volume. If a column contains water, remove water from the column by washing it with ethanol.

Q 4 : How can I wash the COSMOSIL SL-II column?

A 4 : The SL-II column can be washed with tetrahydrofuran, methanol, ethanol, methylene chloride, n-Hexane or n-Heptane.

Q 5 : How can I store the COSMOSIL SL-II column?

A 5 : Store the SL- II column with the shipping screw tighten in order to prevent the solvent in the column from volatilization. In case where a solvent containing halogens is used, replace the solvent in the column with a solvent without halogens such as n-Heptane before storing.

Q 6 : I get peak tailing in my run. What can I do about it?

A 6 : ① In case where the sample contains acidic compounds, add approx. 0.5% of acetic acid to the mobilephase.  
② In case where the sample contains basic compounds, add approx. 0.5% of triethylamine to the mobile phase.

Q 7 : I get no peaks. How can I troubleshoot?

A 7 : First, make sure that there is no problem with the system. If the problems are with the sample, sample solvent or mobile phase, try following check list.  
1. The analyte may not be eluted from a column because the retention of the analyte is too strong. In this case, use a stronger eluent (mobile phase).  
2. The sample contains chelating compounds or basic compounds. They may be adsorbed to the packing materials. In this case, add 0.1% - 1% acid (trifluoroacetic acid or acetic acid) to the mobile phase.





## 8. Inner diameter of column (scale down and scale up)

The figure below shows general parameters for 1.0 mm to 50 mm I.D. COSMOSIL columns : flow rate, equipment, inner diameter of pipe, application, surface ratio (compared with 4.6 mm I.D.) and particle size. It may help to scale up or down from the most commonly used 4.6 mm I.D. column.

Inner diameter (mm I.D.)	1.0	2.0	3.0	4.6	10	20	28	50
Flow rate (ml/min)	0.05	0.2	0.4	1.0	5.0	19	37	70
Detector cell · Injector	for Semi-micro		for Analytical			for Preparative		
Inner diameter of pipe (mm)	0.05	0.1	0.2-0.3			1.0		
Application	LC-MS Solvent saving		Solvent saving with standard system	Standard	Preparative (small scale)	Preparative (medium scale)	Preparative (large scale)	Preparative (super large scale)
Surface ratio with 4.6 mm I.D.	0.05	0.19	0.43	1.00	4.73	18.90	37.05	118.15
Particle size (μm)	3 or 5			5		15 or more		


### Scale down

When scaling down from the most commonly used analytical column (4.6 mm I.D.) to a semi-micro or 3.0 mm I.D. analytical HPLC column (of the same column length), sample loading dose is proportionate to the cross section of column. The 3.0 mm I.D. columns provide high sensitivity and solvent saving without the need to change the existing equipment settings. Semi-micro columns (2.0 mm I.D. and 1.0 mm I.D.) provide higher sensitivity and enable analysis of minor components, but one needs to change the piping of HPLC equipment, the injector and the detector cell for semimicro columns.

Column size	4.6 mm I.D. × 150 mm	3.0 mm I.D. × 150 mm	2.0 mm I.D. × 150 mm	1.0 mm I.D. × 150 mm
Chromatogram				
Flow rate (ml/min)	1.0	0.4	0.2	0.05
Pressure (MPa)	3.4	3.6	3.8	3.6
Injection volume(μl)	1.0	0.4	0.2	0.05
Detector Cell · Injector	for Analytical		for Semi-micro	
Detector sensitivity(AUFS)	0.08		0.04	
Inner diameter of pipe (mm)	0.25		0.10	0.05
	Column Mobile phase Flow rate Temperature Detection	COSMOSIL 5C <sub>18</sub> -MS- II Acetonitrile : Water = 70 : 30 1.0 ml/min 30°C UV 254 nm		Sample 1. Benzene 2. Toluene 3. Ethylbenzene 4. Propylbenzene 5. Butylbenzene 6. Amylbenzene

# V INDEX

Sample name .....	178
Packing material's name .....	188



## Sample name

A	
Acenaphthene	87, 148
Acenaphthylene	148
Acesulfame K	123
Acetaminophen	31, 80, 110
Acetanilide	90, 99, 110
Acetic Acid	18
Acetoaminophen	15
Acetophenone	14, 30, 147
4'-Acetoxyacetanilide	31, 80
17 $\alpha$ -Acetoxyprogesterone	86
Acettrizoic Acid	81
N-Acetylprocainamide	14, 110
Acetylsalicylic Acid (Aspirine)	110
Acid Red	123
Aconitine	120
Acrinol	80
Acrylic Acid	73
Actinase E	139
Adenine	18, 135
Adenosine	18, 100, 135
5'-ADP	135
L-Adrenaline	18, 72
[Ala2]-Leu-Enkephalin	138
[Ala2]-Met-Enkephalin	138
Albiflorin	115
Albumin (Bovine)	44, 53, 139
Albumin (Goat)	140
Albumin (Human)	139
Albumin (Ovalbumin)	140
Albumin (Rabbit)	140
Alcohol Dehydrogenase	139
Alizalin	119
Allantoic Acid	112
Allantoin	112
Allopurinol	129
Allylbenzene	149
Amaranth	123
Amidotrizoic Acid	81
Amikacin Sulfate 2,4,6-trinitrobenzenesufonic Acid Derivative	81
4-Aminoacetophenone	90, 98
Aminoantipyrine	110
4-Aminoantipyrine	15
4-Aminobenzoic Acid	80
4-Amino-n-butyric Acid [GABA]	137
4-Amino-6-chlorobenzene-1,3-disulfonamide	98
4-(Aminomethyl) benzoic Acid	94
m-Aminophenol	147
o-Aminophenol	147
p-Aminophenol	147
p-Aminophenol Hydrochloride	31, 80
Amlodipine Besilate	81
Ammelide	39
Ammeline	39
Amoxicillin	81
5'-AMP	135
dAMP	135
Ampicillin	99, 111
Ampicillin Sodium Salt	81
D-(-)-Amygdalin	114
$\alpha$ -Amylase	139
Amylbenzene	12, 21, 34, 119, 160
4-Androstene-3,17-dione	128
Angiotensin I (Human)	16, 137
Angiotensin II	138
Angiotensin II (Human)	16, 137, 138
Angiotensin II, [Asn1,Val5]	138
Angiotensin II, [Sar1,Ala8]	138
Angiotensin II, [Sar1,Ile8]	138
Angiotensin II, [Sar1,Thr8]	138
Angiotensin II, [Val5]	138
Angiotensin II, Des-Asp1-[Ile8]	138
Aniline	147
Anisole	147, 155
Anthracene	22, 148, 157
Anthrurufin	119
Antipyrine	15, 110
Apigenin	118
Aprotinin	48
Arabinose	41
Arbutin	113
Ascorbic Acid	19, 125, 126
L-(+)-Ascorbic Acid	125
Aspartame	123
Aspartic Acid	137
Asulam	121
Atenolol	81
5'-ATP	135
Atrazine	122
Atropine Sulfate	15, 116, 117
Azathioprine	80
Aztreonam	80
B	
Baclofen	97
Baicalein	23
Baicalin	113, 115
Bamethane Sulfate	102

Barbital	110
Barbital Sodium	117
Beclometasone Dipropionate	86, 102, 103
Bensulide	121
Benzamide	147
Benz[a]anthracene	148, 158
Benzene	14, 21, 30, 31, 34, 144, 146, 147, 155, 157, 160
Benzo[b]fluoranthene	148
Benzoic Acid	36, 80, 106, 123, 146
Benzonitrile	147
Benzo[g,h,i]perylene	148
Benzophenone	86, 88, 106, 107
Benzo[a]pyrene	148
Benz[a]pyrene	148
Benzyl Alcohol	9, 14, 106
Benzyl p-Hydroxybenzoate	98
Benzylpenicillin Potassium	103
2-Benzylpyridine	148
3-Benzylpyridine	148
4-Benzylpyridine	148
Berberine Chloride	103, 113
Berberine Hydrochloride	25
Betahistine Mesilate	103
Betamethasone	103, 129
Betamethasone Sodium Phosphate	103
Betamethasone Valerate	103
Betamethasone 17-Valerate	172
Bezafibrate	102
Bibenzyl	22
1,1'-Binaphthyl	149
D-Biotin [Vitamin H]	125
Biphenyl	97
Bisoprolol	98
Boc-(D)Phe-(L)Phe-OMe	137
Boc-(L)Phe-(L)Phe-OMe	137
Bombesin	138
Bovine Serum	142
BPMC	120
Bradykinin	138
Brilliant Blue FCF	123
Bromacil	122
Bromhexine Hydrochloride	102
Bromobenzene	144, 147, 159
Brucine Dihydrate	116, 117
Bufalin	115
Bufexamac	100
Butylbenzene	12, 21, 34, 160
n-Butylbenzene	147
Butylbenzoate	12
2-tert-Butylhydroquinone	124
3-tert-Butyl-4-hydroxyanisol	124
Butyl p-Hydroxybenzoate	16, 82, 84, 85, 86, 88, 98, 101, 103, 108, 117
iso-Butyl p-Hydroxybenzoate	108
Butyl Salicylate	16
C	
C60	54, 56, 58, 59
C60O	56
C70	54, 56, 58
C76	56, 58
C84	56, 58
C86	57
Caffeine	23, 31, 35, 93, 96, 101, 105, 107, 109, 115, 125
Calcium Folate	104
Camostat Mesilate	84
Campesterol	26
Capsaicin	27, 116, 120
Captopril	83
Carbamazepine	109, 110
Carbazochrome	84
Carbazochrome Sodium Sulfonate	84
Carbazole	33
Carbenicillin	111
Carbidopa	84
Carbonic Anhydrase	139, 141
Carbosulfan	122
Carnitine	125
$\alpha$ -Carotene	118
$\beta$ -Carotene	118, 123, 127
Carvylacetate (cis, trans)	148
CAT	122
Catalase	140
Catechin (C)	23, 35
Catechin gallate(CG)	23, 35
Cefaclor	90
Cefadroxil	90
Cefsulodin Sodium Salt Hydrate	90
Ceftazidime	90, 91
Ceftriaxone Sodium	91
Cephaeline Hydrobromide	116
Cephalothin Sodium Salt	90
Cetirizine Dihydrochloride	89
Chlomethoxnyl	121
Chlormadinone Acetate	86
Chlorobenzene	144, 147, 159
4-Chlorobenzoic Acid	82
p-Chlorobenzoic Acid	102
4-Chlorobenzophenone	86
Chloro IPC	122
m-Chlorophenol	143
o-Chlorophenol	143
p-Chlorophenol(4-)	85, 86, 93
4-(4-Chlorophenyl)-4-hydroxypiperidine	97
Chlorotoluron	122
Chloroxuron	122

Chlorpheniramine Maleate	86, 95
Chlorpropamide	86
Chlortetracycline	111
Cholesterol	26
Choline Oxidase	139
Chondroitinase AC-II digested	132
Chondroitin Sulfate A	132
Chondroitin Sulfate C	132
Chrysazin	119
Chrysene	148, 158
Chrysin	23
$\alpha$ -Chymotrypsinogen A	140, 141
Ciclosporin (Cyclosporin)	88
Cilastatin Sodium Salt	89
Cimetidine	110
Cinchonidine	84
trans-Cinnamic Acid	117
Cinnamyl Alcohol	36
Cinobufagin	115
Citric Acid	18, 125
L-Citrulline	120
Clarithromycin	85
Clindamycin Hydrochloride	85
Clindamycin Phosphate	85
Clobetasol Propionate	86
Clofibrate	85, 86
Cloperastine Hydrochloride	86
Cloxacillin Sodium Monohydrate	85
CNP	121
Conalbumin	48, 140
Concanavalin A	140, 141
Coomassie Brilliant Blue G-250	149
Cortisone	24, 87, 129
Cortisone Acetate	87, 101, 129
Coumarin	119
Creatinine	128
m-Cresol	106
Crotonic Acid	73
Crude fullerenes	58
Cyanazine	122
Cyanocobalamin	87, 105
Cyanocobalamin Derivative	87
Cyanuric Acid	39
$\alpha$ -Cyclodextrin	130
$\beta$ -Cyclodextrin	130
$\gamma$ -Cyclodextrin	130
Cytidine	18, 135
Cytochrome C	44, 138, 140
Cytosine	18, 135
<b>D</b>	
Daidzein	25

Daidzin	25
Dansyl Amino Acids	136
Dantrolene Sodium	91
dATP	18
Daunorubicin Hydrochloride	83, 91
dCMP	135
DCMU	121
DCEPA	121
dCTP	18
Deferoxamine Mesilate	93
Dehydroacetic Acid	123
Deoxyribonuclease I	140
2'-Deoxyuridine	82
Desipramine Hydrochloride	109
Dexamethasone	92, 129
dGMP	135
dGTP	18
Diaphorase	141
Dibenz[a,h]anthracene	148
Dibenzosuberone	149, 157
Dibenzosuberone	149, 157
Dibucaine Hydrochloride	92
2,6-Di-tert-butyl-4-hydroxytoluene	124
Di-n-butyl Phthalate	96
Diclofenac Sodium	88
Diclofenamide	88
Diethylene Glycol	145
Diethyl Terephthalate	91
1,2-Difluorobenzene	144
1,3-Difluorobenzene	144
1,4-Difluorobenzene	144
Digitoxigenin	111
Digitoxin	87, 111
Digoxin	87
Dihydrocoumarin	119
Dihydroquinine	84
7,8-Dihydroxyflavone	23
3,4-Dihydroxyphenylacetic Acid	128
Diketopiperazine	123
Diltiazem Hydrochloride	89
Dimefuron	122
4-Dimethylaminoantipyrine	89
Dimethylaminophenol	97
N,N-Dimethylaniline	146, 147
N,N-Dimethylbenzylamine	9
1,5-Dimethylnaphthalene	28, 156
2,4-Dinitroaniline	147
1,2-Dinitrobenzene	83
1,5-Dinitronaphthalene(1,5-DNN)	28, 156
1,8-Dinitronaphthalene(1,8-DNN)	28, 156
Diphenhydramine Hydrochloride	110
Diphenyl	124
Diphenylimidazole	100

Diphenylmethane	149
Diphenyl Phthalate	95
Dipyridamole	88
Disulfiram	88
DNA	135
DNA-wrapped CNT	60
Dobutamine Hydrochloride	94
Domperidone	95
L-DOPA	18, 72, 128
Dopamine	128
Dopamine Hydrochloride	18, 72, 94
Doxifluridine	93
Doxorubicin Hydrochloride	83, 93
Doxycycline Hydrochloride	93
dTTP	18
<b>E</b>	
E.coli proteins	49
Egg White	142
Elaidic Acid	133
Emetine Hydrochloride	116
Epicatechin(EC)	23, 35
Epicatechin gallate(ECG)	23, 35
Epigallocatechin(EGC)	23, 35
Epigallocatechin gallate(EGCG)	23, 35
Epinephrine	128
Epirubicin Hydrochloride	83
meso-Erythritol	41, 130, 131, 146
Erythrosine	123
Estradiol	128
17- $\alpha$ -Estradiol	128
17- $\beta$ -Estradiol	128
Estradiol Benzoate	83
Estriol	83, 128
Estrone	128
Ethidimuron	122
Ethinylestradiol	109
Ethiofencarb	122
Ethofumesate	122
o-Ethoxybenzamide	95
4-Ethoxyphenol	85
Ethylbenzene	21, 34, 147, 160
p-Ethylbenzoic Acid	146
Ethylene Glycol	145
Ethyl p-Hydroxybenzoate(4-)	16, 84, 87, 88, 89, 95, 100, 102, 147
2-Ethylpyridine	146
Ethyl Salicylate	16, 147
Ethyl Vanillin	104

<b>F</b>	
FAD	100
Famotidine	99, 110
Fast Green FCF	123
Fibrinogen	140
Fisetin	23
Flavin Mononucleotide Sodium Salt	125
Flavone	23
Flopropione	102
Flumequine	112
Fluocinolone Acetonide	100, 101, 129
Fluocinonide	100
Fluoranthene	148
Fluorene	148, 149
Fluorobenzene	144, 147
Fluorodifen	121
Fluorometholone	101, 129
Fluoxymesterone	100
Flurbiprofen	101, 111
FMN	100
Folic Acid	104, 105, 125
Fructooligosaccharides	42
Fructose	41, 42, 130, 131
D-Fructose-6-phosphate	132
Fullerene toluene extract	55
Fumaric Acid	18, 98
Furan	147
Furfural	147
Furfuryl Alcohol	147
<b>G</b>	
Gabexate Mesilate	84
Galactose	41
Gallic Acid	16, 113
Gallocatechin(GC)	23, 35
Gallocatechin gallate(GCG)	23, 35
Gd@C82(I)	57
Geniposide	113, 115
Genistein	25
Genistin	25
Gentisic Acid	16
6-Gingerol	120
Ginkgolide A	120
Ginkgolide B	120
Ginkgolide C	120
Ginsenoside Rb1	114
Ginsenoside Rc	114
Ginsenoside Re	114
Ginsenoside Rg1	114
Gitoxin	87
Glucitol (Sorbitol)	41, 131

Glucose	40, 41, 42, 130, 131, 132
α-D-Glucose-1-phosphate	132
D-Glucose-6-phosphate	132
Glucose Oxidase	48, 139
β-Glucosidase	53
Glucuronic Acid	42
Glutamic Acid	137
L-Glutamic Dehydrogenase	139
Glyceric Acid	146
Glycerol	41, 145
Glycine	38
Glycitein	25
Glycitin	25
Glycolic Acid	18
Glycylglycine	38
Glycyrrhizic Acid	113, 124
Griseofulvin	85
Guaiacol Glycerol Ether	81, 85
Guanine	18, 135
Guanosine	18, 111, 135
<b>H</b>	
Haloperidol	97
Hemoglobin	53
Hemoglobin, Bovine	46
Hesperidin	117
Hibenzic Acid	92
Hinokitiol	119
Hippuric Acid	128
Histidine	73
Homovanillic Acid	18
Honokiol	114
Human Serum	142
Hyaluronic Acid	132
Hydralazine Hydrochloride	111
Hydrochlorothiazide	98
Hydrocortisone	24, 98, 129
Hydrocortisone Acetate	87, 98
Hydrocortisone Succinate	98
Hydroquinidine Hydrochloride	84
Hydroquinone	113
Hydroxocobalamin Acetate	105
p-Hydroxybenzoic Acid	87, 95
3-Hydroxyflavone	23, 119
5-Hydroxyflavone	23, 119
6-Hydroxyflavone	23, 119
7-Hydroxyflavone	23, 119
Hypoxanthine	129
<b>I</b>	
Ibuprofen	86, 111

Iodoxuridine	82
Imazalil	124
Imipramine Hydrochloride	109
Indeno[1,2,3-c,d]pyrene	148
Indigo Carmine	123
Indometacin (Indomethacin)	82, 115
Inositol	41, 131
Insulin	151
Iodobenzene	147
5-Iodouracil	82
Ipratropium Bromide	82
Ipratropium Bromide Derivative	82
Iprodione	121
Isoamyl Benzoate	172
Isoascorbic Acid	125
Isobutyl p-Hydroxybenzoate	81, 105
Isoleucine	39
Isomaltooligosaccharides	130
Isoniazid	82
Isonicotinamide	90
Isonicotinic Acid	82
Isophthalonitrile	26, 156
Isoproc carb	122
Isopropyl p-Hydroxybenzoate	91, 100
Isoproturon	122
Isoxsuprine Hydrochloride	82
<b>J</b>	
Josamycin	89
<b>K</b>	
Kanamycin Sulfate 2,4,6-trinitrobenzenesulfonic Acid Derivative	81
Ketamine Hydrochloride	109
Ketoprofen	87, 111
<b>L</b>	
Lactic Acid	18
β-Lactoglobulin	53
Lactose	43, 131
Latamoxef Sodium (isomer 1)	106, 107
Latamoxef Sodium (isomer 2)	106, 107
Leu-Enkephalin	138
Levallorphan Tartrate	108
Levothyroxine Sodium	109
Lidocaine	107, 109
Lincomycin Hydrochloride	85, 108
Linolenic Acid p-Bromophenacyl Ester	134
γ-Linolenic Acid p-Bromophenacyl Ester	134
Linuron	121, 122
Liothyronine Sodium	107

Lisinopril	107
Luteolin	23
Lypoxytase	51
Lysozyme	44, 140
Lysyl Endopeptidase digested	138

**M**

Magnolol	114
Maleic Acid	18, 86, 92
Malic Acid	18, 120, 126
Maltitol	41, 130, 131
Maltoheptaose	42
Maltohexaose	42
Maltopentaose	42
Maltose	40, 41, 42, 43, 130, 131
Maltotetraitol	131
Maltotetraose	42
Maltotriose	42, 131
Mandelic Acid	128
Maneb	121
Mannitol	131
Mannosamine	42
Mannose	42
MCC	121
MCP	121
MCPB	121
MCPP	121
Mecobalamin	105
Melamine	39
Menatetrenone	106
Mepivacaine Hydrochloride	106
5-Mercapto-1-methyltetrazole	106
6-Mercaptopurine Hydrate	80
Meropenem Trihydrate	106
Met-Enkephalin	138
Metabromuron	122
Metazachlor	122
Metformin Hydrochloride	105
Methabenzthiazuron	122
Methicillin Sodium Salt	111
Methomyl	122
Methotrexate	105
DL-3-Methoxy-4-hydroxymandelic Acid	18
6-Methoxyflavone	23
m-Methoxyphenol	143
o-Methoxyphenol	143
p-Methoxyphenol	143
3-Methoxytyramine	128
Methyl Benzoate	14, 30, 83, 91, 147
N-Methylbenzylamine	9
Methyldopa	84
o-Methyl Hippuric Acid	128

Methyl p-Hydroxybenzoate	16, 85, 93, 94, 96, 97, 99, 103, 105, 108, 113
Methyl Linoleate	133
Methyl Linolenate	133
Methyl Margarate	133
1-Methylnaphthalene	28, 156
Methyl n-Nonadecanoate	133
Methyl Oleate	133
Methyl Prednisolone	100
Methyl Salicylate	16, 147
Methyl Stearate	133
17-Methyltestosterone	128
Meticrane	105
Metoxuron	122
Metronidazole	105
Mexiletine Hydrochloride	104
Milk	139
4-epi-Minocycline	104
Minocycline Hydrochloride	104
Mitomycin C	104
Mizoribine	104
Monolinuron	122
MPMC	122
MTMC	122
Myoglobin	44, 48, 51, 53, 140
Myosin	138
Myricetin	23
<b>N</b>	
Nabumetone	96
NAC	122
Nalidixic Acid	96
Naphazoline Nitrate	95
Naphthalene	28, 148, 156, 157
2-Naphthalenesulfonic Acid	91
2-Naphthalenesulfonic Acid Sodium	83
1-Naphthol	37
2-Naphthol	37
Naringenin	118
Naringin	116, 117, 118
Neostigmine Methylsulfate	97
Neurotensin	138
New Coccine	123
Niacin	125
Nicardipine Hydrochloride	96
Nicergoline	97
Nicotinamide	96, 125
Nicotinic Acid	93, 96, 125
Nifedipine	96
Nitrendipine	97
p-Nitroaniline	147
Nitrobenzene	147, 155

p-Nitrobenzyl Alcohol	36	2-Phenylethyl- $\beta$ -gentiobioside	19
Nitrofen	121	2-Phenylethyl- $\beta$ -lactoside	19
Nitrofurantoin	112	2-Phenylethyl- $\beta$ -maltoside	19
Nitrofurazone	112	2-Phenylethyl- $\beta$ -melibioside	19
p-Nitrophenol	102	o-Phenylphenol	124
L-Noradrenaline	18, 72	Phenytol	109, 110
Norepinephrine	128	Phloxine	123
<b>O</b>		L- $\alpha$ -Phosphatidyl-L-serine	134
Ofloxacin	112	L- $\alpha$ -Phosphatidyl Choline	134
Oleic Acid	133	L- $\alpha$ -Phosphatidylethanolamine	134
Omeprazole	83	L- $\alpha$ -Phosphatidylinositol Sodium Salt	134
Ovalbumin	44	Phthalic Acid	16, 98
Oxalic Acid	146	Phthalonitrile	26, 156
Oxamic Acid	146	Phytonadione	106
Oxine-copper	14	Pimaricin	99
Oxolinic Acid	112	Piperacillin Sodium	99
Oxyfluorfen	121	Piperonyl Butoxide	124
Oxytetracycline	111	Pirimicarb	122
Oxytocin	16, 137	Piromidic Acid	112
Oxytocin Acetate Salt	83	cis-Platin (CDDP)	111
<b>P</b>		Polyethylene Glycol	49
PAC	122	Poly-L-Glutamic Acid Sodium Salt	45
Paeoniflorin	115	Polyoxyethylene Hexadecyl Ether	148
Paeonol	8, 117	Polyoxyethylene Lauryl Ether	148
PA-Glucose Oligomer (DP=3-22)	132	Potassium Biphthalate	94
Palatinit	131	Potassium Clavulanate	84
Palatinose	130	Pravastatin Sodium	100
Palmatine Chloride	103	Prednisolone	24, 95, 101, 129
Palmatine Chloride Hydrate	113	Prednisolone Acetate	101, 129
Palmatine Hydrochloride	25	Prednisone	98
D-Pantothenic Acid	125	Procainamide	14, 110
D-Pantothenic Acid Sodium Salt	125	Procaine Hydrochloride	101
Papaverine Hydrochloride	120	Procaterol Hydrochloride	101
pd(T)12-18	135	Progesterone	83
Penicillin G Potassium Salt	111	Prometryn	122
Perfluorotetradecanoic Acid	146	Propazine	122
Peroxidase	48, 140	Propionic Acid	18
Perylene	149	Propranolol Hydrochloride	102
PHC	122	Propylbenzene	21, 34, 149, 160
Phenacetin	15, 110, 114	n-Propylbenzene	147
Phenanthrene	148	Propyl Benzoate	100
Phenethylamine Hydrochloride	104	Propyl Gallate	124
Phenol	87, 106, 109, 145, 147	Propyl p-Hydroxybenzoate	16, 88, 89, 91, 92, 97, 100, 107, 113, 124, 147
Phenolphthalein	147	Propyl Salicylate	16
[6,6]-Phenyl-C61 Butyric Acid Butyl Ester	59	Prostaglandin D2	129
[6,6]-Phenyl-C61 Butyric Acid Methyl Ester	59	Prostaglandin E2	129
Phenylalanine	73	Prostaglandin F2 $\alpha$	129
Phenyl Benzoate	89	Prostaglandin I 2	129
2-Phenylethyl- $\beta$ -cellobioside	19	Protocatechuic Acid	16
		PTH-Arg	136
		PTH-Asp	136
		PTH-His	136

PTH-Lys	136
PTH-Ser	136
PTH-Trp	136
Puerarin	113
Pullulan	49
Purpurin	119
Pyrene	148
Pyridine	91, 145, 147
Pyridoxine (Vitamin B6)	125
Pyridoxine Hydrochloride	125
Pyruvate Kinase	140

**Q**

Quercetin	23
Quinidine	109
Quinidine Sulfate	84
Quinine Hydrochloride	84
Quinizarin	119

**R**

Raffinose	41, 130
Ranitidine	110
Reserpine	108, 120
Resibufogenin	115
13-cis-Retinoic Acid	127
all-trans-Retinoic Acid	127
all-trans-Retinol Acetate	123
Rhamnose	41, 42
Riboflavin	100
Riboflavin (Vitamin B2)	125
Riboflavin Phosphate	125
Riboflavin Phosphate Sodium Salt	125
Ribonuclease A	44, 140
Rifampicin	108
Ritodrine Hydrochloride	107, 108
threo-Ritodrine Hydrochloride	107
Rose Bengale	123
Roxatidine Acetate Hydrochloride	109
Roxithromycin	109
Rthoxyquin	124

**S**

Saccharin	123
Saikosaponin a	35, 114, 118
Saikosaponin b1	35, 118
Saikosaponin b2	35, 115, 118
Saikosaponin c	35, 118
Saikosaponin d	35, 114, 118
Saikosaponin h	35
Salicylamide	36, 94

Salicylic Acid	36, 87, 106, 110
Santonin	87
Sc2@C76(I)	57
Sc2@C78	57
Sc2@C80(I)	57
Scopolamine n-Butyl Bromide	99
Scopolamine Hydrobromide	15, 99, 117
Sennoside A	115, 116
Sennoside B	115
6-Shogaol	120
Sitosterol	26
Sodium Amoxicillin	84
Sodium Copper Chlorophyll	119
Sodium Copper Chlorophyllin	123
Sodium 5'-Cytidylate	124
Sodium 5'-Guanylate	124
Sodium 5'-Inosinate	124
Sodium Iron Chlorophyllin	123
Sodium Polystyrene Sulfonate	49
Sodium 5'-Uridylate	124
Sorbic Acid	73, 123, 125
Soybean Milk	139
Stearic Acid	133
Stigmasterol	26
cis-Stilbene	144
trans-Stilbene	144
Streptomycin Sulfate	112
Strychnine	117
Styrene oxide	147
Substance P	16, 137, 138
Succinic Acid	18
Sucrose	41, 43, 130, 131
Sulbactam Sodium	89
Sulfachloropyridazine	112
Sulfadimethoxine	112
Sulfaisomidin	112
Sulfamethazine	112
Sulfamethoxazole	112
Sulfamethoxypyridazine	112
Sulfamonomethoxine	112
Sulfathiazole	82, 112
Sunset Yellow FCF	123
Swertiamarin	116

**T**

Tartaric Acid	18
Tartrazine	123
Taurine	137
Terbutryn	122
Terbutylazine	122
Terephthalonitrile	26, 156
m-Terphenyl	143

o-Terphenyl	12, 21, 143
p-Terphenyl	88, 143
Testosterone	128
Testosterone Propionate	102, 128
Tetracycline	111
Theobromine	31, 109
Theophylline	31, 91, 106, 109, 116
Thermolysin	139
Thiabendazole	124
Thiamine Hydrochloride	91, 114
Thioanisole	147
Thiopental Sodium	91
Thiopental Sodium Isomer	91
Thiram	121
Thymidine	18, 135
Thymine	18, 88, 135
Thyroglobulin	48, 140
Ticarcillin Sodium	92
Tiglic Acid	73
Timolol Maleate	92
Tipepidine	92
Tocopherol	93
α-Tocopherol	27, 93, 127
β-Tocopherol	27, 127
γ-Tocopherol	27, 127
δ-Tocopherol	27, 127
Tocopherol Acetate	93
Tocopherol Nicotinate	93
α-Tocotrienol	126
β-Tocotrienol	126
γ-Tocotrienol	126
δ-Tocotrienol	126
Todralazine Hydrochloride	94, 111
Tolazamide	94
Tolbutamide	94
Tolnaftate	95
Toluene	14, 21, 30, 31, 34, 144, 147, 159, 160
o-Toluic Acid	92, 146
m-Tolunitrile	24, 155
o-Tolunitrile	24, 155
p-Tolunitrile	24, 155
N-(m-Toluoyl) glycine	128
N-(o-Toluoyl) glycine	128
N-(p-Toluoyl) glycine	128
Tranexamic Acid	94
Transferrin	51, 139, 142
Trehalose	43, 130
Triamcinolone	94, 129
Triamcinolone Acetonide	95, 101, 129
Trimetazidine Hydrochloride	95
1-(2,3,4-Trimethoxybenzyl) piperazine Dihydrochloride	95
Trimethylene Glycol	145
Triphenylene	12, 21

Tris (hydroxymethyl) aminomethane	146
Triton X-100	29
Trypsin Inhibitor	51
Trypsin Inhibitor, Soybean	141
Trypsinogen	141
Tryptophan	73
Tyrosine	73

## U

Ubidecarenone	106
Uracil	9, 18, 94, 109, 110, 135, 145, 146, 147
Urea	146
Uric Acid	129
Uridine	18, 135, 145

## V

cis-Vaccenic Acid	133
trans-Vaccenic Acid	133
Valerophenone	12
Valine	39
Vancomycin Hydrochloride	97
Vanillylnonanamide	116
N-Vanillylnonanamide	27
[Arg8]-Vasopressin	83
Vincristine Sulfate	99
Vinyl Acetate	117
2-Vinylpyridine	103
1-Vinyl-2-pyrrolidone	117
Vitamin A	127
Vitamin A Acetate	126
Vitamin B1	125
Vitamin B1 Hydrochloride	125
Vitamin B2	125
Vitamin B6	125
Vitamin B12	125
Vitamin C	125
Vitamin D2	126
Vitamin D3	126
Vitamin E	126
Vitamin E Acetate	126
Vitamin K1	126
Vitamin K3	126

## X

Xanthene	92
Xanthine	129
XMC	122
p-Xylene	36
Xylitol	41, 130, 131
Xylooligosaccharides	130

Xylose	42
--------	----

## Z

Zidovudine	88
------------	----

## Packing material's name

<b>C<sub>18</sub>-MS-II</b>	Benzyl p-Hydroxybenzoate	98
Acesulfame K	Benzylpenicillin Potassium	103
Acetaminophen	Berberine Hydrochloride	25
Acetanilide	Betahistine Mesilate	103
Acetoaminophen	Betamethasone	103, 129
Acetophenone	Betamethasone Sodium Phosphate	103
4'-Acetoxyacetanilide	Betamethasone Valerate	103
17 $\alpha$ -Acetoxyprogesterone	Betamethasone 17-Valerate	172
Acetrizoic Acid	D-Biotin [Vitamin H]	125
N-Acetylprocainamide	Boc-(D)Phe-(L)Phe-OMe	137
Acid Red	Boc-(L)Phe-(L)Phe-OMe	137
Aconitine	BPMC	122
Acrylic Acid	Brilliant Blue FCF	123
Adenosine	Bromacil	122
Alizalin	Bromhexine Hydrochloride	102
Amaranth	Bromobenzene	144, 147, 159
Amidotrizoic Acid	Butylbenzene	12, 21, 34, 160
4-Aminoantipyrine	n-Butylbenzene	147
p-Aminophenol Hydrochloride	Butylbenzoate	12
Amlodipine Besilate	2-tert-Butylhydroquinone	124
Ampicillin	3-tert-Butyl-4-hydroxyanisol	124
Ampicillin Sodium Salt	Butyl p-Hydroxybenzoate	82, 84, 85, 86, 88, 98, 101, 103, 117
D-(-)-Amygdalin	C60	54
Amylbenzene	C70	54
Aniline	Caffeine	23, 31, 101, 105, 107, 115
Anisole	Calcium Folate	104
Anthracene	Campesterol	26
Anthrarufin	Capsaicin	27, 120
Antipyrine	Captopril	83
Apigenin	Carbamazepine	109, 110
L(+)-Ascorbic Acid	Carbenicillin	111
Aspartame	Carbosulfan	122
Asulam	$\alpha$ -Carotene	118
Atrazine	$\beta$ -Carotene	118
Atropine Sulfate	CAT	122
Baicalein	Catechin (C)	23
Baicalin	Catechin gallate(CG)	23
Bamethane Sulfate	Cefaclor	90
Barbital	Cefadroxil	90
Beclometasone Dipropionate	Ceftazidime	91
Bensulide	Ceftriaxone Sodium	91
Benzamide	Cephaeline Hydrobromide	116
Benz[a]anthracene	Cephalothin Sodium Salt	90
Benzene	Chlormethoxynil	121
Benzoic Acid	Chlormadinone Acetate	86
Benzonitrile	Chlorobenzene	144, 147, 159
Benzophenone	4-Chlorobenzoic Acid	82
Benzyl Alcohol	Chloro IPC	122
	m-Chlorophenol	143
	o-Chlorophenol	143

p-Chlorophenol	143	Diphenyl	124
Chlorotoluron	122	Disulfiram	88
Chloroxuron	122	Dobutamine Hydrochloride	94
Chlortetracycline	111	Dopamine Hydrochloride	94
Cholesterol	26	Doxorubicin Hydrochloride	93
Chrysazin	119	Doxycycline Hydrochloride	93
Chrysene	158	Elaidic Acid	133
Chrysin	23	Emetine Hydrochloride	116
Cimetidine	110	Epicatechin(EC)	23
Clarithromycin	85	Epicatechin gallate(EGC)	23
Clindamycin Hydrochloride	85	Epigallocatechin(EGC)	23
Clobetasol Propionate	86	Epigallocatechin gallate(EGCG)	23
Cloxacillin Sodium Monohydrate	85	Erythrosine	123
CNP	121	17- $\alpha$ -Estradiol	128
Coomassie Brilliant Blue G-250	149	17- $\beta$ -Estradiol	128
Cortisone	24, 87, 129	Estradiol Benzoate	83
Cortisone Acetate	87, 101, 129	Estriol	83
Coumarin	119	Ethidimuron	122
Creatinine	128	Ethiofencarb	122
m-Cresol	106	Ethofumesate	122
Crotonic Acid	73	Ethylbenzene	21, 34, 147, 160
Cyanazine	122	p-Ethylbenzoic Acid	146
Daidzein	25	Ethyl p-Hydroxybenzoate	84, 87, 89, 100, 147
Daidzin	25	Ethyl Salicylate	147
Dansyl Amino Acids	136	FAD	100
Daunorubicin Hydrochloride	91	Famotidine	110
DCMU	121	Fast Green FCF	123
DCPA	121	Fisetin	23
Dehydroacetic Acid	123	Flavin Mononucleotide Sodium Salt	125
Desipramine Hydrochloride	109	Flavone	23
Dexamethasone	129	Flumequine	112
Dibenzosuberone	157	Fluocinolone Acetonide	129
Dibenzosuberone	157	Fluocinonide	100
2,6-Di-tert-butyl-4-hydroxytoluene	124	Fluorobenzene	144, 147
Di-n-butyl Phthalate	96	Fluorodifen	121
Diclofenamide	88	Fluorometholone	101, 129
Diethyl Terephthalate	91	Flurbiprofen	111
1,2-Difluorobenzene	144	FMN	100
1,3-Difluorobenzene	144	Folic Acid	104, 105, 125
1,4-Difluorobenzene	144	Fullerene toluene extract	55
Dihydrocoumarin	119	Furan	147
7,8-Dihydroxyflavone	23	Furfural	147
Diketopiperazine	123	Furfuryl Alcohol	147
Diltiazem Hydrochloride	89	Gabexate Mesilate	84
Dimefuron	122	Gallocatechin(GC)	23
Dimethylaminophenol	97	Gallocatechin gallate(GCG)	23
N,N-Dimethylaniline	147	Geniposide	115
N,N-Dimethylbenzylamine	9	Genistein	25
1,5-Dimethylnaphthalene	28, 156	Genistin	25
2,4-Dinitroaniline	147	6-Gingerol	120
1,5-Dinitronaphthalene(1,5-DNN)	28, 156	Ginkgolide A	120
1,8-Dinitronaphthalene(1,8-DNN)	28, 156	Ginkgolide B	120
Diphenhydramine Hydrochloride	110	Ginkgolide C	120

Ginsenoside Rb1	114	Methicillin Sodium Salt	111
Ginsenoside Rc	114	Methomyl	122
Glycitein	25	Methotrexate	105
Glycitin	25	6-Methoxyflavone	23
Griseofulvin	85	m-Methoxyphenol	143
Guaiacol Glycerol Ether	81, 85	o-Methoxyphenol	143
Hibenzic Acid	92	p-Methoxyphenol	143
Hippuric Acid	128	Methyl Benzoate	14, 30, 83, 147
Histidine	73	N-Methylbenzylamine	9
Hydralazine Hydrochloride	111	o-Methyl Hippuric Acid	128
Hydrocortisone	24, 98, 129	Methyl p-Hydroxybenzoate	94, 103, 105
Hydrocortisone Acetate	87, 98	Methyl Linoleate	133
Hydrocortisone Succinate	98	Methyl Linolenate	133
p-Hydroxybenzoic Acid	95	Methyl Margarate	133
3-Hydroxyflavone	23, 119	1-Methylnaphthalene	28, 156
5-Hydroxyflavone	23, 119	Methyl n-Nonadecanoate	133
6-Hydroxyflavone	23, 119	Methyl Oleate	133
7-Hydroxyflavone	23, 119	Methyl Salicylate	147
Ibuprofen	111	Methyl Stearate	133
Imazalil	124	Meticrane	105
Imipramine Hydrochloride	109	Metoxuron	122
Indigo Carmine	123	Metronidazole	105
Indometacin (Indomethacin)	82	Monolinuron	122
Iodobenzene	147	MPMC	122
Iprodione	121	MTMC	122
Isoamyl Benzoate	172	Myricetin	23
Isobutyl p-Hydroxybenzoate	81	NAC	122
Isoniazid	82	Naphthalene	28, 156, 157
Isonicotinic Acid	82	2-Naphthalenesulfonic Acid	91
Isophthalonitrile	26, 156	Naringenin	118
Isoprocarb	122	Naringin	116, 118
Isopropyl p-Hydroxybenzoate	91	Neostigmine Methylsulfate	97
Isoproturon	122	New Coccine	123
Isoxsuprine Hydrochloride	82	Niacin	125
Ketamine Hydrochloride	109	Nicardipine Hydrochloride	96
Ketoprofen	111	Nicergoline	97
Latamoxef Sodium (isomer 1)	106	Nicotinamide	96, 125
Latamoxef Sodium (isomer 2)	106	Nicotinic Acid	93, 96, 125
Lidocaine	109	Nifedipine	96
Linuron	121, 122	Nitrendipine	96
Lisinopril	107	p-Nitroaniline	147
Luteolin	23	Nitrobenzene	147, 155
Mandelic Acid	128	Nitrofen	121
Maneb	121	Ofloxacin	112
MCC	121	Oleic Acid	133
MCP	121	Oxine-copper	14
MCPB	121	Oxolinic Acid	112
MCPP	121	Oxyfluorfen	121
Menatetrenone	106	Oxytetracycline	111
5-Mercapto-1-methyltetrazole	106	Oxytocin Acetate Salt	83
Metabromuron	122	PAC	122
Metazachlor	122	Paeonol	8, 117
Methabenzthiazuron	122	Palmatine Hydrochloride	25

D-Pantothenic Acid Sodium Salt	125	Papaverine Hydrochloride	120
Penicillin G Potassium Salt	111	Penicillin G Potassium Salt	111
PHC	122	Phenacetin	15
Phenacetin	15	Phenol	106, 145, 147
Phenol	106, 145, 147	Phenolphthalein	147
Phenolphthalein	147	Phenylalanine	73
Phenylalanine	73	Phenyl Benzoate	89
Phenyl Benzoate	89	o-Phenylphenol	124
o-Phenylphenol	124	Phenytoin	109, 110
Phenytoin	109, 110	Phloxine	123
Phloxine	123	L- $\alpha$ -Phosphatidyl-L-serine	134
L- $\alpha$ -Phosphatidyl-L-serine	134	L- $\alpha$ -Phosphatidyl Choline	134
L- $\alpha$ -Phosphatidyl Choline	134	L- $\alpha$ -Phosphatidylethanolamine	134
L- $\alpha$ -Phosphatidylethanolamine	134	L- $\alpha$ -Phosphatidylinositol Sodium Salt	134
L- $\alpha$ -Phosphatidylinositol Sodium Salt	134	Phthalonitrile	26, 156
Phthalonitrile	26, 156	Phytonadione	106
Phytonadione	106	Pimaricin	99
Pimaricin	99	Piperacillin Sodium	99
Piperacillin Sodium	99	Pirimicarb	122
Pirimicarb	122	Piromidic Acid	112
Piromidic Acid	112	Potassium Biphthalate	94
Potassium Biphthalate	94	Pravastatin Sodium	100
Pravastatin Sodium	100	Prednisolone	24, 95, 101, 129
Prednisolone	24, 95, 101, 129	Prednisolone Acetate	101, 129
Prednisolone Acetate	101, 129	Procaïnamide	14
Procaïnamide	14	Procaine Hydrochloride	101
Procaine Hydrochloride	101	Progesterone	83
Progesterone	83	Prometryn	122
Prometryn	122	Propazine	122
Propazine	122	Propranolol Hydrochloride	102
Propranolol Hydrochloride	102	Propylbenzene	21, 34, 160
Propylbenzene	21, 34, 160	n-Propylbenzene	147
n-Propylbenzene	147	Propyl Benzoate	100
Propyl Benzoate	100	Propyl Gallate	124
Propyl Gallate	124	Propyl p-Hydroxybenzoate	91, 92, 97, 147
Propyl p-Hydroxybenzoate	91, 92, 97, 147	Prostaglandin D2	129
Prostaglandin D2	129	Prostaglandin E2	129
Prostaglandin E2	129	Prostaglandin F2 $\alpha$	129
Prostaglandin F2 $\alpha$	129	Prostaglandin I2	129
Prostaglandin I2	129	Puerarin	113
Puerarin	113	Purpurin	119
Purpurin	119	Pyridine	91, 145, 147
Pyridine	91, 145, 147	Pyridoxine Hydrochloride	125
Pyridoxine Hydrochloride	125	Quercetin	23
Quercetin	23	Quinidine	109
Quinidine	109	Quinizarin	119
Quinizarin	119	Ranitidine	110
Ranitidine	110	Reserpine	120
Reserpine	120	13-cis-Retinoic Acid	127
13-cis-Retinoic Acid	127	all-trans-Retinoic Acid	127
all-trans-Retinoic Acid	127	Riboflavin	100

Riboflavin Phosphate Sodium Salt	125	Riboflavin Phosphate Sodium Salt	125
Rose Bengale	123	Roxatidine Acetate Hydrochloride	109
Roxatidine Acetate Hydrochloride	109	Rthoxyquin	124
Rthoxyquin	124	Saccharin	123
Saccharin	123	Saikosaponin a	114, 118
Saikosaponin a	114, 118	Saikosaponin b1	118
Saikosaponin b1	118	Saikosaponin b2	115, 118
Saikosaponin b2	115, 118	Saikosaponin c	118
Saikosaponin c	118	Saikosaponin d	114, 118
Saikosaponin d	114, 118	Salicylamide	97
Salicylamide	97	Salicylic Acid	106
Salicylic Acid	106	Santonin	87
Santonin	87	Scopolamine Hydrobromide	15
Scopolamine Hydrobromide	15	Sennoside A	115, 116
Sennoside A	115, 116	Sennoside B	115
Sennoside B	115	6-Shogaol	120
6-Shogaol	120	Sitosterol	26
Sitosterol	26	Sodium Copper Chlorophyll	119
Sodium Copper Chlorophyll	119	Sorbic Acid	73, 123
Sorbic Acid	73, 123	Stearic Acid	133
Stearic Acid	133	Stigmasterol	26
Stigmasterol	26	cis-Stilbene	144
cis-Stilbene	144	trans-Stilbene	144
trans-Stilbene	144	Streptomycin Sulfate	112
Streptomycin Sulfate	112	Styrene oxide	147
Styrene oxide	147	Sulbactam Sodium	89
Sulbactam Sodium	89	Sulfachloropyridazine	112
Sulfachloropyridazine	112	Sulfadimethoxine	112
Sulfadimethoxine	112	Sulfaisomidin	112
Sulfaisomidin	112	Sulfamethazine	112
Sulfamethazine	112	Sulfamethoxazole	112
Sulfamethoxazole	112	Sulfamethoxypyridazine	112
Sulfamethoxypyridazine	112	Sulfamonomethoxine	112
Sulfamonomethoxine	112	Sulfathiazole	112
Sulfathiazole	112	Sunset Yellow FCF	123
Sunset Yellow FCF	123	Swertiamarin	116
Swertiamarin	116	Tartrazine	123
Tartrazine	123	Terbutryn	122
Terbutryn	122	Terbutylazine	122
Terbutylazine	122	Terephthalonitrile	26, 156
Terephthalonitrile	26, 156	m-Terphenyl	143
m-Terphenyl	143	o-Terphenyl	12, 21, 143
o-Terphenyl	12, 21, 143	p-Terphenyl	143
p-Terphenyl	143	Testosterone Propionate	102
Testosterone Propionate	102	Tetracycline	111
Tetracycline	111	Theobromine	31
Theobromine	31	Theophylline	31, 106, 116
Theophylline	31, 106, 116	Thiabendazole	124
Thiabendazole	124	Thioanisole	147
Thioanisole	147	Thiopental Sodium	91
Thiopental Sodium	91	Thiopental Sodium Isomer	91
Thiopental Sodium Isomer	91	Thiram	121
Thiram	121		

Thymine	88
Tiglic Acid	73
Tipepidine	92
Tocopherol	93
α-Tocopherol	27, 93
β-Tocopherol	27
γ-Tocopherol	27
δ-Tocopherol	27
Tocopherol Acetate	93
Tocopherol Nicotinate	93
α-Tocotrienol	126
β-Tocotrienol	126
γ-Tocotrienol	126
δ-Tocotrienol	126
Todalazine Hydrochloride	94, 111
Toluene	14, 21, 30, 31, 34, 144, 147, 159, 160
o-Toluic Acid	146
m-Tolunitrile	24, 155
o-Tolunitrile	24, 155
p-Tolunitrile	24, 155
Triamcinolone	94, 129
Triamcinolone Acetonide	94, 129
Trimetazidine Hydrochloride	95
1-(2,3,4-Trimethoxybenzyl) piperazine Dihydrochloride	95
Triphenylene	21
Triton X-100	29
Tryptophan	73
Tyrosine	73
Uracil	9, 94, 110, 146, 147
cis-Vaccenic Acid	133
trans-Vaccenic Acid	133
Valerophenone	12
Vancomycin Hydrochloride	97
N-Vanillylnonanamide	27
[Arg8]-Vasopressin	83
2-Vinylpyridine	103
Vitamin A Acetate	126
Vitamin B1 Hydrochloride	125
Vitamin B2 [Riboflavin]	125
Vitamin B6	125
Vitamin B12	125
Vitamin C	125
Vitamin D2	126
Vitamin D3	126
Vitamin E	126
Vitamin E Acetate	126
Vitamin K1	126
Vitamin K3	126
Xanthene	92
XMC	122
Zidovudine	88

C<sub>18</sub>-AR-II

Acenaphthene	148
Acenaphthylene	148
Acetaminophen	110
Acetanilide	90
N-Acetylprocainamide	110
Acetylsalicylic Acid (Aspyrine)	110
Acid Red	123
Albiflorin	115
Allopurinol	129
Amaranth	123
4-Aminoacetophenone	98
Aminoantipyrine	110
4-Aminobenzoic Acid	80
4-Amino-6-chlorobenzene-1,3-disulfonamide	98
Amylbenzene	12
4-Androstene-3,17-dione	128
Angiotensin I (Human)	16
Angiotensin II (Human)	16, 138
Angiotensin II, [Asn1,Val5]	138
Angiotensin II, [Sar1,Ala8]	138
Angiotensin II, [Sar1,Ile8]	138
Angiotensin II, [Sar1,Thr8]	138
Angiotensin II, [Val5]	138
Angiotensin II, Des-Asp1-[Ile8]	138
Anthracene	22, 148
Antipyrine	110
Atenolol	81
Azathioprine	80
Aztreonam	80
Baclofen	97
Baicalin	113
Benz[a]anthracene	148
Benzo[b]fluoranthene	148
Benzoic Acid	80
Benzo[g,h,i]perylene	148
Benzophenone	86, 106, 107
Benzo[a]pyrene	148
Benz[a]pyrene	148
Benzyl Alcohol	106
Bezafibrate	102
Bibenzyl	22
Biphenyl	97
Brilliant Blue FCF	123
Bufalin	115
Bufexamac	100
Butylbenzene	12
Butylbenzoate	12
Butyl p-Hydroxybenzoate	16, 84, 108
iso-Butyl p-Hydroxybenzoate	108
Butyl Salicylate	16
Caffeine	93, 96, 109

Camostat Mesilate	84
Carbazole	33
Ceftazidime	90
p-Chlorobenzoic Acid	102
4-Chlorobenzophenone	86
4-Chlorophenol	86
4-(4-Chlorophenyl)-4-hydroxypiperidine	97
Chlorpheniramine Maleate	95
Chlorpropamide	86
Chrysene	148
Cilastatin Sodium Salt	89
trans-Cinnamic Acid	117
Cinobufagin	115
Clofibrate	86
Cloperastine Hydrochloride	86
dAMP	135
dCMP	135
dGMP	135
Cytochrome C	138
Deferoxamine Mesilate	93
Dibenz[a,h]anthracene	148
Dibucaine Hydrochloride	92
Diclofenac Sodium	88
Digitoxigenin	111
Digitoxin	111
3,4-Dihydroxyphenylacetic Acid	128
N,N-Dimethylaniline	146
Diphenylimidazole	100
L-DOPA	128
Dopamine	128
Doxifluridine	93
Epinephrine	128
Erythrosine	123
Estradiol	128
Estriol	128
Estrone	128
Ethinylestradiol	109
o-Ethoxybenzamide	95
Ethyl p-Hydroxybenzoate(4-)	16, 88, 100, 102
2-Ethylpyridine	146
Ethyl Salicylate	16
Famotidine	99
Fast Green FCF	123
Flopropione	102
Fluocinolone Acetonide	100
Fluoranthene	148
Fluorene	148
Flurbiprofen	101
Gallic Acid	16
Geniposide	113
Gentisic Acid	16
Glycyrrhizic Acid	113, 124
Haloperidol	97

Hesperidin	117
Hibenzic Acid	92
Hinokitiol	119
Honokiol	114
Hydrochlorothiazide	98
p-Hydroxybenzoic Acid	87
Hypoxanthine	129
Ibuprofen	86
Indeno[1,2,3-c,d]pyrene	148
Indigo Carmine	123
Indometacin (Indomethacin)	115
Isobutyl p-Hydroxybenzoate	105
Isopropyl p-Hydroxybenzoate	100
Josamycin	89
Ketoprofen	87
Levallorphan Tartrate	108
Levothyroxine Sodium	109
Lidocaine	107
Lysyl Endopeptidase digested	138
Magnolol	114
Mepivacaine Hydrochloride	106
6-Mercaptopurine Hydrate	80
Meropenem Trihydrate	106
Metformin Hydrochloride	105
3-Methoxytyramine	128
Methyl Benzoate	91
Methyl p-Hydroxybenzoate	16, 93, 96, 97, 99, 108, 113
Methyl Salicylate	16
17-Methyltestosterone	128
Nabumetone	96
Nalidixic Acid	96
Naphazoline Nitrate	95
Naphthalene	148
Naringin	117
New Coccine	123
Nicotinic Acid	96
Nitrofurantoin	112
Nitrofurazone	112
p-Nitrophenol	102
Norepinephrine	128
Oxytocin	16
Paeoniflorin	115
Perfluorotetradecanoic Acid	146
Phenacetin	110
Phenanthrene	148
Phenol	87, 109
Phloxine	123
Phthalic Acid	16, 98
Piperonyl Butoxide	124
Poly-L-Glutamic Acid Sodium Salt	45
Procainamide	110
Procaterol Hydrochloride	101
Propyl p-Hydroxybenzoate	16, 88, 100, 113, 124

Propyl Salicylate	16
Protocatechuic Acid	16
PTH-Arg	136
PTH-Asp	136
PTH-His	136
PTH-Lys	136
PTH-Ser	136
PTH-Trp	136
Pyrene	148
Reserpine	108
Resibufogenin	115
Ritodrine Hydrochloride	107, 108
threo-Ritodrine Hydrochloride	107
Rose Bengale	123
Salicylic Acid	87, 110
Sodium Copper Chlorophyllin	123
Sodium 5'-Cytidylate	124
Sodium 5'-Guanylate	124
Sodium 5'-Inosinate	124
Sodium Iron Chlorophyllin	123
Sodium 5'-Uridylate	124
Substance P	16
Sunset Yellow FCF	123
Tartrazine	123
o-Terphenyl	12
Testosterone	128
Testosterone Propionate	128
Theobromine	109
Theophylline	109
Thiamine Hydrochloride	91
Tipepidine	92
Uracil	109
Uric Acid	129
Valerophenone	12
Xanthine	129

C<sub>18</sub>-PAQ

Acenaphthene	87
Acetic Acid	18
Acrinol	80
Adenine	18, 135
Adenosine	18, 135
5'-ADP	135
L-Adrenaline	18, 72
Amikacin Sulfate 2,4,6-trinitrobenzenesulfonic Acid Derivative	81
4-Aminoacetophenone	91
4-Amino-n-butyric Acid [GABA]	137
4-(Aminomethyl) benzoic Acid	94
m-Aminophenol	147
o-Aminophenol	147
p-Aminophenol	147
Amoxicillin	81

5'-AMP	135
Amylbenzene	12
Arbutin	113
Ascorbic Acid	19
Aspartic Acid	137
5'-ATP	135
Atropine Sulfate	116, 117
Barbital Sodium	117
Berberine Chloride	103, 113
Brucine Dihydrate	116, 117
Butylbenzene	12
Butylbenzoate	12
Caffeine	107, 125
Carbazochrome	84
Carbazochrome Sodium Sulfonate	84
Carbidopa	84
Carnitine	125
Cefaclor	90
Cefsulodin Sodium Salt Hydrate	90
Chlorpheniramine Maleate	86
Ciclosporin (Cyclosporin)	88
Cinchonidine	84
Citric Acid	18, 125
Cyanocobalamin	105
Cytidine	18
Cytosine	18
dATP	18
dCTP	18
2'-Deoxyuridine	82
dGTP	18
Digitoxin	87
Digoxin	87
Dihydroquinine	84
Diphenyl Phthalate	95
L-DOPA	18, 72
Dopamine Hydrochloride	18, 72
dTTP	18
Fumaric Acid	18
Gallic Acid	113
Ginsenoside Re	114
Ginsenoside Rg1	114
Gitoxin	87
Glutamic Acid	137
Glycine	38
Glycolic Acid	18
Glycylglycine	38
Guanine	18
Guanosine	18
Homovanillic Acid	18
Hydroquinidine Hydrochloride	84
Hydroquinone	113
Hydroxocobalamin Acetate	105
Idoxuridine	92

5-Iodouracil	92
Isonicotinamide	90
Kanamycin Sulfate 2,4,6-trinitrobenzenesulfonic Acid Derivative	81
Lactic Acid	18
Latamoxef Sodium (isomer 1)	107
Latamoxef Sodium (isomer 2)	107
Liothyronine Sodium	107
Lisinopril	107
Maleic Acid	18, 86
Malic Acid	18
Mecobalamin	105
DL-3-Methoxy-4-hydroxymandelic Acid	18
Methyldopa	84
Mizoribine	104
Nicotinamide	125
L-Noradrenaline	18, 72
PA-Glucose Oligomer (DP=3-22)	132
Palmitine Chloride	103
Palmitine Chloride Hydrate	113
Phenacetin	114
2-Phenylethyl-β-cellobioside	19
2-Phenylethyl-β-gentiobioside	19
2-Phenylethyl-β-lactoside	19
2-Phenylethyl-β-maltoside	19
2-Phenylethyl-β-melibioside	19
Potassium Clavulanate	84
Propionic Acid	18
Propyl p-Hydroxybenzoate	107
Quinidine Sulfate	84
Quinine Hydrochloride	84
Riboflavin Phosphate	125
Roxithromycin	109
Scopolamine Hydrobromide	117
Senenoside A	116
Sodium Amoxicillin	84
Strychnine	117
Succinic Acid	18
Sulfathiazole	82
Tartaric Acid	18
o-Terphenyl	12
Thiamine Hydrochloride	114
Thymidine	18
Thymine	18
Tolnaftate	95
Tranexamic Acid	94
Triphenylene	12
Ubidecarenone	106
Uracil	18
Uridine	18
Valerophenone	12
Vitamin B1	125
Vitamin B6	125

## Cholester

Amylbenzene	21, 160
Angiotensin I (Human)	137
Angiotensin II (Human)	137
Anthracene	22
Apigenin	118
Baicalein	23
Benz[a]anthracene	158
Benzene	21, 160
Bibenzyl	22
Bromobenzene	144
Butylbenzene	21, 160
Caffeine	23, 35
α-Carotene	118
β-Carotene	118
Catechin (C)	23, 35
Catechin gallate(CG)	23, 35
Chlorobenzene	144
m-Chlorophenol	143
o-Chlorophenol	143
p-Chlorophenol	143
Chondroitinase AC-II digested	132
Chondroitin Sulfate A	132
Chondroitin Sulfate C	132
Chrysene	158
Chrysin	23
7,8-Dihydroxyflavone	23
Elaidic Acid	133
Epicatechin(EC)	23, 35
Epicatechin gallate(EGC)	23, 35
Epigallocatechin(EGC)	23, 35
Epigallocatechin gallate(EGCG)	23, 35
Ethylbenzene	21, 160
Fisetin	23
Flavone	23
Gallocatechin(GC)	23, 35
Gallocatechin gallate(GCG)	23, 35
3-Hydroxyflavone	23
5-Hydroxyflavone	23
6-Hydroxyflavone	23
7-Hydroxyflavone	23
Linolenic Acid p-Bromophenacyl Ester	134
γ-Linolenic Acid p-Bromophenacyl Ester	134
Luteolin	23
6-Methoxyflavone	23
m-Methoxyphenol	143
o-Methoxyphenol	143
p-Methoxyphenol	143
Methyl Linoleate	133
Methyl Linolenate	133
Methyl Margarate	133
Methyl n-Nonadecanoate	133

Methyl Oleate	133
Methyl Stearate	133
Myricetin	23
Naringenin	118
Naringin	118
Oleic Acid	133
Oxytocin	137
Phenol	145
Propylbenzene	21, 160
Pyridine	145
Quercetin	23
13-cis-Retinoic Acid	127
all-trans-Retinoic Acid	127
Saikosaponin a	35, 118
Saikosaponin b1	35, 118
Saikosaponin b2	35, 118
Saikosaponin c	35, 118
Saikosaponin d	35, 118
Saikosaponin h	35
Stearic Acid	133
cis-Stilbene	144
trans-Stilbene	144
Substance P	137
m-Terphenyl	143
o-Terphenyl	21, 143
p-Terphenyl	143
Toluene	21, 144, 160
Triphenylene	21
cis-Vaccenic Acid	133
trans-Vaccenic Acid	133
Vitamin D2	126
Vitamin D3	126
<b><math>\pi</math>NAP</b>	
Adenine	135
Adenosine	135
Anisole	155
Benzene	144, 155
Berberine Hydrochloride	25
Betamethasone	129
Campesterol	26
Capsaicin	27, 120
Cholesterol	26
Cortisone	24, 129
Cortisone Acetate	129
Cytidine	129
Cytosine	129
Daidzein	25
Daidzin	25
Dexamethasone	129
1,2-Difluorobenzene	144
1,3-Difluorobenzene	144

1,4-Difluorobenzene	144
Fluocinolone Acetonide	129
Fluorobenzene	144
Fluorometholone	129
Genistein	25
Genistin	25
Glycitein	25
Glycitin	25
Guanine	135
Guanosine	135
Hydrocortisone	24, 129
Isophthalonitrile	26
Nitrobenzene	155
Palmitate Hydrochloride	25
Phthalonitrile	26
Prednisolone	24, 129
Prednisolone Acetate	129
Sitosterol	26
Stigmasterol	26
Terephthalonitrile	26
Thymidine	135
Thymine	135
$\alpha$ -Tocopherol	27
$\beta$ -Tocopherol	27
$\gamma$ -Tocopherol	27
$\delta$ -Tocopherol	27
$\alpha$ -Tocotrienol	126
$\beta$ -Tocotrienol	126
$\gamma$ -Tocotrienol	126
$\delta$ -Tocotrienol	126
m-Tolunitrile	24
o-Tolunitrile	24
p-Tolunitrile	24
Triamcinolone	129
Triamcinolone Acetonide	129
Uracil	135
Uridine	135
N-VanillylNonanamide	27

**PYE**

Amylbenzene	160
Anisole	155
Anthracene	157
Benz[a]anthracene	158
Benzene	155, 157, 160
2-Benzylpyridine	148
3-Benzylpyridine	148
4-Benzylpyridine	148
Bromobenzene	159
Butylbenzene	160
C60	54
C70	54

Campesterol	26
Capsaicin	27
Carvylacetate (cis, trans)	148
Chlorobenzene	159
Cholesterol	26
Chrysene	158
Dibenzosuberone	157
Dibenzosuberone	157
1,5-Dimethylnaphthalene	28, 156
1,5-Dinitronaphthalene(1,5-DNN)	28, 156
1,8-Dinitronaphthalene(1,8-DNN)	28, 156
Ethylbenzene	160
Isophthalonitrile	26, 156
1-Methylnaphthalene	28, 156
Naphthalene	28, 156, 157
Nitrobenzene	155
Phthalonitrile	26, 156
Propylbenzene	160
Sitosterol	26
Stigmasterol	26
Terephthalonitrile	26, 156
$\alpha$ -Tocopherol	27
$\beta$ -Tocopherol	27
$\gamma$ -Tocopherol	27
$\delta$ -Tocopherol	27
Toluene	159, 160
m-Tolunitrile	155
o-Tolunitrile	155
p-Tolunitrile	155
N-(m-Toluoyl) glycine	128
N-(o-Toluoyl) glycine	128
N-(p-Toluoyl) glycine	128
N-VanillylNonanamide	27

**NPE**

Amylbenzene	160
Benzene	144, 160
Butylbenzene	160
C60	54, 59
C70	54
1,2-Difluorobenzene	144
1,3-Difluorobenzene	144
1,4-Difluorobenzene	144
1,5-Dimethylnaphthalene	28, 156
1,5-Dinitronaphthalene(1,5-DNN)	28, 156
1,8-Dinitronaphthalene(1,8-DNN)	28, 156
Ethylbenzene	160
Fluorobenzene	144
Isophthalonitrile	156
1-Methylnaphthalene	28, 156
Naphthalene	28, 156
[6,6]-Phenyl-C61 Butyric Acid Butyl Ester	59

[6,6]-Phenyl-C61 Butyric Acid Methyl Ester	59
Phthalonitrile	156
Propylbenzene	160
Terephthalonitrile	156
Toluene	160
m-Tolunitrile	155
o-Tolunitrile	155
p-Tolunitrile	155

**PBB-R**

Allylbenzene	149
Amylbenzene	160
Anthracene	157
Benzene	157, 160
1,1'-Binaphthyl	149
Bromobenzene	159
Butylbenzene	160
Chlorobenzene	159
Dibenzosuberone	149, 157
Dibenzosuberone	149, 157
Diphenylmethane	149
Ethylbenzene	160
Fluorene	149
Naphthalene	157
Perylene	149
Polyoxyethylene Hexadecyl Ether	148
Polyoxyethylene Lauryl Ether	148
Propylbenzene	149, 160
Toluene	159, 160
Triton X-100	29

**CN-MS**

Acetaminophen	31
Acetophenone	30
4'-Acetoxyacetanilide	31
p-Aminophenol Hydrochloride	31
Benzene	30, 31
Caffeine	31
4-Chlorophenol	85
Clofibrate	85
4-Ethoxyphenol	85
Methyl Benzoate	30
Theobromine	31
Theophylline	31
Toluene	30, 31

**C<sub>22</sub>-AR-II**

Acetophenone	30
Benzene	30
Methyl Benzoate	30

Toluene	30
<b>C<sub>8</sub>-MS</b>	
Acetophenone	30
Benzene	30
Bisoprolol	98
Butyl p-Hydroxybenzoate	108
Cetirizine Dihydrochloride	89
Clindamycin Phosphate	85
Cyanocobalamin	87
Cyanocobalamin Derivative	87
1,2-Dinitrobenzene	83
Dipyridamole	88
Domperidone	95
Ethyl p-Hydroxybenzoate	95
Fumaric Acid	98
Ipratropium Bromide	82
Ipratropium Bromide Derivative	82
Lincomycin Hydrochloride	85, 108
Methyl Benzoate	30
Methyl p-Hydroxybenzoate	85
Mexiletine Hydrochloride	104
4-epi-Minocycline	104
Minocycline Hydrochloride	104
Omeprazole	83
Phenethylamine Hydrochloride	104
Propyl p-Hydroxybenzoate	89
Rifampicin	108
Scopolamine n-Butyl Bromide	99
Scopolamine Hydrobromide	99
p-Terphenyl	88
Toluene	30
Vincristine Sulfate	99
Vinyl Acetate	117
1-Vinyl-2-pyrrolidone	117

**C<sub>4</sub>-MS**

Acetophenone	30
Benzene	30
Methyl Benzoate	30
Toluene	30

**TMS-MS**

Acetophenone	30
Benzene	30
Daunorubicin Hydrochloride	83
Doxorubicin Hydrochloride	83
Epirubicin Hydrochloride	83
Methyl Benzoate	30
2-Naphthalenesulfonic Acid Sodium	83

Ticarcillin Sodium	92
Toluene	30
o-Toluic Acid	92

**PE-MS**

Acetophenone	30
Anisole	155
Benz[a]anthracene	158
Benzene	30, 155
Berberine Hydrochloride	25
Capsaicin	116, 120
Chrysene	158
Cortisone	24
Daidzein	25
Daidzin	25
Dexamethasone	92
Ethyl Vanillin	104
Genistein	25
Genistin	25
Glycitein	25
Glycitin	25
Hydrocortisone	24
Maleic Acid	92
Methyl Benzoate	30
Mitomycin C	104
Nitrobenzene	155
Palmatine Hydrochloride	25
Prednisolone	24
Timolol Maleate	92
Toluene	30
m-Tolunitrile	24, 155
o-Tolunitrile	24, 155
p-Tolunitrile	24, 155
Vanillylnonanamide	116

**SL-II**

Acetaminophen	110
Acetanilide	110
Benzoic Acid	36
β-Carotene	123, 127
Cetirizine Dihydrochloride	89
Cinnamyl Alcohol	36
Dantrolene Sodium	91
4-Dimethylaminoantipyrine	89
Fluocinolone Acetonide	101
Fluoxymesterone	100
Hydrocortisone	98
Methyl Prednisolone	100
1-Naphthol	37
2-Naphthol	37
p-Nitrobenzyl Alcohol	36

Phenacetin	110
L-α-Phosphatidyl Choline	134
Prednisone	98
all-trans-Retinol Acetate	123
Salicylamide	36
Salicylic Acid	36
Theophylline	91
α-Tocopherol	127
β-Tocopherol	127
γ-Tocopherol	127
δ-Tocopherol	127
Tolazamide	94
Tolbutamide	94
Triamcinolone Acetonide	101
Vitamin A	127
p-Xylene	36

**HILIC**

Allantoic Acid	112
Allantoin	112
Ammelide	39
Ammeline	39
Angiotensin II, [Asn1,Val5]	138
Angiotensin II, [Sar1,Ile8]	138
Angiotensin II, [Val5]	138
Angiotensin II, Des-Asp1-[Ile8]	138
Ascorbic Acid	125, 126
L(+)-Ascorbic Acid	125
L-Citrulline	120
Cyanuric Acid	39
Diethylene Glycol	145
meso-Erythritol	146
Ethylene Glycol	145
D-Fructose-6-phosphate	132
Glucose	132
α-D-Glucose-1-phosphate	132
D-Glucose-6-phosphate	132
Glyceric Acid	146
Glycerol	145
Glycine	38
Glycylglycine	38
Isoascorbic Acid	125
Isoleucine	39
Leucine	39
Malic Acid	120, 126
Melamine	39
Nicotinamide	125
Nicotinic Acid	125
Oxalic Acid	146
Oxamic Acid	146
D-Pantothenic Acid	125
Pyridoxine (Vitamin B6)	125

Riboflavin (Vitamin B2)	125
Sorbic Acid	125
Taurine	137
Trimethylene Glycol	145
Tris (hydroxymethyl) aminomethane	146
Uracil	145
Urea	146
Uridine	145
Valine	39

**Sugar-D**

Arabinose	41
α-Cyclodextrin	130
β-Cyclodextrin	130
γ-Cyclodextrin	130
meso-Erythritol	41, 130, 131
Fructooligosaccharides	42
Fructose	41, 42, 130, 131
Galactose	41
Glucitol (Sorbitol)	41, 131
Glucose	40, 41, 42, 130, 131
Glucuronic Acid	42
Glycerol	41
Guanosine	111
Inositol	41, 131
Isomaltulooligosaccharides	130
Lactose	43, 131
Maltitol	41, 130, 131
Maltoheptaose	42
Maltohexaose	42
Maltopentaose	42
Maltose	40, 41, 42, 43, 130, 131
Maltotetraitol	131
Maltotetraose	42
Maltotriitol	131
Maltotriose	42, 131
Mannitol	131
Mannosamine	42
Mannose	42
Palatinit	131
Palatinose	130
cis-Platin (CDDP)	111
Raffinose	41, 130
Rhamnose	41, 42
Sucrose	41, 43, 130, 131
Trehalose	43, 130
Xylitol	41, 130, 131
Xylooligosaccharides	130
Xylose	42

NH<sub>2</sub>-MS

Arabinose	41
Fructose	42
Galactose	41
Glucose	40, 42
Lactose	43
Maltitol	41
Maltose	40, 43
Mannose	42
Rhamnose	42
Sucrose	43
Trehalose	43
Xylose	42

## Protein-R

Actinase E	139
Albumin (Bovine)	44, 139
Albumin (Goat)	140
Albumin (Human)	139
Albumin (Ovalbumin)	140
Albumin (Rabbit)	140
Alcohol Dehydrogenase	139
α-Amylase	139
Angiotensin I (Human)	137
Angiotensin II (Human)	137
Carbonic Anhydrase	139
Catalase	140
Choline Oxidase	139
α-Chymotrypsinogen A	140
Conalbumin	140
Concanavalin A	140, 141
Cytochrome C	44, 140
Deoxyribonuclease I	140
Diaphorase	141
Fibrinogen	140
Glucose Oxidase	139
L-Glutamic Dehydrogenase	139
Lysozyme	44, 140
Milk	139
Myoglobin	44, 140
Myosin	138
Ovalbumin	44
Oxytocin	137
pd(T)12-18	135
Peroxidase	140
Poly-L-Glutamic Acid Sodium Salt	45
Pyruvate Kinase	140
Ribonuclease A	44, 140
Soybean Milk	139
Substance P	137
Thermolysin	139

Thyroglobulin	140
Transferrin	139

C<sub>18</sub>-AR-300

[Ala2]-Leu-Enkephalin	138
[Ala2]-Met-Enkephalin	138
Angiotensin II	138
Bombesin	138
Bradykinin	138
Carbonic Anhydrase	141
α-Chymotrypsinogen A	141
Hemoglobin, Bovine	46
Insulin	141
Leu-Enkephalin	138
Met-Enkephalin	138
Neurotensin	138
Substance P	138
Transferrin	141
Trypsin Inhibitor, Soybean	141
Trypsinogen	141

C<sub>8</sub>-AR-300

Carbonic Anhydrase	141
α-Chymotrypsinogen A	141
Hemoglobin, Bovine	46
Insulin	141
Transferrin	141
Trypsin Inhibitor, Soybean	141
Trypsinogen	141

C<sub>4</sub>-AR-300

Carbonic Anhydrase	141
α-Chymotrypsinogen A	141
Hemoglobin, Bovine	46
Insulin	141
Transferrin	141
Trypsin Inhibitor, Soybean	141
Trypsinogen	141

## Ph-AR-300

Carbonic Anhydrase	141
α-Chymotrypsinogen A	141
Hemoglobin, Bovine	46
Insulin	141
Transferrin	141
Trypsin Inhibitor, Soybean	141
Trypsinogen	141

## Diol-120-II

Aprotinin	48
Bovine Serum	142
Conalbumin	48
Egg White	142
Glucose Oxidase	48
Human Serum	142
Myoglobin	48
Peroxidase	48
Polyethylene Glycol	49
Thyroglobulin	48

## Diol-300-II

Aprotinin	48
Bovine serum	142
Conalbumin	48
E.coli proteins	49
Egg White	142
Glucose Oxidase	48
Human Serum	142
Myoglobin	48
Peroxidase	48
Pullulan	49
Sodium Polystyrene Sulfonate	49
Thyroglobulin	48

## DEAE

Lyxoytase	51
Myoglobin	51
Transferrin	51
Trypsin Inhibitor	51

## QA

Lyxoytase	51
-----------	----

## CM

Lyxoytase	51
-----------	----

## SP

Lyxoytase	51
-----------	----

## HIC

Albumin (Bovine)	53
β-Glucosidase	53
Hemoglobin	53
β-Lactoglobulin	53

Myoglobin	53
-----------	----

## Buckyprep

C60	54, 56, 59
C60O	56
C70	54, 56
C76	56
C84	56
C86	57
Crude fullerenes	58
Fullerene toluene extract	55
Gd@C82(I)	57
[6,6]-Phenyl-C61 Butyric Acid Butyl Ester	59
[6,6]-Phenyl-C61 Butyric Acid Methyl Ester	59
Sc2@C76(I)	57
Sc2@C78	57
Sc2@C80(I)	57

## Buckyprep-M

C60	54
C70	54
C86	57
Gd@C82(I)	57
Sc2@C76(I)	57
Sc2@C78	57
Sc2@C80(I)	57

## PBB

C60	54, 58, 59
C70	54, 58
C76	58
C84	58
Crude fullerenes	58
[6,6]-Phenyl-C61 Butyric Acid Butyl Ester	59
[6,6]-Phenyl-C61 Butyric Acid Methyl Ester	59

## PBB-R

Allylbenzene	149
1,1'-Binaphthyl	149
Dibenzosuberenone	149
Dibenzosuberone	149
Diphenylmethane	149
Fluorene	149
Perylene	149
Polyoxyethylene Hexadecyl Ether	148
Polyoxyethylene Lauryl Ether	148
Propylbenzene	149

I. COSMOSIL HPLC column	<p><b>CNT-300</b></p> <p>DNA-wrapped CNT ..... 60</p>
II. Liquid chromatography related product	<p><b>CNT-1000</b></p> <p>DNA ..... 135</p> <p>DNA-wrapped CNT ..... 60</p> <p>Hyaluronic Acid ..... 132</p> <p><b>CNT-2000</b></p> <p>DNA ..... 135</p> <p>DNA-wrapped CNT ..... 60</p>
III. Application data	
IV. Technical note	
V. Index	

**Warranties and Disclaimers:**

Nacalai Tesque warrants that its products shall conform to the description of such products as provided by Nacalai Tesque through its catalog, analytical data or other literature. Nacalai Tesque makes no other warranty, express or implied, as to the fitness of these products for any particular purpose. Nacalai Tesque shall not in any event be liable for incidental or consequential damages that may result from any use or failure of the products.

For more information on products and pricing, please contact your local distributor.

---

**NACALAI TESQUE, INC.**

Nijo Karasuma, Nakagyo-ku  
Kyoto 604-0855 JAPAN  
TEL :+81-(0)75-251-1730  
FAX :+81-(0)75-251-1763  
Web site :<http://www.nacalai.com>  
E-mail :[info.intl@nacalai.com](mailto:info.intl@nacalai.com)